

TECHNICAL BULLETIN

PowerPlex[®] 16 and PowerPlex[®] ES Monoplex Systems

Instructions for Use of Products

DC6551, DC6561, DC6571, DC6581, DC6591, DC6601, DC6611,
DC6621, DC6631, DC6641, DC6651, DC6661, DC6671, DC6681,
DC6691 and DC6751



PowerPlex® 16 and PowerPlex® ES Monoplex Systems

All technical literature is available at: www.promega.com/protocols/
 Visit the web site to verify that you are using the most current version of this Technical Bulletin.
 E-mail Promega Technical Services if you have questions on use of this system: techserv@promega.com

1. Description.....	1
2. Product Components and Storage Conditions	2
3. Amplification Protocols	3
4. Detection Methods	4
4.A. Detection Using the ABI PRISM® 310, 3100 or 3100- <i>Avant</i> Genetic Analyzer and the Applied Biosystems 3130 or 3130 <i>xl</i> Genetic Analyzer.....	4
4.B. Detection Using the Hitachi FMBIO® II Fluorescence Imaging System.....	5
5. Summary of Changes	6

1. Description

The PowerPlex® 16 and PowerPlex® ES Monoplex Systems contain primer pairs that have the same sequence as those used in the PowerPlex® 16 (Cat.# DC6531), PowerPlex® 16 BIO (Cat.# DC6541), PowerPlex® ESX 16 Fast System (Cat.# DC1611), PowerPlex® ESI 16 Fast System (Cat.# DC1621), PowerPlex® ESX 17 Fast System (Cat.# DC1711) and PowerPlex® ESI 17 Fast System (Cat.# DC1721). Refer to the *PowerPlex® 16 System Technical Manual #TMD012*, for locus-specific information, allele size range, allelic ladder components and genotypes of commonly used standard DNA templates. Refer to the *PowerPlex® ESX 17 Fast System Technical Manual #TMD040* for SE33 locus-specific information. Technical Manuals are available upon request or at: www.promega.com/protocols/

The PowerPlex® 16 and PowerPlex® ES Monoplex Systems were developed for human identification applications including forensic analysis, relationship testing and research use.



2. Product Components and Storage Conditions

PRODUCT	SIZE	CAT.#
PowerPlex® 16 Monoplex System, Penta E (Fluorescein) ^(a,b)	100 reactions	DC6591
PowerPlex® 16 Monoplex System, Penta D (JOE) ^(a,b)	100 reactions	DC6651
PowerPlex® ES Monoplex System, SE33 (JOE) ^(a)	100 reactions	DC6751

Not For Medical Diagnostic Use. Each system includes:

- 150µl Internal Lane Standard
- 300µl Gold ST★R Buffer
- 250µl 10X Primer Pair Mix
- 70µl Allelic Ladder

PRODUCT*	SIZE	CAT.#
PowerPlex® 16 Monoplex System, D3S1358 (Fluorescein) ^(a)	100 reactions	DC6551
PowerPlex® 16 Monoplex System, TH01 (Fluorescein) ^(a)	100 reactions	DC6561
PowerPlex® 16 Monoplex System, D21S11 (Fluorescein) ^(a)	100 reactions	DC6571
PowerPlex® 16 Monoplex System, D18S51 (Fluorescein) ^(a)	100 reactions	DC6581
PowerPlex® 16 Monoplex System, D5S818 (JOE) ^(a)	100 reactions	DC6601
PowerPlex® 16 Monoplex System, D13S317 (JOE) ^(a)	100 reactions	DC6611
PowerPlex® 16 Monoplex System, D7S820 (JOE) ^(a)	100 reactions	DC6621
PowerPlex® 16 Monoplex System, D16S539 (JOE) ^(a)	100 reactions	DC6631
PowerPlex® 16 Monoplex System, CSF1PO (JOE) ^(a)	100 reactions	DC6641
PowerPlex® 16 Monoplex System, vWA (TMR) ^(a)	100 reactions	DC6661
PowerPlex® 16 Monoplex System, D8S1179 (TMR) ^(a)	100 reactions	DC6671
PowerPlex® 16 Monoplex System, TPOX (TMR) ^(a)	100 reactions	DC6681
PowerPlex® 16 Monoplex System, FGA (TMR) ^(a)	100 reactions	DC6691

Not For Medical Diagnostic Use. Each system includes:

- 150µl Internal Lane Standard
- 300µl Gold ST★R Buffer
- 250µl 10X Primer Pair Mix

Storage Conditions: Store all components at –20°C. The fluorescent 10X Primer Pair Mix is light-sensitive; therefore minimize light exposure and store in the dark.

*Items listed are available from Promega as “special order” items. A minimum order is not required, but items may not be immediately available for shipment. Contact Promega Customer Service for more information.

3. Amplification Protocols

Follow the protocols in the *PowerPlex® 16 System Technical Manual #TMD012* for the PowerPlex® 16 Monoplex Systems and in the *PowerPlex® ESX 17 Fast System Technical Manual #TMD040* for the PowerPlex® ES Monoplex System. However, follow the directions below for the amount of AmpliTaq Gold® DNA polymerase to use per 25µl reaction. We highly recommend the use of gloves and aerosol-resistant pipette tips.

1. Prepare the amplification mix as directed in Tables 1 and 2.
2. Add template DNA as directed in Tables 1 and 2. For optimal amplification results, we recommend 0.5–2ng of DNA per reaction.

Table 1. Amplification Mix for a Single Sample.

PCR Component	Volume Per Sample
Nuclease-Free Water	17.4µl
Gold ST★R 10X Buffer	2.5µl
10X Primer Pair	2.5µl
AmpliTaq Gold® DNA polymerase ¹	0.1µl (0.5 units)
Total amplification mix volume	22.5µl
Template DNA ² (to be added)	2.5µl
Total Reaction Volume	25.0µl

¹Assumes the AmpliTaq Gold® DNA polymerase is 5u/µl. If the enzyme concentration is different, adjust the volume of enzyme used accordingly.

²Assumes the template DNA volume is 2.5µl.

Table 2. Preparation of Amplification Mix for Multiple Samples.

PCR Component	Volume Per Sample	× Number of Reactions	= Final Volume (µl)
Nuclease-Free Water	µl		
Gold ST★R 10X Buffer	2.5µl		
10X Primer Pair	2.5µl		
AmpliTaq Gold®			
DNA polymerase ¹	0.1µl (0.5 units)		
Total amplification mix volume²	µl		
Template DNA (to be added)	µl		
Total Reaction Volume	25.0µl		

¹Assumes the AmpliTaq Gold® DNA polymerase is 5u/µl. If the enzyme concentration is different, adjust the volume of enzyme used accordingly.

²The volume of template DNA plus the volume of amplification mix should equal 25.0µl.



4. Detection Methods

Note: Allelic ladders are only included with the PowerPlex® 16 Monoplex System, Penta E (Fluorescein), PowerPlex® 16 Monoplex System, Penta D (JOE) and PowerPlex® ES Monoplex System, SE33 (JOE). Allelic ladders for the other systems are available by custom order. Allelic ladder options are included in Tables 3 and 4.

4.A. Detection Using the ABI PRISM® 310, 3100 or 3100-Avant Genetic Analyzer and the Applied Biosystems 3130 or 3130xl Genetic Analyzer

Specific directions for detection using the ABI PRISM® 310, 3100 or 3100-Avant Genetic Analyzers and the Applied Biosystems 3130 or 3130xl are given in the *PowerPlex® 16 System Technical Manual* #TMD012 and the *PowerPlex® ESX 17 Fast System Technical Manual* #TMD040. Please follow these protocols when using the PowerPlex® Monoplex Systems. For additional questions, contact Promega Technical Services.

Table 3. Allelic Ladder Options for the ABI PRISM® 310, 3100 and 3100-Avant Genetic Analyzers and the Applied Biosystems 3130 or 3130xl Genetic Analyzer.

PowerPlex® 16/ES Monoplex Loci	Allelic Ladder Options ^{1,2}
D18S51, D21S11, TH01, D3S1358, FGA, D8S1179, vWA	PowerPlex® 16
Penta E	Penta E or PowerPlex® 16
Penta D	Penta D or PowerPlex® 16
CSF1PO, D16S539, D7S820, D13S317, D5S818, TPOX	PowerPlex® 16
SE33	SE33

¹The PowerTyper™ Macros (Cat.# DG3470) can be used for data analysis.

²All allelic ladders are available by custom order.

4.B. Detection Using the Hitachi FMBIO® II Fluorescence Imaging System

For detection using the Hitachi FMBIO® II Fluorescence Imaging System, see the *PowerPlex® 16 BIO System Technical Manual #TMD016*. The filter set used for detection of the PowerPlex® 16 and PowerPlex® ES Monoplex loci is shown in Table 4.

Table 4. Allelic Ladder Options for the Hitachi FMBIO® II Fluorescence Imaging System. The filter set used for detection of each PowerPlex® 16 or PowerPlex® ES Monoplex locus is shown in bold in the “Filters Needed for Detection” Column.

PowerPlex® 16/ES Monoplex Locus	Allelic Ladder Options ¹	Filters Needed for Detection (nm)	PowerPlex® Filter Set
D18S51, D21S11, TH01, D3S1358	PowerPlex® 2.1	505 , 585, 650	1.1/2.1
	PowerPlex® 16 BIO ²	505 , 577, 598, 665	16 BIO
CSF1PO, D16S539, D7S820, D13S317, D5S818	PowerPlex® 16 BIO ²	505, 577 , 598, 665	16 BIO
	JOE Ladder Mix	585 , 650 577 , 665	1.1/2.1 16 BIO
FGA, TPOX, D8S1179, vWA	PowerPlex® 2.1	505, 585 , 650	1.1/2.1
		505, 577 , 665	16 BIO ³
SE33	SE33 Allelic Ladder	585 , 650	1.1/2.1
		577 , 665	16 BIO
Penta E	PowerPlex® 2.1	505 , 585, 650	1.1/2.1
	PowerPlex® BIO	505 , 577 598, 665	16 BIO
	Penta E Allelic Ladder	505 , 650	1.1/2.1
		505 , 665	16 BIO
Penta D	PowerPlex® BIO	505, 577 598, 665	16 BIO
	JOE Ladder Mix	585 , 650	1.1/2.1
		577 , 665	16 BIO
	Penta D Allelic Ladder	585 , 650	1.1/2.1
577 , 665		16 BIO	

¹All allelic ladders are available by custom order.

²The Matrix 16 BIO is required for use with the PowerPlex® 16 BIO Allelic Ladder Mix. The Internal Lane Standard 600, labeled with CXR, can be detected with the 665nm filter.

³When using the PowerPlex® 16 BIO filter set for TMR-labeled loci, do not include the data from the 598nm filter in the FMBIO® Analysis Project or do not scan using the 598nm filter. Both TMR and CXR will be detected with the 598nm filter making color separation impossible. The TMR-labeled loci will be detected in the 577nm scan (JOE channel).



5. Summary of Changes

The following changes were made to the 8/15 revision of this document:

1. Replaced mention of discontinued products with current products.
2. Removed expired disclaimers.
3. Moved document into new format.

^(a) The purchase of this product does not convey a license to use AmpliTaq Gold[®] DNA polymerase. You should purchase AmpliTaq Gold[®] DNA polymerase licensed for the forensic and human identity field directly from your authorized enzyme supplier.

^(b) U.S. Pat. No. 6,238,863, Chinese Pat. No. ZL99802696.4, European Pat. No. 1058727, Japanese Pat. No. 4494630 and other patents pending.

© 2002, 2003, 2007, 2011, 2015 Promega Corporation. All Rights Reserved.

PowerPlex is a registered trademark of Promega Corporation. PowerTyper is a trademark of Promega Corporation.

ABI PRISM is a registered trademark of Applied Biosystems Corporation. AmpliTaq Gold is a registered trademark of Roche Molecular Systems, Inc. FMBIO is a registered trademark of Hitachi Software Engineering Company, Ltd

Products may be covered by pending or issued patents or may have certain limitations. Please visit our Web site for more information.

All prices and specifications are subject to change without prior notice.

Product claims are subject to change. Please contact Promega Technical Services or access the Promega online catalog for the most up-to-date information on Promega products.