# Mupid-exu Electrophoresis System

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### Electrophoresis Tank

Overall dimension	183mm(W) × 56mm(H) × 164mm(L)
Material characteristic	UV transparent (50% at 254nm, 80% at 312nm)
Solution volume	approx. 270 - 320ml
Quantity	1

### Safety Lid

Overall dimension	197mm(W) × 38mm(H) × 169mm(L)
Material characteristic	UV non-transparent
Quantity	1

### Intelligent Power Supply

Overall dimension	75mm(W) × 62mm(H) × 170mm(L)
Weight	410 g
Input Voltage	AC100 - 240V, 50/60Hz
Output Voltage	135V, 100V, 50V, 25V, 70V, 35V, 18V Constant peak voltage of 140V, and duty-control
Timer	Timer operation 0 - 99 min, and continuous mode Temporary shutdown supported
Safety Switch	Micro-switches in the Power Supply Without safety lid, no output of Power
Memory Function	Automatic memory (the last used V and T)
Quantity	1

### Gel Tray

Gel tray(small)	130mm(W) × 13mm(H) × 59.5mm(L)
Gel tray(large)	130mm(W) × 13mm(H) × 122mm(L)
Quantity	small: 2, large: 1

### Comb

Multiple number of wells	13 - 26 Wells 13 wells: 9mm spacing 26 wells: 4.5mm spacing
Quantity	4 (supports making 13 or 26 wells)

### Gel Casting Stand


Overall dimension	149mm(W) × 20mm(H) × 125mm(L)
Quantity	1

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### 1. Mupid-EX<sup>U</sup> Construction and Functions

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
### 1-2 Gel Casting Set

Made of UV transparent plastic, all Mupid-EX casting trays allow visualization of samples under UV light without removing the tray. This saves time and minimizes user exposure during handling of gels that may contain hazardous chemicals such as ethidium bromide.

#### Gel Tray (Small)

Small gels are made using this tray. Convenient transmittance of UV illumination for observation using a transilluminator. Maximum separation length is approx. 5.5 cm with this tray. The width of this tray suits the usage of a multi-channel pipette.

#### Gel Tray (Large)

Large gels are made with this tray. Convenient transmittance of UV illumination for observation using a transilluminator. Maximum separation length is ca. 11cm with this tray. The width of this tray suits the usage of a multi-pipette.


#### Combs

Set these combs in the gel maker stand to form wells. Double-sided combs are provided: 13 teeth (5 mm wide wells), and 26 teeth (2 mm wide wells). The spacing of these wells is compatible with a multi-channel pipette.

#### shape

The 26-toothed comb is designed for excellent band resolution.

Direction of electrophoresis (in the case of DNA)  
With the 26-toothed comb, electrophoresis is the direction as shown in this figure.

#### Gel Casting Stand

Set gel trays in this stand to cast gels.

#### Center Partition

When using the large gel tray, remove this center partition.

#### Comb Setting Sills

Up to four combs can be set in the stand for large gels, use the small round marks as the setting position (See page 6).

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