



Total-Glutathione (GSH) Assay Kit

(Cat/No.:BC031 Size:96T/48T Microplate method)

1. Principle of Measurement

The total-GSH or oxidized GSH (Glutathione disulfide, GSSG for short) content of tissues or body fluid is determined through the reversible reaction with 5,5'-dithiobis(2-nitrobenzoic acid) (DTNB).

2. Composition and Preparation(The kit is valid for 6 months)

	Composition	48T	96T	Storage
Reagent I	Substrate I Powder	1 Bottles	2 Bottles	4 °C
	Buffer Solution	1 Bottle×5ml	2 Bottles×5ml	4 °C
Preparation of Reagent I Solution : Dilute each one bottle powder with 5ml buffer solution and dissolve thoroughly. The solution should be preserved at 4 °C				
Reagent II	Substrate II Stock	25μl×1 tube	25μl×2 tubes	-20 °C
	Diluent	500μl×tube	500μl×2 tubes	-20 °C
Preparation of Reagent II Solution : Dilute the stock with diluent at the ratio 1:19. Prepare the solution prior to use with the amount needed.				
Reagent III	Powder	Powder×3 tubes	Powder×5 tubes	-20 °C
	Dilution	1ml×3 tubes	1ml×5 tubes	4 °C
Preparation of Reagent III Solution : Dissolve the powder with 1ml diluent. Prepare the solution prior to use with the amount needed.				
Reagent IV	Powder	Powder×3 bottles	Powder×5 bottles	4 °C
Preparation of Reagent IV : Dissolve each bottle of powder with boiled double distilled water (DDW) to 10ml and cool the solution. The solution can be preserved at 4 °C for 3 days.				
Reagent V	Stock Solution	1 Bottle×15μl	1 Bottle×30μl	4 °C
	Reagent V Solvent	1 Bottle×150μl	1 Bottle×300μl	4 °C
Preparation of Reagent V : Prior to use, dilute the stock solution with solvent at the ratio 1:9.				
Reagent VI	Solution	1 Bottle×250μl	1 Bottle×500μl	4 °C
Note : Reagent VI is viscid and be patient and careful when extracting the solution.				
GSSG Standard	6.13mg	Powder×1 tube	Powder×1 tube	4 °C
Preparation of 1mM GSSG Stock Solution : Dissolve the powder with 10ml DDW and the solution is stable at -20 °C for a month.				
Preparation of 50μM GSSG Standard Solution: Dilute the 1mM solution to 20 times the initial volume. Prepare				

before use.

GSH Standard	3.07mg	Powder×1 tube	Powder×1 tube	4 °C
Preparation of 1mM GSH Stock Solution : Dissolve the powder with 10ml DDW and the solution is stable at -20 °C for a month.				



Preparation of **50 μ M GSH Standard Solution**: Dilute the 1mM solution to 20 times the initial volume. Prepare before use.

3. Required Instruments and Reagents

Adjustable 405nm wavelength microplate reader and 96-well plate (one plate included), 37°C water bath or incubator, vortex mixer, stopwatch, double-distilled water (or distilled water).

4. Sample Pretreatment

I. Whole Blood Sample

- a. Take blood and add Heparin or EDTA for anticoagulation.
- b. Extract 100 μ l blood sample and add 400 μ l freshly prepared reagent IV solution with the dilute coefficient of 5. Vortex for 30 s and then set aside the sample at 4°C for 5 min.
- c. Centrifuge at 3,500 rpm for 10 min and extract the supernatant at 4°C. The supernatant should be placed under -20°C for preservation.

II. Red Blood Cell Sample

- a. Take blood and add Heparin or EDTA for anticoagulation.
- b. Centrifuge at 2,000 rpm for 10 min and carefully remove the upper plasma and the white blood cell layer on the surface of red blood cell layer.
- c. Extract 100 μ l red blood cell and add 400 μ l freshly prepared reagent IV solution with the dilute coefficient of 5. Vortex for 30 s and then set aside the sample at 4°C for 5 min.
- d. Centrifuge at 3,500 rpm for 10 min and extract the supernatant at 4°C. The supernatant should be placed under -20°C for preservation.

III. Serum or Plasma Sample

- a. Take blood and add Heparin or EDTA for anticoagulation.
- b. Centrifuge at 2,000 rpm for 10 min and carefully extract the upper plasma layer.
- c. Extract 100 μ l plasma and add 400 μ l freshly prepared reagent IV solution with the dilute coefficient of 5. Vortex for 30 s and then set aside the sample at 4°C for 5 min.
- d. Centrifuge at 3,500 rpm for 10 min and extract the supernatant at 4°C. The supernatant should be placed under -20°C for preservation.

IV. Tissue Sample

- a. Rinse the freshly prepared tissues with saline and remove the excess water on tissues.
- b. Weigh the tissues and add reagent IV with the ratio of 1:4(g/ml) and homogenize in an ice water bath.
- c. Centrifuge the homogenate at 3,500rpm for 10 min and extract the supernatant at 4°C. The supernatant should be placed under -20°C for preservation.

5. Procedures of Measurement

I. T-GSH Measurement

Compositions	Standard	Blank
50 μ M GSH Standard(μ l)	10	
Sample(μ l)		10
Reagent I(μ l)	100	100
Reagent II(μ l)	10	10
Mix and after mixing, set aside at room temperature for 2 min		
Reagent III(μ l)	50	50

Record the timing right after the addition of reagent III. Mix and extract the solution into microplate. At 405nm, record the absorbance at 30 s (A_1) and set aside for 5 min at room temperature. At 5 min 30 s, record the absorbance (A_2) at the same wavelength.

II. GSSG Measurement

a. Pre-Treatment

Compositions	Standard	Sample
50 μ M GSSG Standard(μ l)	100	
Sample(μ l)		100
Reagent V(μ l)	2	2
Reagent VI(μ l)	5	5

Vortex for 1min and warm the solution at 37°C for 30 min for further use.

b. GSSG Measurement

Compositions	Standard	Sample
Standard Pretreatment(μ l)	10	
Sample Pretreatment(μ l)		10
Reagent I(μ l)	100	100
Reagent II(μ l)	10	10
Mix and after mixing, set aside at room temperature for 2 min		
Reagent III(μ l)	50	50

Record the timing right after the addition of reagent III. Mix and extract the solution into microplate. At 405nm, record the absorbance at 30 s (A_1) and set aside for 5 min at room temperature. At 5 min 30 s, record the absorbance (A_2) at the same wavelength.

Note: Pre-measurement can be done for 1-2 samples and in case high absorbance results obtained, dilute with reagent IV so that the T-GSH or GSSG concentration lays within the concentration range required for measurement.

6. Standard Curve Establishment

I. GSH Standard Curve

- Dissolve GSH in DDW to the concentration of 1mM as the standard solution and standard solution can be preserved at -20°C for a month.
- Dilute the standard solution to 0mM (Blank), 0.0125mM, 0.025mM, 0.05mM and 0.1mM at 4°C for further use.

c. Procedures


Compositions (μl)	GSH Test Tube
GSH Solution with Different Concentrations(μl)	10
Reagent I(μl)	100
Reagent II(μl)	10
Mix and set aside the solution at room temperature for 2 min	
Reagent III(μl)	50

Record the timing right after the addition of reagent III. Mix and extract the solution into microplate. At 405nm, record the absorbance at 30 s (A_1) and set aside for 5 min at room temperature. At 5 min 30 s, record the absorbance (A_2) at the same wavelength.

d. Results

GSH Concentration	A_1	A_2	ΔA
0μmol/L	0.1686	0.1686	0.0000
12.5μmol/L	0.1851	0.2173	0.0322
25μmol/L	0.2005	0.2588	0.0583
50μmol/L	0.2306	0.3374	0.1068
100μmol/L	0.2722	0.4786	0.2064

e. Standard Curve



7. GSSG Standard Curve

- Dissolve GSSG in DDW to the concentration of 1mM as the standard solution and standard solution can be preserved at -20°C for a month.
- Dilute the standard solution to 0mM (Blank), 0.0125mM, 0.025mM, 0.05mM and 0.1mM at 4°C for further use.
- Procedures

i. Pre-treatment

Compositions	GSSG Test Tube
GSSG Solution with Different Concentrations (μl)	100
Reagent V (μl)	2
Reagent VI (μl)	5

Vortex for 1min and warm the solution at 37°C for 30 min for further use.

ii. GSSG Measurement


Compositions	GSSG Test Tube
Pre-treated Solution with Different Concentrations (μl)	10
Reagent I (μl)	100
Reagent II (μl)	10
Mix and then set aside at room temperature for 2min	
Reagent III	50

Record the timing right after the addition of reagent III. Mix and extract the solution into microplate. At 405nm, record the absorbance at 30 s (A_1) and set aside for 5 min at room temperature. At 5 min 30 s, record the absorbance (A_2) at the same wavelength.

- Results

GSH Concentration	A ₁	A ₂	ΔA
0μmol/L	0.1640	0.1640	0.0000
12.5μmol/L	0.1720	0.2302	0.0582
25μmol/L	0.1744	0.2799	0.1055
50μmol/L	0.1799	0.3781	0.1982
100μmol/L	0.1906	0.5679	0.3773

e. GSSG Standard Curve



8. Calculation Formula

I. Formula

$$\text{Total — GSH Content } (\mu\text{mol/L}) = \frac{\Delta A_{\text{Sample}}(A_2 - A_1)}{\Delta A_{\text{Standard}}(A_2 - A_1)} \times \frac{C_{\text{Standard}}}{50\mu\text{M}} \times \text{Coefficient of Dilution (pretreatment Caused)}$$

$$\text{GSSG Content } (\mu\text{mol/L}) = \frac{\Delta A_{\text{Sample}}(A_2 - A_1)}{\Delta A_{\text{Standard}}(A_2 - A_1)} \times \frac{C_{\text{Standard}}}{50\mu\text{M}} \times \text{Coefficient of Dilution (Pretreatment Caused)}$$

$$\text{Reduced GSH Content} = \text{Total — GSH Content} - 2 \times \text{GSSG Content}$$

II. Standard Curve

- Find the corresponding T-GSH content with the standard curve established and multiply by the coefficient of dilution according to the steps of pre-treatment in order to achieve the T-GSH of the sample.
- Find the corresponding GSSG content with the standard curve established and multiply by the coefficient of dilution according to the steps of pre-treatment in order to achieve the GSSG of the sample.
- Reduced GSH Content = Total — GSH Content — 2 × GSSG Content

9. Note

- I. For the fast rate of GSH metabolism, please treat the sample with reagent IV as soon as possible in order to lower the loss of GSH in sample.
- II. The supernatant which is obtained via the treatment of samples with reagent IV can be preserved at -20°C for 6 months
- III. The 5 minutes reaction time should be recorded precisely.

10. Significance of Measurement

Researchers show increasingly interest in Glutathione for its widespread existence and multifunction in vivo. It directly or indirectly functions in various biological phenomena and thus involved in numerous research topics such as mechanism of enzyme catalysis, synthesis of macromolecules like protein and DNA, metabolism, radiation, neoplasm, immunization, environmental toxins and ageing. Quite a lot of topics involving the GSH roles in cell defense against active oxygen or radicals. Among the samples measured, the GSH and GSSG content change in serum sample is a key indicator of oxidative stress implicated in vivo.