

## Certificate of Analysis

### ImProm-II™ Reverse Transcriptase

Cat. #	Size
A3801	10 reactions
A3802	100 reactions
A3803	500 reactions

**Description:** ImProm-II™ Reverse Transcriptase<sup>(a)</sup> is provided in quantities sufficient for 10–500 first-strand cDNA synthesis reactions of 20µl each. **Additional components are required to perform the reactions.** Contains one of the following:

#### Cat. # A3801

Part No.	Component	Size
M314B	ImProm-II™ Reverse Transcriptase, 1µl/reaction	10µl
M289B	ImProm-II™ 5X Reaction Buffer	100µl
A351B	MgCl <sub>2</sub> , 25mM	750µl

#### Cat. # A3802

Part No.	Component	Size
M314A	ImProm-II™ Reverse Transcriptase, 1µl/reaction	100µl
M289A	ImProm-II™ 5X Reaction Buffer	600µl
A351H	MgCl <sub>2</sub> , 25mM	1.2ml

#### Cat. # A3803

Part No.	Component	Size
M314C	ImProm-II™ Reverse Transcriptase, 1µl/reaction	500µl
M289C	ImProm-II™ 5X Reaction Buffer	2 × 1,250µl
A351H	MgCl <sub>2</sub> , 25mM	3 × 1.2ml

**Storage Conditions:** See the Product Information Label for storage recommendations. Avoid multiple freeze-thaw cycles and exposure to frequent temperature changes. These fluctuations can greatly alter product stability.

**Usage Note:** Additional components are required to perform RT-PCR. These include: appropriate primers, nuclease-free water, dNTPs and *Taq* DNA polymerase. Recombinant RNasin® Ribonuclease Inhibitor is optional but recommended.

Part# 9PIA380

Revised 10/16



AF9PIA380 1016A380



**Promega**

#### Promega Corporation

2800 Woods Hollow Road	
Madison, WI 53711-5399	USA
Telephone	608-274-4330
Toll Free	800-356-9526
Fax	608-277-2516
Internet	<a href="http://www.promega.com">www.promega.com</a>

## Quality Control Assays

### Activity Assays

**First-Strand cDNA Synthesis:** Two hundred units of enzyme are used to produce cDNA from 1µg of 1.2kb and 6.5kb control RNAs in separate reactions, using [<sup>32</sup>P] dCTP as a tracer. The minimum specification is 120ng of first-strand cDNA made from 1µg of RNA. The cDNA product must be >90% full length as determined by gel electrophoresis and autoradiography.

**Amplification:** When 0.25 zeptomoles (approximately 100 copies) of 1.2kb Kanamycin Positive Control RNA (Cat.# C1381) is reverse transcribed at 42°C and amplified, the result is a clear, discrete 323bp DNA product as visualized on an agarose gel by ethidium bromide staining.

#### PRODUCT USE LIMITATIONS, WARRANTY DISCLAIMER

Promega manufactures products for a number of intended uses. Please refer to the product label for the intended use statements for specific products. Promega products contain chemicals which may be harmful if misused. Due care should be exercised with all Promega products to prevent direct human contact.

Each Promega product is shipped with documentation stating specifications and other technical information. Promega products are warranted to meet or exceed the stated specifications. Promega's sole obligation and the customer's sole remedy is limited to replacement of products free of charge in the event products fail to perform as warranted. Promega makes no other warranty of any kind whatsoever, and SPECIFICALLY DISCLAIMS AND EXCLUDES ALL OTHER WARRANTIES OF ANY KIND OR NATURE WHATSOEVER, DIRECTLY OR INDIRECTLY, EXPRESS OR IMPLIED, INCLUDING, WITHOUT LIMITATION, AS TO THE SUITABILITY, PRODUCTIVITY, DURABILITY, FITNESS FOR A PARTICULAR PURPOSE OR USE, MERCHANTABILITY, CONDITION, OR ANY OTHER MATTER WITH RESPECT TO PROMEGA PRODUCTS. In no event shall Promega be liable for claims for any other damages, whether direct, incidental, foreseeable, consequential, or special (including but not limited to loss of use, revenue or profit), whether based upon warranty, contract, tort (including negligence) or strict liability arising in connection with the sale or the failure of Promega products to perform in accordance with the stated specifications.

© 2001–2016 Promega Corporation. All Rights Reserved.

RNasin is a registered trademark of Promega Corporation. ImProm-II is a trademark of Promega Corporation.

Products may be covered by pending or issued patents or may have certain limitations. Please visit our Web site for more information.

All prices and specifications are subject to change without prior notice.

Product claims are subject to change. Please contact Promega Technical Services or access the Promega online catalog for the most up-to-date information on Promega products.

Part# 9PIA380  
Printed in USA. Revised 10/16.

Signed by:

R. Wheeler, Quality Assurance

## I. Standard Reverse Transcription Protocol

### A. Prepare RNA Target and Primer

1. Use sterile, nuclease-free, thin-walled tubes, prechilled on ice.
2. For each 20µl reverse transcription (RT) reaction, combine:

Component	Volume
RNA template	up to 1µg
primer	<u>20pmol or 0.5µg</u>
nuclease-free water to a final volume of	5µl

3. Incubate at 70°C for 5 minutes.
4. Quick-chill at 4°C for 5 minutes and hold on ice.

#### Notes:

- RT reactions of 20µl process up to 1µg of total RNA, mRNA or synthetic transcript, and allow detection of 100 copies of a specific RNA target.
- Primer amounts recommended for efficient cDNA synthesis are 0.5µg of oligo(dT) or random primers, or 20pmol of gene-specific primers per 20µl reaction.
- For best results, heat denature the primers and target in a controlled-temperature heat block; then chill immediately in an ice water bath for at least 5 minutes. Hold on ice until addition to RT reaction.

### B. Prepare Reverse Transcription Mix

1. For each 20µl RT reaction, combine

Component	Final Concentration	Volume
nuclease-free water		Xµl
ImProm-II™ 5X Reaction Buffer	1X	4µl
MgCl <sub>2</sub> , 25mM*	3mM*	2.4µl
dNTP mix* (10mM each dNTP)	0.5mM*	1µl
Recombinant RNasin® Ribonuclease Inhibitor	1u/µl	20u
Vortex the mixture. Add:		
ImProm-II™ Reverse Transcriptase		<u>1µl</u>
Final Volume RT Mix per 20µl reaction		15µl

\*See notes below

2. Vortex gently to mix.
3. Dispense 15µl aliquots into reaction tubes.

#### Notes:

- Combine indicated reagents in the order specified in a sterile, nuclease-free reaction tube.
- MgCl<sub>2</sub> may be optimized between 1.5–8.0mM.
- Final dNTP concentration is 0.5mM.
- Use of Recombinant RNasin® Ribonuclease Inhibitor is optional but recommended.
- We recommend keeping the RT reaction mix chilled on ice prior to incubation.

### C. Add Template + Primers to the Reaction Mix

1. For each individual reaction, add 5µl of the appropriate Template + Primer mix to the 15µl reverse transcription or RT-PCR mix. If necessary, overlay with nuclease-free mineral oil. The final volume for each individual reaction will be 20µl.

#### Notes:

- Add the RNA Template + Primer mix immediately prior to incubation.
- Keep tubes closed whenever possible to avoid contamination, and be careful to avoid cross contamination.

## D. Reverse Transcription

1. Anneal at 25°C for 5 minutes.
2. Extend the first strand for 60 minutes at 42°C. The extension temperature may be optimized between 37–55°C.
3. Heat-inactivate the ImProm-II™ Reverse Transcriptase by incubating at 70°C for 15 minutes.
4. Analyze cDNA, proceed with PCR or store frozen.

#### Notes:

- Place the reaction tubes in a controlled-temperature heating block.
- Annealing conditions may require optimization. See Technical Manual #TM236 for details regarding optimization.
- First-strand cDNA synthesis occurs during the extension step. The extension temperature may require optimization. See Technical Manual #TM236 for details regarding optimization.

## E. PCR Amplification (Optional)

1. The cDNA may be amplified by adding the products of the heat inactivated reverse transcription reaction directly to the PCR mix and proceeding with thermal cycling.
2. Prepare PCR Mix, minus the cDNA sample by combining the amplification reagents in a sterile 1.5ml microcentrifuge tube on ice. Combine the components in the order listed, vortex gently to mix, and keep on ice.

Component	Final Concentration	Volume per 100µl reaction (20µl RT reaction)
nuclease-free water		Xµl
10X thermophilic polymerase reaction buffer (without MgCl <sub>2</sub> )*	1X	8.0µl
MgCl <sub>2</sub> , 25mM	2mM*	5.6µl
dNTP mix, 10mM	0.2mM	1µl
upstream primer	1µM	Xµl
downstream primer	1µM	Xµl
Taq DNA polymerase	<u>5.0u</u>	<u>1.0µl</u>
PCR Mix volume per reaction		80µl
Volume of RT reaction added		<u>20.0µl</u>
Total PCR Volume		100.0µl

\*For experimental systems, Mg<sup>2+</sup> should be optimized between 1.5–2.5mM.

3. Dispense the appropriate volume of PCR mix into each chilled reaction tube.
4. Add the appropriate aliquot of the reverse transcription reaction to the PCR mix.
5. If necessary, overlay the reaction with nuclease-free mineral oil and proceed with amplification program.

#### Notes:

- In the table of PCR component additions shown, the amounts of 10X Reaction Buffer, 25mM MgCl<sub>2</sub> and dNTP mix combined in the PCR amplifications supplement the specified volume of reverse transcription reaction that is added. If the MgCl<sub>2</sub> concentration of the RT reaction is changed, these proportions will require modifications.
- **The volume of the PCR amplification may be scaled for volumes less than 100µl. Carryover concentrations of MgCl<sub>2</sub>, dNTPs, buffer and primers must be considered when calculating the volumes of the PCR mix components.**