

TECHNICAL MANUAL

ADCC Reporter Bioassay, Core Kit

Instructions for Use of Products
G7010 and G7018



ADCC Reporter Bioassay, Core Kit

All technical literature is available at: www.promega.com/protocols/
 Visit the web site to verify that you are using the most current version of this Technical Manual.
 E-mail Promega Technical Services if you have questions on use of this system: techserv@promega.com

1.	Description.....	1
2.	Product Components and Storage Conditions.....	6
3.	General Considerations	7
4.	Assay Protocol for ADCC Reporter Bioassay, Core Kit	7
4.A.	Before You Begin	7
4.B.	Preparation of Components, Reagents and Bioassay Starting Materials.....	9
4.C.	Plate Layout Design	10
4.D.	Preparing Antibody Serial Dilutions.....	10
4.E.	Plating Target Cells.....	12
4.F.	Adding Antibody to Target Cells in Assay Plates	13
4.G.	Plating ADCC Bioassay Effector Cells	13
4.H.	Adding Bio-Glo™ Luciferase Assay Reagent	14
4.I.	Data Analysis	14
5.	Troubleshooting.....	15
6.	References.....	17
7.	Appendix: Representative Assay Results	18
8.	Related Products.....	21
9.	Summary of Changes	23

1. Description

Antibody-dependent cell-mediated cytotoxicity (ADCC) is a mechanism of action of antibodies through which virus-infected or other diseased cells are targeted for destruction by components of the cell-mediated immune system, such as natural killer cells. The ADCC Reporter Bioassay^(a-d) is a bioluminescent reporter assay for quantifying biological activity on pathway activation by therapeutic antibody drugs in an ADCC mechanism of action (MOA) assay. The assay combines a simple, add-mix-read format, effector cells provided in a frozen, thaw-and-use format, and an optimized protocol for a bioassay that has low variability and high accuracy. Moreover the bioassay can be performed in a single day. These performance characteristics make the bioassay suitable for applications across antibody drug research, development and manufactured lot release. The thaw-and-use cells provided in the ADCC Reporter Bioassay kits are generated under highly controlled conditions that drive low assay variability run to run, while providing the convenience of an assay reagent that eliminates the need to propagate and prepare cells each time.

In addition to the ADCC Reporter Bioassay, Core Kit (Cat.# G7010, G7018), we offer two versions of a Complete Kit: ADCC Reporter Bioassay, Complete Kit (WIL2-S), Cat.# G7014, and ADCC Reporter Bioassay, Complete Kit (Raji), Cat.# G7015, in which we provide the target cells and Control Ab, Anti-CD20. We also offer the effector cells for banking and propagation, the ADCC Reporter Bioassay, Cell Propagation Model (Cat.# G7102), under a unique purchase agreement.

ADCC is a desirable mechanism for killing target cancer cells using antibody-based drugs. The antibody binds to target antigens on the cell surface. When the Fc effector portion of target-bound antibodies also binds to FcγRIIIa receptors on the cell surface of effector cells (natural killer cells predominantly), multiple cross-linking of the two cell types occurs, leading to pathway activation of ADCC MOA (1). Killing of target cells is an endpoint of this pathway activation and is used in classic ADCC bioassays, which use donor peripheral blood mononuclear cells (PBMCs) or the natural killer (NK) cell subpopulation as effector cells (2). These cells can be highly variable in response, are tedious to prepare and can result in high background readings.

The ADCC Reporter Bioassay uses an alternative readout at an earlier point in ADCC MOA pathway activation: the activation of gene transcription through the NFAT (nuclear factor of activated T-cells) pathway in the effector cell (3,4). In addition, the ADCC Reporter Bioassay uses engineered Jurkat cells stably expressing the FcγRIIIa receptor, V158 (high affinity) variant, and an NFAT response element driving expression of firefly luciferase as effector cells. Antibody biological activity in ADCC MOA is quantified through the luciferase produced as a result of NFAT pathway activation; luciferase activity in the effector cell is quantified with luminescence readout (Figure 1). Signal is high, and assay background is low.




Figure 1. Representation of the ADCC Reporter Bioassay. Readout is luminescence signal from NFAT response element driving expression of firefly luciferase.

The ADCC Reporter Bioassay exhibits the clear specificity desired for a bioassay, as shown in Figure 2. A good assay response is only obtained when target cells with the correct surface antigen, the correct specific antibody, and effector cells expressing FcγRIIIa are present. If any one of these is missing, there is no response.




Figure 2. Specificity of the ADCC Reporter Bioassay. Serial dilutions of rituximab (anti-CD20 chimeric monoclonal antibody drug), trastuzumab (anti-Her2 humanized monoclonal antibody drug), or assay medium control (no antibody) were incubated for 6 hours of induction at 37°C with engineered Jurkat effector cells (ADCC Bioassay Effector Cells), with or without ADCC Bioassay Target Cells (WIL2-S), as indicated. Luciferase activity was quantified using Bio-Glo™ Reagent^(e). Data were fitted to a 4PL curve using GraphPad Prism® software.

The ADCC Reporter Bioassay has performance characteristics suitable for many applications of a bioassay used across antibody drug discovery, development and manufacture: It is stability-indicating and has the precision and accuracy suitable for a lot-release bioassay (Figure 3). Additionally the assay can be used to quantify effects of glycosylation differences on Fc effector function of antibodies in ADCC MOA (Figure 4), which would be useful for ADCC efficiency variant analysis, for example (5). Benchmarking studies demonstrate the ADCC Reporter Bioassay provides antibody activity ranking equivalent to classic ADCC bioassay using PBMCs and LDH release as a measure of target cell death (Figure 5).

Parameter		Results	
		WIL2-S Target Cells	Raji Target Cells
Accuracy	% Expected Relative Potency	% Recovery	% Recovery
	50	97.7	101.0
	75	88.5	101.2
	125	98.4	93.8
	150	98.4	96.7
Repeatability (%CV)	(100% Reference)	5.0	9.1
Intermediate Precision (%CV)		7.3	3.0
Linearity (r^2)		0.995	0.997
Linearity ($y = mx + b$)		$y = 1.016x - 0.052$	$y = 0.922x + 5$

Figure 3. Bioassay characterization. The ADCC Reporter Bioassay was characterized in studies that evaluated accuracy, repeatability, intermediate precision and linearity across the 50–150% relative potency range. Dilution ranges for an anti-CD20 antibody were selected to ensure good coverage of upper and lower asymptotes and sufficient points in the intermediate dose-range for accurate slope and EC_{50} determinations. A series of relative potency samples, of 50%, 75%, 125% and 150% theoretical relative potency, were evaluated as triplicate dilution series of antibody dose on each of 3 different days. The effector-to-target cell ratio (E:T ratio) was 6:1. The ADCC Reporter Bioassay was characterized using ADCC Bioassay Target Cells (WIL2-S) and ADCC Bioassay Target Cells (Raji). Data were fitted to a 4PL curve using GraphPad Prism® software, and relative potencies were calculated after parallelism determination using SAS Institutes, Inc. JMP® software. Relative potencies were calculated using the 100% reference sample run as a triplicate dilution series in the same assay plate as the test sample.




Figure 4. Detection of antibody glycosylation. Rituximab-blended samples containing mixtures of fully deglycosylated and fully N-glycosylated antibody were assayed in the ADCC Reporter Bioassay against a 100% relative activity reference sample of fully N-glycosylated rituximab in the same plate. Target cells were ADCC Bioassay Target Cells (WIL2-S), and the E:T ratio was 6:1. Biological activity was expressed relative to the 100% control run in the same assay plate and plotted against the % of N-glycosylation present. Linear regression analysis was performed to determine correlation.




Figure 5. Antibody biological activity correlated with antibody glycosylation in an ADCC reporter assay and a classic PBMC-based ADCC assay. Blended antibody samples containing 10%, 20%, 30%, 40% or 50% of untreated trastuzumab (fully glycosylated) were prepared by mixing untreated and PNGase F-treated trastuzumab in appropriate proportions. Each of the blended trastuzumab samples was then assayed in triplicate against an untreated trastuzumab (reference sample) in the same plate using an ADCC reporter assay or a PBMC-based ADCC assay. Freshly cultured HER2⁺ SK-BR-3 cells were used as target cells in both assays. For the ADCC reporter assay, the E:T ratio was 15:1 and induction time was 6 hours. Results are mean \pm SD of three independent experiments using the same batch of frozen, thaw-and-use Jurkat effector cells. For the PBMC-based ADCC assay, the effector cells were PBMCs from the donors with Fc γ R3IIa V/V genotype for aa 158 and the E:T ratio was 50:1. After an overnight incubation, target cell lysis was measured by detecting the release of lactate dehydrogenase from lysed SK-BR-3 cells using CytoTox 96[®] Non-Radioactive Cytotoxicity Assay. Results are mean \pm SD of three independent experiments using PBMCs isolated from different donors. Percent ADCC biological activity is defined as the ratio of the EC₅₀ of untreated trastuzumab (N-glycosylated) to the EC₅₀ of the blended trastuzumab mixture in the same assay plate.

The ADCC Reporter Bioassay Core Kit is suitable for use with your biologic Ab and target cell line of interest. It offers the greatest flexibility for different target cell types that you provide from continuous cell culture. Suspension or adherent target cells can be used with the engineered Jurkat effector cells (Figures 9–11 in the Appendix, Section 7). To obtain the optimal results in routine application of the assay using your target system, we recommend optimizing the E:T ratio, dose-range of antibody to generate a full dose-response, and possibly response induction time. We recommend that any first-time user read the protocol for this kit in its entirety before using it. Since the assay is different from a classic ADCC bioassay, we recommend a first-time user also purchase an ADCC Reporter Bioassay, Complete Kit (WIL2-S) or ADCC Reporter Bioassay, Complete Kit (Raji), as a starter kit and use it in initial experiments to gain familiarity with format and protocol.

2. Product Components and Storage Conditions

Note: The ADCC Reporter Bioassay components are shipped separately because of temperature requirements. The ADCC Bioassay Effector Cells are shipped on dry ice. The Bio-Glo™ Luciferase Assay System⁽⁶⁾ and Low IgG Serum are shipped on dry ice, separately from the cells. The RPMI 1640 Medium is shipped at ambient temperature.

PRODUCT	SIZE	CAT.#
ADCC Reporter Bioassay, Core Kit	1 kit	G7010

Each system contains sufficient reagents for 120 assays using the inner 60 wells of two 96-well plates. Includes:

- 1 vial ADCC Bioassay Effector Cells (0.65ml)
- 4ml Low IgG Serum
- 10ml Bio-Glo™ Luciferase Assay Buffer
- 1 vial Bio-Glo™ Luciferase Assay Substrate (lyophilized)
- 36ml RPMI 1640 Medium

PRODUCT	SIZE	CAT.#
ADCC Reporter Bioassay, Core Kit 5X	5 kits	G7018

Each system contains sufficient reagents for 600 assays using the inner 60 wells of ten 96-well plates. Includes:

- 5 vials ADCC Bioassay Effector Cells (0.65ml)
- 5 × 4ml Low IgG Serum
- 5 × 10ml Bio-Glo™ Luciferase Assay Buffer
- 5 vials Bio-Glo™ Luciferase Assay Substrate (lyophilized)
- 5 × 36ml RPMI 1640 Medium

Storage Conditions: Upon arrival, immediately transfer the vials of ADCC Bioassay Effector Cells for long-term storage below –140°C (freezer or liquid nitrogen vapor phase). The cells are sensitive, and care should be taken when handling. For safety reasons do not store cell vials submerged in liquid nitrogen. Low IgG Serum should be stored at –20°C. Avoid multiple freeze-thaw cycles of the serum. Bio-Glo™ Luciferase Assay Buffer and Bio-Glo™ Luciferase Assay Substrate should be stored at –20°C. For optimal performance, reconstituted Bio-Glo™ Luciferase Reagent should be used on the day of preparation. However, once reconstituted, Bio-Glo™ Luciferase Assay Reagent can be stored –20°C for up to 6 weeks. RPMI 1640 Medium should be stored at 4°C protected from fluorescent light.

3. General Considerations

The ADCC Reporter Bioassay differs from classic ADCC bioassays in a number of ways. **Please read through the entire protocol for this kit to become familiar with the assay, the components and the protocol in general before beginning.** The ADCC Bioassay Effector Cells, when thawed and diluted as instructed, will be at the proper concentration for the bioassay. The effector:target (E:T) cell ratio, the antibody dose range, assay buffer and incubation times may differ from those used in a classic ADCC bioassay with PBMCs or natural killer cells as effector cells. **We recommend that you evaluate these parameters with your target cells and select the best conditions for your target system.**

The ADCC Bioassay Effector Cells are provided in a frozen, thaw-and-use format and are ready to be used without any culturing procedures. Although they are ready-to-use, the cells are sensitive, and care should be taken not to over-mix the cell reagents, as is commonly done with other cell preparations that tend to clump or settle. Follow the protocol instructions, warming the cells until just thawed (Section 4.G).

Because the ADCC Reporter Bioassay produces a bioluminescent readout, the assay requires a sensitive luminometer or luminescence plate reader for the detection of signal. See Related Products, Section 8, for a list of GloMax[®] Detection Systems available from Promega. The bioassay produces a strong signal; therefore, an integration time of 0.5sec/well should be sufficient. If your luminometer/plate reader requires gain adjustment for luminescence, use the well with the highest Ab concentration. Finally, if you have the ability to select the multiwell plate type in your reader's software and that multiwell plate is not listed in the software, a generic 96-well plate selection will suffice. We recommend white, flat-bottom 96-well assay plates (Corning Cat.# 3917).

4. Assay Protocol for ADCC Reporter Bioassay, Core Kit

4.A. Before You Begin

Materials to Be Supplied by the User

- user-defined target cells
- sterile clear 96-well, V-bottom plate with lid (Linbro Cat.# 76-223-05 or equivalent) for preparing antibody dilutions
- white, flat-bottom 96-well assay plates (Corning Cat.# 3917 or equivalent)
- pipettes (single channel and 12-channel)
- sterile 15ml and 50ml conical tubes
- sterile reagent reservoirs (Corning Cat.# 4870 or equivalent)
- 37°C CO₂ incubator
- 37°C water bath
- plate reader with glow luminescence read capability or luminometer (e.g., GloMax[®] Discover System)
- reference antibody
- test antibody

Overview

One recommended protocol, described here as an example protocol, is designed to test two antibody samples in a single assay run. Each test antibody and a reference antibody are run as triplicate 10-point dilution series in a single 96-well assay plate using the same target cells, for a total of two plates. Other protocols and plate layouts are possible and may need to be optimized for your specific target antibody and target cells.




Figure 6. Schematic protocol for the ADCC Reporter Bioassay.

4.B. Preparation of Components, Reagents and Bioassay Starting Materials

1. **Bio-Glo™ Luciferase Assay Reagent:** Prepare Bio-Glo™ Luciferase Assay Reagent 4–6 hours before use on the day of assay. First, thaw one bottle of Bio-Glo™ Luciferase Assay Buffer. Equilibrate Bio-Glo™ Luciferase Assay Buffer and Bio-Glo™ Luciferase Assay Substrate to ambient temperature. Transfer the buffer into the amber bottle containing Substrate. Mix by inversion until the Substrate is thoroughly dissolved. Store reconstituted Bio-Glo™ Luciferase Assay Reagent at ambient temperature (22–25°C), protected from light. Approximate stability of Bio-Glo™ Reagent after reconstitution: 18% loss of luminescence over 24 hours at ambient temperature.
2. **ADCC Assay Buffer:** On the day of assay, thaw the Low IgG Serum in a 37°C water bath. In a 50ml conical tube, add 1.4ml Low IgG Serum to 33.6ml of RPMI 1640 Medium to make 35ml of ADCC Assay Buffer for two assay plates. Mix well and warm to 37°C prior to use.
3. **Starting dilutions (Dilu1, 3X final concentration) for Reference antibody and two Test antibodies:** Decide the starting concentration (1X) for reference antibody and two test antibody samples based on previous testing results in conventional ADCC cytotoxicity assay if available. If the working concentration of test antibody is unknown, use 1µg/ml as starting concentration, and adjust later based on the assay results; this starting concentration has worked for both rituximab and trastuzumab in the ADCC Reporter Bioassay.

Make 400µl of starting dilution for reference antibody (dilu1, 3X final concentration) and make 200µl of starting dilution for each of the test antibodies (dilu1, 3X final concentration). Use ADCC Assay Buffer to prepare and dilute antibodies in 1.5ml microcentrifuge tubes. Store the tubes containing the antibody starting dilutions appropriately before making antibody serial dilutions.

4.C. Plate Layout Design

We recommend orienting samples within an assay plate in a nonclustered fashion to help minimize any well positional effects on the response. For the protocol we describe here, use the plate layouts in Figure 7 as a guide.

The protocol uses serial replicate dilutions (n = 3) of reference antibody and each of the test antibodies to generate two 10-point dose response curves in each plate.

Recommended Plate Layout Design													
	1	2	3	4	5	6	7	8	9	10	11	12	
A	B	B	B	B	B	B	B	B	B	B	B	B	Assay Buffer (B)
B	B	no Ab	dilu9	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1	B	Reference Ab
C	B	no Ab	dilu9	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1	B	Test Ab
D	B	no Ab	dilu9	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1	B	Test Ab
E	B	no Ab	dilu9	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1	B	Reference Ab
F	B	no Ab	dilu9	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1	B	Reference Ab
G	B	no Ab	dilu9	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1	B	Test Ab
H	B	B	B	B	B	B	B	B	B	B	B	B	Assay Buffer (B)

Figure 7. Example plate layout showing nonclustered sample locations of reference antibody dilution series, a single test antibody dilution series, and “ADCC Assay Buffer” control, color coded for location.

4.D. Preparing Antibody Serial Dilutions

Preparing Serial Dilutions from a Single Antibody Dilution Stock to Generate Triplicates

Note: Alternatively, you can make three independent antibody dilution stocks to generate triplicates of each dose-response curve.

In order to establish a full dose-response range for any antibody to be tested in the ADCC Reporter Bioassay, we suggest that you first determine the starting concentrations and serial dilution schemes optimal for the antibody. For your reference, when tested in ADCC Reporter Bioassay, the starting concentrations and serial dilution schemes are 1×10^{-6} g/ml, fourfold serial dilution for rituximab, and 1×10^{-6} g/ml, threefold serial dilution for trastuzumab. These provide full dose response curves in both cases.

4.D. Preparing Antibody Serial Dilutions (continued)

1. Obtain a sterile clear V-bottom 96-well plate for preparing antibody serial dilutions. For threefold serial dilutions, perform the dilutions described in Steps 2–8 below. You will need 400µl of reference antibody at 3X the highest antibody concentration in your dose-response curve. You will need 200µl of each test antibody at 3X the highest antibody concentration in each of their dose-response curves. Adjust all volumes accordingly for other dilution schemes.
2. Add 150µl of reference antibody starting dilution (dilu1, 3X final concentration) to both well A11 and well B11.
3. Add 150µl of Test Antibody 1 and 150µl of Test Antibody 2 starting dilution (dilu1, 3X final concentration) to well E11 and well G11, respectively (see Figure 8).
4. Add 100µl of ADCC Assay Buffer to other wells in these four rows, from column 10 to column 2.
5. Using a multichannel pipette, transfer 50µl from the antibody starting dilutions in column 11 into column 10. Mix well by pipetting. Avoid creating bubbles.
6. Repeat equivalent threefold serial dilutions across columns from right to left until you reach column 3.
Note: Wells A2, B2, E2 and G2 contain 100µl of ADCC Assay Buffer as a “no-antibody” control.
7. Place the plate with antibody dilutions on the bench during preparation of target cells at the next step. Cover with a lid.

Recommended Plate Layout: Antibody Dilutions Prepared from a Single Antibody Stock													
	1	2	3	4	5	6	7	8	9	10	11	12	
A		no Ab	dilu9	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1		Reference Ab
B		no Ab	dilu9	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1		Reference Ab
C													
D													Reference Ab
E		no Ab	dilu9	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1		Test Ab 1
F													
G		no Ab	dilu9	dilu8	dilu7	dilu6	dilu5	dilu4	dilu3	dilu2	dilu1		Test Ab 2
H													

Figure 8. Example plate layout showing serial dilutions of antibodies. Reference and test antibodies for serial dilutions from a single antibody stock to generate triplicates.

4.E. Plating Target Cells

Recommendations for Plating Target Cells



Follow institutional guidelines for handling, including use of personal protective equipment (PPE), and waste disposal for hazardous material.

To prepare the target cells for use with the ADCC Reporter Bioassay Core Kit, the cells need to be cultured using standard practices to maintain viability of cells and density in a range satisfactory for good performance in a conventional ADCC assay. Several suspension target cell lines and adherent target cell lines grown in continuous culture have been tested in the ADCC Reporter Bioassay and demonstrated good results (see Figure 11 for examples).

As a possible alternative that may fit your needs, we have identified appropriate cell growth and freezing conditions that allow specific target cells to be used directly in bioassay without cell culture after thaw. Two suspension B-cell lines (WIL2-S and Raji) have been prepared in frozen, thaw-and-use formats and have demonstrated good results in ADCC Reporter Bioassay. These materials are currently available from Promega in ADCC Reporter Bioassay, Complete Kits (Cat.# G7014 and G7015, respectively). If your target is expressed on these cells you have the option to use one of the frozen, thaw-and-use target cell lines as provided in these Complete kits.

For assay optimization, try different effector-to-target cell (E:T) ratios in the range of 2.5:1 to 25:1. Keep the cell density of ADCC Bioassay Effector Cells constant, and vary the cell density of target cells. As a reference, we use 75,000 cells per well for ADCC Bioassay Effector Cells and an E:T ratio of 6:1 when working with ADCC Bioassay Target Cells and an anti-CD20 antibody.

Preparation of suspension target cell lines from continuous culture: On the day of assay, first estimate the target cell numbers needed. Harvest enough target cells (two to three times the required cell number) by centrifugation at $130\text{--}200 \times g$ for 10 minutes, wash cells once with 10ml of D-PBS, and recentrifuge. Resuspend cells in 6–9ml of ADCC Assay Buffer (prewarmed to 37°C) at a cell density that results in the required target cell number in 25µl in each assay well. Transfer the cell suspension to a sterile reagent reservoir. Using a multichannel pipette immediately add 25µl of cells to the inner 60 wells in both white 96-well assay plates. Dispense 75µl of ADCC Assay Buffer into outermost wells, labeled “B” in Figure 7. Cover plates with lids, and keep the plates on the bench before adding antibody dilutions and ADCC Bioassay Effector Cells.

Preparation of adherent target cell lines from continuous culture: Twenty to twenty four hours before the assay, remove cells from propagation flasks by trypsinization (or other standard procedure), and prepare them for plating in 96-well plates using fresh culture medium. Resuspend cells at an appropriate cell density, so that there will be the appropriate cell number required for each well in the ADCC bioassay when you plate 100µl cells per well. Transfer the cells to a sterile reagent reservoir and immediately add 100µl of cells to the inner 60 wells of white 96-well assay plates using a multichannel pipette. Dispense 100µl of culture medium into those outermost wells, labeled “B” in Figure 7, of both assay plates. Place lids on the plates and incubate overnight in a CO₂ incubator at 37°C.

On the morning of the assay, use a multichannel pipette to remove 95µl of culture medium from each of the wells. Add 25µl per well of ADCC Assay Buffer (prewarmed to 37°C) to the inner 60 wells of both assay plates. Always allow the pipette tips to touch the wall of the well, and add buffer gently to the wells to minimize disruption of cells. Dispense 75µl of ADCC Assay Buffer into those outermost wells, labeled “B” in Figure 7.

Note: This ADCC Reporter Bioassay protocol recommends using ADCC Assay Buffer, which contains 4% Low IgG Serum (FBS, low IgG). If you experience any cell viability or assay performance issues with the recommended ADCC Assay Buffer, we suggest that you test several other serum concentrations (in the range of 1–10%) to determine the optimal serum concentration for your test antibody and target cells.

4.F. Adding Antibody to Target Cells in Assay Plates

1. Using a multichannel pipette, add 25µl per well of antibody dilution series from the antibody dilution plates you prepared in Section 4.D to the white, 96-well assay plates already containing target cells, according to the plate layout in Figure 7.
2. Cover plates with lids and keep the plates on the bench before adding ADCC Bioassay Effector Cells at the next step.

4.G. Plating ADCC Bioassay Effector Cells

1. Label a sterile 15ml conical tube, “ADCC Bioassay Effector Cells”. Add 3.6ml of ADCC Assay Buffer (prewarmed to 37°C) to the tube.
 2. Remove 1 vial of ADCC Bioassay Effector Cells from –140°C freezer storage or vapor phase of liquid nitrogen to dry ice for transport to the bench on day of use. Thaw vial in a 37°C water bath until cells are **just** thawed (about 2–3 minutes). While thawing, gently agitate and visually inspect. Do not invert.
-  **Note:** The recommended thawing protocol above is important to the performance of the cells. No further handling is required or recommended.
3. Gently mix the cell suspension by pipetting 1 or 2 times. Transfer 630µl of cells to the 15ml “ADCC Bioassay Effector Cells” tube containing 3.6ml of ADCC Assay Buffer. Mix well by gently inverting the tube 2 times.
 4. Transfer cell suspension to a sterile reagent reservoir. Immediately, using a multichannel pipette, add 25µl of cells to the inner 60 wells of the 96-well assay plates already containing target cells and antibody.
 5. Cover plates with lids, and incubate the plates for 6 hours at 37°C in a humidified CO₂ incubator. Do not stack plates within the incubator.

4.H. Adding Bio-Glo™ Luciferase Assay Reagent

1. Remove assay plates from the 37°C incubator and equilibrate to ambient temperature (22–25°C) on the bench for 15 minutes.
2. Using a manual multichannel pipette, add 75µl of Bio-Glo™ Luciferase Assay Reagent to all the inner 60 wells of the assay plates; avoid creating any bubbles.
Note: Bio-Glo™ Reagent should be at ambient temperature.
3. Add 75µl of Bio-Glo™ Luciferase Assay Reagent to well B1, C1 and D1 in each assay plate to determine plate background.
4. Incubate at ambient temperature for 5–30 minutes.
5. Measure luminescence using a plate reader with glow-type luminescence read capabilities.

4.I. Data Analysis

1. Determine Plate Background by calculating the average RLU from wells B1, C1 and D1.
2. Calculate Fold of Induction = $\text{RLU (induced-background)} / \text{RLU (no antibody control-background)}$
Note: When calculating Fold of Induction, if the sample RLUs are 100 times or higher than the plate background RLU, there is no need to subtract plate background from sample RLU.
3. Graph data as RLU versus Log_{10} [antibody] and Fold of Induction versus Log_{10} [antibody]. Fit curves and determine EC_{50} of antibody response using appropriate curve fitting software (such as GraphPad Prism® software).

5. Troubleshooting

For questions not addressed here, please contact your local Promega Branch Office or Distributor. Contact information available at: www.promega.com. E-mail: techserv@promega.com

Symptoms

Possible Causes and Comments

High background

As a bioluminescent assay, the ADCC Reporter Bioassay generally gives low assay background and high signal response. There are multiple possible causes for high background such as a matrix effect from assay buffer or antibody stock solution, signal crosstalk from neighboring wells due to use of unsuitable assay plates or improper settings for the detection instrument. See also the “matrix effect” comment below.

Poor or low luminescence measurements (RLU readout)

Choose a sensitive instrument designed for luminescence detection. Instruments primarily designed for fluorescence detection are not recommended. If you must use an instrument primarily designed for fluorescence detection, ensure no filters are used. Luminometers measure and report luminescence as relative values, and actual numbers will vary between instruments. See Section 3 for more recommendations on how to set up the luminometer.

An insufficient number of effector cells could lead to low RLU. Handle and plate the effector cells appropriately according to the instructions in this protocol to ensure that there are sufficient viable effector cells per well in the assay.

Low activity of Bio-Glo™ Luciferase Assay Reagent also leads to low RLU. Store and handle Bio-Glo™ Luciferase Assay Reagent appropriately according to the instructions in the protocol.

Weak ADCC response (see section below).

Possible issues with matrix effect

IgG, complement (or other components from serum) supernatant of phage display or hybridoma culture could nonspecifically impact antibody binding to the FcγRIIIa receptor or affect the NFAT-RE signaling pathway directly, and cause a matrix effect. Use Low IgG Serum or perform further dilution of antibody starting preparation to minimize any matrix effect. The use of heat-inactivated or Low IgG Serum for growth of target cells also helps.

5. Troubleshooting (continued)

Symptoms

Weak ADCC response

Possible Causes and Comments

Optimize the E:T ratio while keeping the effector cell number constant at 75,000 cells per well. Since the readout of the ADCC Reporter Bioassay is from the effector cells, improvement of the response can be achieved by increasing the number of target cells per well.

Make sure to use the optimal concentration range for the antibody, which can provide a full dose response with complete upper and lower asymptotes. Note that any EC_{50} of antibody in the ADCC Reporter Bioassay will not necessarily be the same as those from other ADCC bioassays, thus some adjustment on the antibody starting concentration and serial dilution schemes may be needed to achieve the maximal response in ADCC Reporter Bioassay.

Optimize assay incubation time within a range of 6–24 hours, and choose the incubation time that gives optimal ADCC response.

Optimize the composition of ADCC Assay Buffer by varying the concentration of Low IgG FBS in a range of 1–10%, and choose the serum concentration that gives the optimal ADCC response.

Will I see the same ranking of therapeutic Abs in the Promega ADCC Reporter Bioassay as in a classic ADCC bioassay?

The ADCC Reporter Bioassay and classic ADCC bioassays show the same expected relative potency differences for Ab variants known to differ in ADCC efficiencies. This has been observed in several different studies using antibodies that differ in glycosylation (including fucosylation) and amino acid sequence.

Note: Assays should be performed using conditions that can differentiate activities in the ranges expected. These conditions may not be the same for both assays.

Symptoms

EC₅₀ for Ab varies between classic ADCC bioassay and Promega ADCC Reporter Bioassay

Possible Causes and Comments

EC₅₀ refers to the concentration of the substance (mAb in this assay) that gives 50% of the maximal biological response. The EC₅₀ value is determined not only by the binding affinity of the antibody but also by the assay conditions used in that particular assay such as the E:T ratio, incubation time and assay buffer in the case of ADCC bioassays. The EC₅₀ value of any antibody can differ dramatically between different assays and is not an intrinsic property of the antibody. It is normal if the EC₅₀ value for an Ab differs between ADCC Reporter Bioassay and other ADCC bioassays.

Some of the cells in the vial provided with the kit were dead

There will be some dead cells in the vial upon thawing but we designed and tested the bioassay such that this will not affect performance of the bioassay as long as instructions for handling cells are followed carefully. Carefully follow the instructions on gentle thawing and handling of the cells as outlined in the protocol (Section 4.G). Use the thawed cells immediately in the assay.

6. References

1. Hogarth, P.M. and Pietersz, G.A. (2012) Fc receptor-targeted therapies for the treatment of inflammation, cancer and beyond. *Nat. Rev. Drug Discov.* **11**, 311–31.
2. Chung, S. *et al.* (2012) Quantitative evaluation of fucose reducing effects in a humanized antibody on Fcγ receptor binding and antibody-dependent cell-mediated cytotoxicity activities. *mAbs* **4**, 326–40.
3. Parekh, B.S. *et al.* (2012) Development and validation of an antibody-dependent cell-mediated cytotoxicity-reporter gene assay. *mAbs* **4**, 310–8.
4. Chung, S. *et al.* (2014) Characterization of in vitro antibody-dependent cell-mediated cytotoxicity activity of therapeutic antibodies — Impact of effector cells. *J. Immunol. Methods* **407**, 63–75.
5. Cheng, Z.J. *et al.* (2014) Development of a robust reporter-based ADCC assay with frozen, thaw-and-use cells to measure Fc effector function of therapeutic antibodies. *J. Immunol. Methods* **414**, 69–81.

7. Appendix: Representative Assay Results

The following data were generated using the ADCC Bioassay Effector Cells with the ADCC Bioassay Target Cells (Figures 9 and 10) in the thaw-and-use format or the ADCC Bioassay Effector Cells with adherent target cells grown in continuous culture prior to bioassay (Figure 11).




Figure 9. ADCC Reporter Bioassay response to rituximab (Panel A; trade name: RITUXAN® or Panel B; Control Ab, Anti-CD20). ADCC Bioassay Target Cells (WIL2-S) were incubated with a series of concentrations of rituximab or Control Ab, Anti-CD20, followed by addition of ADCC Bioassay Effector Cells. The E:T ratio was 6:1. After 6 hours of induction at 37°C, Bio-Glo™ Luciferase Assay Reagent was added and luminescence was determined using a GloMax®-Multi+ Luminometer. Data were fitted to a 4PL curve using GraphPad Prism® software. The EC₅₀ response determined was 5.5ng/ml for rituximab and 3.2ng/ml for the Control Ab, Anti-CD20.




Figure 10. ADCC Reporter Bioassay response to rituximab (Panel A; trade name: RITUXAN®), or Control Ab, Anti-CD20 (Panel B). ADCC Bioassay Target Cells (Raji) were incubated with a series of concentrations of rituximab or Control Ab, Anti-CD20, followed by addition of ADCC Bioassay Effector Cells. The E:T ratio was 6:1. After 6 hours of induction at 37°C, Bio-Glo™ Luciferase Assay Reagent was added and luminescence was determined using a GloMax®-Multi+ Luminometer. Data were fitted to a 4PL curve using GraphPad Prism® software. The EC₅₀ response determined was 30.4ng/ml for rituximab and 18.3ng/ml for Control Ab, Anti-CD20.




Figure 11. ADCC Reporter Bioassay response to trastuzumab (trade name: HERCEPTIN[®], Panel A) or cetuximab (trade name: ERBITUX[®], Panel B) using ADCC Bioassay Effector Cells and human breast cancer cell line SK-BR-3 (Panel A) or human epithelial carcinoma cell line A431 (Panel B) as target cells. Cultured SK-BR-3 cells or A431 cells were plated at the density of 5,000 cells per well in complete culture medium overnight before bioassay. On the day of bioassay, the medium was removed carefully, and then the series of concentrations of antibody were added to the cells, followed by addition of ADCC Bioassay Effector Cells. The E:T ratio was 15:1. After 6 hours of induction at 37°C, Bio-Glo[™] Luciferase Assay Reagent was added and luminescence determined using a GloMax[®]-Multi+ Luminometer. The data were fitted to a 4PL curve using GraphPad Prism[®] software. The EC₅₀ of response using trastuzumab/SK-BR-3 target cells was 10.0ng/ml. The EC₅₀ of response using cetuximab/A431 target cells was 15.6ng/ml.

8. Related Products

T Cell Activation Bioassays

Product	Size	Cat.#
T Cell Activation Bioassay (NFAT)	1 each	J1621
T Cell Activation Bioassay (NFAT) 5X	5 each	J1625
T Cell Activation Bioassay (NFAT), Propagation Model	1 each	J1601
T Cell Activation Bioassay (IL-2)	1 each	J1651
T Cell Activation Bioassay (IL-2) 5X	5 each	J1655
T Cell Activation Bioassay (IL-2), Propagation Model	1 each	J1631

Not for Medical Diagnostic Use.

Immune Checkpoint Bioassays

Product	Size	Cat.#
4-1BB Bioassay	1 each	JA2351
4-1BB Bioassay 5X	5 each	JA2355
4-1BB Bioassay, Propagation Model	1 each	J2332
FcγRIIb CHO-K1 Cells	1 each	JA2251
FcγRIIb CHO-K1 Cells 5X	5 each	JA2255
FcγRIIb CHO-K1 Cells, Propagation Model	1 each	J2232
PD-1/PD-L1 Blockade Bioassay	1 each	J1250
PD-1/PD-L1 Blockade Bioassay 5X	5 each	J1255
PD-L1 Negative Cells	1 each	J1191
CTLA-4 Blockade Bioassay	1 each	JA3001
CTLA-4 Blockade Bioassay 5X	5 each	JA3005
TIGIT Negative Cells	1 each	J1921
PD-1+TIGIT Combination Bioassay	1 each	J2211
PD-1+TIGIT Combination Bioassay, 5X	5 each	J2215
Control Ab, Anti-CTLA-4	100µg	JA1020
Control Ab, Anti-PD-1	100µg	J1201
Control Ab, Anti-TIGIT	100µg	J2051
Control Ab, Anti-4-1BB	50µg	K1161
Control Ab, Anti-CD20	5µg	GA1130
Control Ab, Anti-LAG-3	100µg	K1150

Not for Medical Diagnostic Use.

Additional kit formats are available.



Immune Checkpoint Bioassays (continued)

Product	Size	Cat.#
LAG-3/MHCII Blockade Bioassay	1 each	JA1111
LAG-3/MHCII Blockade Bioassay 5X	5 each	JA1115
LAG-3/MHCII Blockade Bioassay, Propagation Model	1 each	JA2355

Not for Medical Diagnostic Use.

Fc Effector Bioassays

Product	Size	Cat.#
ADCC Reporter Bioassay, Complete Kit (Raji)*	1 each	G7015
ADCC Reporter Bioassay, Target Kit (Raji)*	1 each	G7016
ADCC Reporter Bioassay, Core Kit*	1 each	G7010
ADCC Reporter Bioassay, F Variant, Core Kit**	1 each	G9790
Fc γ RIIa-H ADCP Reporter Bioassay, Complete Kit**	1 each	G9901
Fc γ RIIa-H ADCP Reporter Bioassay, Core Kit**	1 each	G9991

*For Research Use Only. Not for use in diagnostic procedures.

**Not for Medical Diagnostic Use.

Additional kit formats are available.

Note: Additional Bioassays are available from Promega Custom Assay Services. To view and order products from Custom Assay Services, see the Early Access listings at:

www.promega.com/applications/biologics-drug-discovery/functional-bioassays/target-pathway-assays/

or email: CAS@promega.com

Detection Reagent

Product	Size	Cat.#
Bio-Glo™ Luciferase Assay System	10ml	G7941
	100ml	G7940

Not for Medical Diagnostic Use.

Luminometers

Product	Size	Cat.#
GloMax® Navigator System	1 each	GM2000
GloMax® Discover System	1 each	GM3000
GloMax® Explorer System	1 each	GM3500

For Research Use Only. Not for use in diagnostic procedures.

9. Summary of Changes

The following changes were made to the 8/19 revision of this document:

1. Figures 9 and 10 were updated.
2. Section 6, References 4 and 5 were updated.
3. Section 8, Related Products was updated.

^(b)NOT FOR MEDICAL DIAGNOSTIC USE. FOR IN VITRO USE ONLY. BY USE OF THIS PRODUCT, RECIPIENT AGREES TO BE BOUND BY THE TERMS OF THIS LIMITED USE STATEMENT. If the recipient is not willing to accept the conditions of this limited use statement, and the product is unused, Promega will accept return of the unused product and provide the recipient with a full refund.

This product may not be further sold or transferred by the recipient and may be used only by the recipient, and then only for (1) research use, (2) drug discovery and development of biologic drugs, (3) quality assurance testing of biologic drugs, and (4) product release assays for biologic drugs. No other commercial use is allowed. "Commercial use" means any and all uses of this product by recipient for monetary or other consideration, including providing a service, information or data to unaffiliated third parties, and resale of this product for any use. Recipient has no right to propagate, modify, derivatize, genetically engineer or otherwise create variations of the Effector Cells or the luciferase gene stably transfected within the Effector Cells. In addition, recipient must use Bio-Glo™ Luciferase Assay System purchased from Promega Corporation for all determinations of luminescence activity of this product, or contact Promega to obtain a license for use of this product with reagents other than Promega's. PROMEGA MAKES NO REPRESENTATIONS OR WARRANTIES OF ANY KIND, EITHER EXPRESSED OR IMPLIED, INCLUDING AS TO MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE WITH REGARDS TO THIS PRODUCT. The terms of this agreement shall be governed under the laws of the State of Wisconsin, USA.

^(b)U.S. Pat. No. 8,008,006 and European Pat. No. 1341808.

^(d)Patent Pending.

^(d)Licensed from Lonza Cologne GmbH under U.S. Pat. Nos. 7,700,357, 8,192,990 and 8,003,389, European Pat. Nos. 1297119, 1522587, 1607484 and 1741778 and other pending and issued patents.

^(e)U.S. Pat. No. 10,077,244 and other patents.

© 2012–2019 Promega Corporation. All Rights Reserved.

CytoTox 96, GloMax and Instinct are registered trademarks of Promega Corporation. Bio-Glo is a trademark of Promega Corporation.

ERBITUX is a registered trademark of ImClone LLC, a wholly owned subsidiary of Eli Lilly and Company. GraphPad Prism is a registered trademark of GraphPad Software, Inc. JMP is a registered trademark of SAS Institute, Inc. HERCEPTIN and RITUXAN are registered trademarks of Genentech USA, Inc.

Products may be covered by pending or issued patents or may have certain limitations. Please visit our Web site for more information.

All prices and specifications are subject to change without prior notice.

Product claims are subject to change. Please contact Promega Technical Services or access the Promega online catalog for the most up-to-date information on Promega products.