

## 500 μL, Model 750 RN SYR,Large Removable NDL, 22 ga,2 in, point style 2

Specification Sheet (Part/REF # 80830)



Split/splitless injection consists of a heated chamber with a glass liner into which the sample is injected through the septum. A microsyringe is used to inject the sample through a rubber septum into a flash vaporizer chamber at the head of the column.

## **Product Specifications**

| Volume                | 500 μL                  |
|-----------------------|-------------------------|
| Syringe Type          | Microliter              |
|                       |                         |
| Syringe Series        | 700                     |
| Termination           | Removable Needle (RN)   |
| Needle                | Large RN                |
| Gauge                 | 22 gauge                |
| Point Style           | 2                       |
| Needle Length         | 2 inch (51 mm)          |
| Fluid Path            | SST, borosilicate glass |
| Autoclavable          | Yes                     |
| Gas Sterilizable      | Yes                     |
| Pressure              | 1000 psig (68.9 bar)    |
| Maximum Temperature   | 115°C / 239°F           |
| Minimum Temperature   | 10°C / 50°F             |
| Barrel Inner Diameter | 0.135 inches (3.26 mm)  |
|                       |                         |



US: +1-800-648-5950 EU: +40-356-635-055



| Specifications are subject to change without notice | Sp                      | ec. Version B |
|---|-------------------------|---------------|
| Sales Quantity                                      | 1 syringe per pack      |               |
| Std Needle Dead Volume                              | 6.81 μL                 |               |
| Needle Outer Diameter                               | 0.0283 inches (0.72 mm) |               |
| Needle Inner Diameter                               | 0.0161 inches (0.41 mm) |               |
| Barrel Outer Diameter                               | 0.305 inches (7.75 mm)  |               |
|   |                         |               |

https://www.hamiltoncompany.com/laboratory-products/syringes/80830



US: +1-800-648-5950 EU: +40-356-635-055