

Extension Tube

Variable Types Available for Different Use

Specifications

- Precisely designed parts enabling easy connection with no leakage Many types satisfying various requiremnts.
 Sterilization: ETO

Product Range

Model	Capacity (mL)	Tube dimension			Remarks	Application
		I.D. (mm)	O.D. (mm)	Length (mm)	Remarks	Application
EX-S	1.0	2.9	4.9	15	Larger tube diameter. With Luter slip connector.	Mainly dialysis use. Connection with AVF, PSV or retention cannula, etc.
EX-M	2.8	2.9	4.9	45		
EX-L	3.7	2.9	4.9	60		
EX-S(T)	1.1	2.9	4.9	15	Larger tube diameter. With Luter slip connector and 3-way cock.	
EX-M(T)	2.9	2.9	4.9	45		
EX-M(T)	3.8	2.9	4.9	60		

Package : 50 pieces \times 10 = 500 pieces