Sigma-Aldrich.

SAFETY DATA SHEET according to Regulation (EC) No. 1907/2006

Version 8.1 Revision Date 17.01.2024 Print Date 23.11.2024 GENERIC EU MSDS - NO COUNTRY SPECIFIC DATA - NO OEL DATA

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1	Product identifiers Product name	:	Lipid Peroxidation (MDA) Assay Kit
	Product Number Brand REACH No.	:	MAK568 Sigma A registration number is not available for this substance as the substance or its uses are exempted from registration, the annual tonnage does not require a registration or the registration is envisaged for a later registration deadline.
1.2	Relevant identified us	es	of the substance or mixture and uses advised against

Identified uses : Laboratory chemicals, Manufacture of substances

1.3 Details of the supplier of the safety data sheet ÷

Company

Emergency telephone 1.4

Emergency Phone # +(44)-870-8200418 (CHEMTREC (GB)) : +(353)-19014670 (CHEMTREC Ireland) 001-803-017-9114 (CHEMTREC India)

This is a summary SDS for a kit. For the full SDS for each of the components listed in section 16 please visit our website.

SECTION 2: Hazards identification 2.1 Classification of the substance or mixture Flammable liquids, (Category 2) H225: Highly flammable liquid and vapor. Skin corrosion, (Sub-category H314: Causes severe skin burns and eye 1C) damage. Serious eye damage, (Category H318: Causes serious eye damage. 1) Long-term (chronic) aquatic H411: Toxic to aquatic life with long lasting

Sigma- MAK568

The life science business of Merck operates as MilliporeSigma in the US and Canada



Page 1 of 11

effects.

2.2 Label elements

Labelling according Regulation (EC) No 1272/2008 Pictogram

Pictogram					
Signal Word	Danger				
Hazard Statements H225 H314 H411	Highly flammable liquid and vapor. Causes severe skin burns and eye damage. Toxic to aquatic life with long lasting effects.				
Precautionary Statements					
P210	Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.				
P233	Keep container tightly closed.				
P273	Avoid release to the environment.				
P280	Wear protective gloves/ protective clothing/ eye protection/ face protection.				
P303 + P361 + P353	IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.				
P305 + P351 + P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.				
Supplemental Hazard Statements	none				
Reduced Labeling (<= 12	Reduced Labeling (<= 125 ml)				

apeim 9(

Pictogram



Signal Word	Danger
Hazard Statements H314	Causes severe skin burns and eye damage.
Precautionary Statements P280	Wear protective gloves/ protective clothing/ eye protection/ face
P303 + P361 + P353	protection. IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.
P305 + P351 + P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
Supplemental Hazard Statements	none

2.3 Other hazards - none

Sigma- MAK568

Page 2 of 11

The life science business of Merck operates as MilliporeSigma in the US and Canada



SECTION 3: Composition/information on ingredients

Refer to component SDS

SECTION 4: First aid measures

Refer to component SDS

SECTION 5: Firefighting measures

Refer to component SDS

SECTION 6: Accidental release measures

Refer to component SDS

SECTION 7: Handling and storage

7.1 Precautions for safe handling

Advice on protection against fire and explosion

Keep away from open flames, hot surfaces and sources of ignition. Take precautionary measures against static discharge.

Hygiene measures

Immediately change contaminated clothing. Apply preventive skin protection. Wash hands and face after working with substance. For precautions see section 2.2.

7.2 Conditions for safe storage, including any incompatibilities

Storage conditions

Keep container tightly closed in a dry and well-ventilated place. Keep away from heat and sources of ignition.

Storage stabilityRecommended storage temperature -25 - -10 °C

Storage class

Storage class (TRGS 510): 3: Flammable liquids

7.3 Specific end use(s)

Apart from the uses mentioned in section 1.2 no other specific uses are stipulated

SECTION 8: Exposure controls/personal protection

Refer to component SDS

Sigma- MAK568

Page 3 of 11

The life science business of Merck operates as MilliporeSigma in the US and Canada



SECTION 9: Physical and chemical properties

Refer to component SDS

SECTION 10: Stability and reactivity

Refer to component SDS

SECTION 11: Toxicological information

Refer to component SDS

SECTION 12: Ecological information

Refer to component SDS

SECTION 13: Disposal considerations

Refer to component SDS

SECTION 14: Transport informat	ion	
14.1 UN number ADR/RID: 3316	IMDG: 3316	IATA: 3316
14.2 UN proper shipping name ADR/RID:CHEMICAL KITIMDG:CHEMICAL KITIATA:Chemical kit		
14.3 Transport hazard class(es) ADR/RID: 9	IMDG: 9	IATA: 9
14.4 Packaging group ADR/RID: II	IMDG: II	IATA: II
14.5 Environmental hazards ADR/RID: no	IMDG Marine pollutant: no	IATA: no
14.6 Special precautions for use	r	
Further information :	No data available	

Sigma- MAK568

The life science business of Merck operates as MilliporeSigma in the US and Canada



Page 4 of 11

SECTION 15: Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

This material safety data sheet complies with the requirements of Regulation (EC) No. 1907/2006.

Refer to component SDS

15.2 Chemical Safety Assessment

For this product a chemical safety assessment was not carried out

SECTION 16: Other information

Kit Components:			
Lysis buffer	Sigma	MAK568A	Eye Irrit. 2; Aquatic Chronic 3; H319, H412
Phosphotungstic Acid Solution	Sigma	MAK568B	Skin Corr. 1C; Eye Dam. 1; Aquatic Chronic 3; H314, H412
BHT Stock/BHT (X100)	Sigma	MAK568C	Flam. Liq. 2; Eye Irrit. 2; Aquatic Chronic 2; H225, H319, H411
TBA MDA standard	Sigma Sigma	MAK568D MAK568E	Flam. Liq. 3; H226
	-		• •

Sigma- MAK568

The life science business of Merck operates as MilliporeSigma in the US and Canada

Page 5 of 11

Full text of other abbreviations

ADN - European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways; ADR - Agreement concerning the International Carriage of Dangerous Goods by Road; AIIC - Australian Inventory of Industrial Chemicals; ASTM -American Society for the Testing of Materials; bw - Body weight; CMR - Carcinogen, Mutagen or Reproductive Toxicant; DIN - Standard of the German Institute for Standardisation; DSL - Domestic Substances List (Canada); ECx - Concentration associated with x% response; ELx - Loading rate associated with x% response; EmS -Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx -Concentration associated with x% growth rate response; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan); ISO - International Organisation for Standardization; KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; n.o.s. -Not Otherwise Specified; NO(A)EC - No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NZIOC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention; PBT - Persistent, Bioaccumulative and Toxic substance; PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; RID - Regulations concerning the International Carriage of Dangerous Goods by Rail; SADT - Self-Accelerating Decomposition Temperature; SDS -Safety Data Sheet; TCSI - Taiwan Chemical Substance Inventory; TECI - Thailand Existing Chemicals Inventory; TSCA - Toxic Substances Control Act (United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative

Classification of the m	Classification procedure:	
Flam. Liq.2	H225	Based on product data or assessment
Skin Corr.1C	H314	Calculation method
Eye Dam.1	H318	Calculation method
Aquatic Chronic2	H411	Calculation method

Further information

The information is believed to be correct but is not exhaustive and will be used solely as a guideline, which is based on current knowledge of the chemical substance or mixture and is applicable to appropriate safety precautions for the product. It does not represent any guarantee of the properties of the product. Sigma-Aldrich Corporation and its Affiliates shall not be held liable for any damage resulting from handling or from contact

Sigma- MAK568

The life science business of Merck operates as MilliporeSigma in the US and Canada



Page 6 of 11

with the above product. See www.sigma-aldrich.com and/or the reverse side of invoice or packing slip for additional terms and conditions of sale.

Copyright 2020 Sigma-Aldrich Co. LLC. License granted to make unlimited paper copies for internal use only.

The branding on the header and/or footer of this document may temporarily not visually match the product purchased as we transition our branding. However, all of the information in the document regarding the product remains unchanged and matches the product ordered. For further information please contact mlsbranding@sial.com.

Sigma- MAK568

The life science business of Merck operates as MilliporeSigma in the US and Canada



Page 7 of 11



Kit Components:

Lysis buffer					
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Sigma - MAK568A		Warning	Causes serious eye irritation. Harmful to aquatic life with long lasting effects.	Wash skin thoroughly after handling. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Dispose of contents/ container to an approved waste disposal plant. Avoid release to the environment. If eye irritation persists: Get medical advice/ attention. Wear eye protection/ face protection.	

Phosphotungstic Acid Solution						
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental	
identifier		Word	Statements	Statements	information	
Sigma -	No. of the second se	Danger	Causes severe	Avoid release to		
MAK568B			skin burns and	the		
	▼		eye damage.	environment. IF		
			Harmful to	ON SKIN (or		
			aquatic life with	hair): Take off		
			long lasting	immediately all		
			effects.	contaminated		
				clothing. Rinse		
				skin with water.		

Sigma- MAK568

The life science business of Merck operates as MilliporeSigma in the US and Canada



Page 8 of 11

Wear protective gloves/ protective clothing/ eye protection/ face protection. IF INHALED: Remove person to fresh air and
a POISON CENTER/ doctor. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Wash contaminated clothing before reuse.

BHT Stock/E		1			
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Sigma -		Danger	Highly	Keep away	
MAK568C	X3		flammable liquid	from heat, hot	
	l Å		and vapor.	surfaces,	
			Causes serious	sparks, open	
			eye irritation.	flames and	
			Toxic to aquatic	other ignition	
			life with long	sources. No	
			lasting effects.	smoking. IF IN	
			_	EYES: Rinse	
				cautiously with	
				water for	
				several minutes.	
				Remove contact	
				lenses, if	
				present and	
				easy to do.	
				Continue	
				rinsing. Keep	
				container tightly	

Sigma- MAK568

The life science business of Merck operates as MilliporeSigma in the US and Canada



Page 9 of 11


closed. Ground and bond container and receiving equipment. Use explosion- proof electrical/ ventilating/ lighting/ equipment. Avoid release to the
environment.

ТВА	ТВА						
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information		
Sigma - MAK568D					No hazard pictogram, no signal word, no hazard statement(s), no precautionary statement(s) required		

MDA standard							
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental		
identifier		Word	Statements	Statements	information		
Sigma - MAK568E		Warning	Flammable liquid and vapor.	Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Keep container tightly closed. Ground and bond container and receiving equipment. Use explosion- proof electrical/ ventilating/ lighting/ equipment. Use non- sparking tools. Take action to			

Sigma- MAK568

The life science business of Merck operates as MilliporeSigma in the US and Canada



	prevent static	
	discharges.	

Sigma- MAK568

The life science business of Merck operates as MilliporeSigma in the US and Canada

Page 11 of 11

