

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Version 8.5

Revision Date 16.04.2025

Print Date 22.04.2025

GENERIC EU MSDS - NO COUNTRY SPECIFIC DATA - NO OEL DATA

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1 Product identifiers

Product name : MILLIPLEX® map Human Aging Panel 1

Product Number : HAGE1MAG-20K

Catalogue No. : 6E8328

Brand : Millipore

REACH No. : This product is a mixture. REACH Registration Number see section 3.

1.2 Relevant identified uses of the substance or mixture and uses advised against

Identified uses : Biochemical research/analysis

1.3 Details of the supplier of the safety data sheet

Company :

1.4 Emergency telephone

Emergency Phone # : +(44)-870-8200418 (CHEMTREC (GB))
 +(353)-19014670 (CHEMTREC Ireland)
 001-803-017-9114 (CHEMTREC India)

This is a summary SDS for a kit. For the full SDS for each of the components listed in section 16 please visit our website.

SECTION 2: Hazards identification

2.1 Classification of the substance or mixture

Acute toxicity, (Category 4)	H302: Harmful if swallowed.
Acute toxicity, (Category 4)	H332: Harmful if inhaled.
Acute toxicity, (Category 4)	H312: Harmful in contact with skin.
Serious eye damage, (Category 1)	H318: Causes serious eye damage.
Skin sensitization, (Category 1)	H317: May cause an allergic skin reaction.
Specific target organ toxicity - repeated exposure, (Category 2), Respiratory Tract	H373: May cause damage to organs through prolonged or repeated exposure.



Long-term (chronic) aquatic hazard, (Category 2)

H411: Toxic to aquatic life with long lasting effects.

2.2 Label elements

Labelling according Regulation (EC) No 1272/2008

Pictogram



Signal Word

Danger

Hazard Statements

H302 + H312 + H332

Harmful if swallowed, in contact with skin or if inhaled.

H317

May cause an allergic skin reaction.

H318

Causes serious eye damage.

H373

May cause damage to organs (Respiratory Tract) through prolonged or repeated exposure.

H411

Toxic to aquatic life with long lasting effects.

Precautionary Statements

P273

Avoid release to the environment.

P280

Wear protective gloves/ protective clothing/ eye protection/ face protection.

P302 + P352 + P312

IF ON SKIN: Wash with plenty of water. Call a POISON CENTER/ doctor if you feel unwell.

P304 + P340 + P312

IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a POISON CENTER/ doctor if you feel unwell.

P305 + P351 + P338

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P314

Get medical advice/ attention if you feel unwell.

Supplemental Hazard Statements

none

Reduced Labeling (<= 125 ml)

Pictogram



Signal Word

Danger

Hazard Statements

H317

May cause an allergic skin reaction.

H318

Causes serious eye damage.

Precautionary Statements

P305 + P351 + P338

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Supplemental Hazard Statements

none

2.3 Other hazards

This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher.

Ecological information:

Millipore- HAGE1MAG-20K

Page 2 of 13

The life science business of Merck operates as MilliporeSigma in the US and Canada



The substance/mixture does not contain components considered to have endocrine disrupting properties according to REACH Article 57(f) or Commission Delegated regulation (EU) 2017/2100 or Commission Regulation (EU) 2018/605 at levels of 0.1% or higher.

Toxicological information:

The substance/mixture does not contain components considered to have endocrine disrupting properties according to REACH Article 57(f) or Commission Delegated regulation (EU) 2017/2100 or Commission Regulation (EU) 2018/605 at levels of 0.1% or higher.

SECTION 3: Composition/information on ingredients

Refer to component SDS

SECTION 4: First aid measures

Refer to component SDS

SECTION 5: Firefighting measures

Refer to component SDS

SECTION 6: Accidental release measures

Refer to component SDS

SECTION 7: Handling and storage

7.1 Precautions for safe handling

Advice on safe handling

Work under hood. Do not inhale substance/mixture. Avoid generation of vapours/aerosols.

Hygiene measures

Immediately change contaminated clothing. Apply preventive skin protection. Wash hands and face after working with substance.

For precautions see section 2.2.

7.2 Conditions for safe storage, including any incompatibilities

Storage conditions

Tightly closed.

Recommended storage temperature see product label.

Storage class

Storage class (TRGS 510): 10: Combustible liquids

7.3 Specific end use(s)

Apart from the uses mentioned in section 1.2 no other specific uses are stipulated

SECTION 8: Exposure controls/personal protection

Refer to component SDS



SECTION 15: Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

This material safety data sheet complies with the requirements of Regulation (EC) No. 1907/2006.

Refer to component SDS

15.2 Chemical Safety Assessment

For this product a chemical safety assessment was not carried out

SECTION 16: Other information

Kit Components:

10X Wash Buffer - MILLIPLEX	Millipore	L-WB	Skin Sens. 1; Aquatic Chronic 2; H317, H411
Milliplex Assay Buffer	Millipore	L-AB	
Bead Diluent MILLIPLEX	Millipore	LBD-4	
Streptavidin-Phycoerythrin	Millipore	L-SAPE4	Eye Irrit. 2; H319
Serum Matrix	Millipore	HAGE1-SM	Aquatic Chronic 3; H412
Human Aging Panel 1 Standard	Millipore	HAGE1-8020	Acute Tox. 4; Eye Dam. 1; STOT RE 2; Aquatic Chronic 3; H302 + H312 + H332, H318, H373, H412
Human Aging Panel 1 Quality Controls 1 & 2	Millipore	HAGE1-6020	Acute Tox. 4; Eye Dam. 1; STOT RE 2; Aquatic Chronic 3; H302 + H312 + H332, H318, H373, H412
Human Aging Panel 1 Detection Antibodies	Millipore	HAGE1-1020	Eye Irrit. 2; H319
15-CTACK Beads Magnetic	Millipore	HACTCK-MAG	
18-FGF21 Beads Magnetic	Millipore	HAFGF21-MAG	
20-GDF11 Beads Magnetic	Millipore	HAGDF11-MAG	
22-GDF15 Beads Magnetic	Millipore	HAGDF15-MAG	
26-GnRH Beads Magnetic	Millipore	HAGNRH-MAG	
34-anti IL-6 MAGNETIC Beads MILLIPLEX	Millipore	HIL6-MAG	
35-IL-10 Beads Magnetic	Millipore	HAIL10-MAG	
36-IL-18 Beads Magnetic	Millipore	HAIL18-MAG	
37-Jag1 Beads Magnetic	Millipore	HAJAG1-MAG	
39-anti Leptin MAGNETIC Beads MILL	Millipore	HLPTN-MAG	



Full text of other abbreviations

ADN - European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways; ADR - Agreement concerning the International Carriage of Dangerous Goods by Road; AIIC - Australian Inventory of Industrial Chemicals; ASTM - American Society for the Testing of Materials; bw - Body weight; CMR - Carcinogen, Mutagen or Reproductive Toxicant; DIN - Standard of the German Institute for Standardisation; DSL - Domestic Substances List (Canada); ECx - Concentration associated with x% response; ELx - Loading rate associated with x% response; EmS - Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx - Concentration associated with x% growth rate response; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan); ISO - International Organisation for Standardization; KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; n.o.s. - Not Otherwise Specified; NO(A)EC - No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention; PBT - Persistent, Bioaccumulative and Toxic substance; PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; RID - Regulations concerning the International Carriage of Dangerous Goods by Rail; SADT - Self-Accelerating Decomposition Temperature; SDS - Safety Data Sheet; TCSI - Taiwan Chemical Substance Inventory; TECI - Thailand Existing Chemicals Inventory; TSCA - Toxic Substances Control Act (United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative

Classification of the mixture

Acute Tox.4	H302
Acute Tox.4	H332
Acute Tox.4	H312
Eye Dam.1	H318
Skin Sens.1	H317
STOT RE2	H373
Aquatic Chronic2	H411

Classification procedure:

Calculation method
Calculation method
Calculation method
Calculation method
Calculation method
Calculation method
Calculation method



Further information


The information is believed to be correct but is not exhaustive and will be used solely as a guideline, which is based on current knowledge of the chemical substance or mixture and is applicable to appropriate safety precautions for the product. It does not represent any guarantee of the properties of the product. Sigma-Aldrich Corporation and its Affiliates shall not be held liable for any damage resulting from handling or from contact with the above product. See www.sigma-aldrich.com and/or the reverse side of invoice or packing slip for additional terms and conditions of sale.

Copyright 2020 Sigma-Aldrich Co. LLC. License granted to make unlimited paper copies for internal use only.

The branding on the header and/or footer of this document may temporarily not visually match the product purchased as we transition our branding. However, all of the information in the document regarding the product remains unchanged and matches the product ordered. For further information please contact mlsbranding@sial.com.




Kit Components:

10X Wash Buffer - MILLIPLEX					
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information
Millipore - L-WB		Warning	<p>May cause an allergic skin reaction.</p> <p>Toxic to aquatic life with long lasting effects.</p>	<p>Avoid breathing mist or vapors.</p> <p>IF ON SKIN:</p> <p>Wash with plenty of water.</p> <p>Contaminated work clothing should not be allowed out of the workplace.</p> <p>If skin irritation or rash occurs: Get medical advice/ attention.</p> <p>Avoid release to the environment.</p> <p>Wear protective gloves.</p>	


Milliplex Assay Buffer					
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information
Millipore - L-AB					No hazard pictogram, no signal word, no hazard statement(s), no precautionary statement(s) required.

Bead Diluent MILLIPLEX					
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information
Millipore - LBD-4					No hazard pictogram, no signal word, no hazard statement(s), no precautionary statement(s) required.



Streptavidin-Phycoerythrin					
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information
Millipore - L-SAPE4		Warning	Causes serious eye irritation.	Wash skin thoroughly after handling. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Wear eye protection/ face protection. If eye irritation persists: Get medical advice/ attention.	


Serum Matrix					
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information
Millipore - HAGE1-SM			Harmful to aquatic life with long lasting effects.	Avoid release to the environment. Dispose of contents/ container to an approved waste disposal plant.	

Human Aging Panel 1 Standard					
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information
Millipore - HAGE1-8020		Danger	Harmful if swallowed, in contact with skin or if inhaled. Causes serious eye damage. May cause damage to organs (Respiratory Tract) through prolonged or repeated exposure. Harmful to aquatic life with long lasting	Avoid release to the environment. IF ON SKIN: Wash with plenty of water. Call a POISON CENTER/ doctor if you feel unwell. Wear protective gloves/ protective clothing/ eye protection/ face protection. IF INHALED:	



			effects.	Remove person to fresh air and keep comfortable for breathing. Call a POISON CENTER/ doctor if you feel unwell. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Get medical advice/ attention if you feel unwell.	
--	--	--	----------	--	--


Human Aging Panel 1 Quality Controls 1 & 2

Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information
Millipore - HAGE1-6020		Danger	<p>Harmful if swallowed, in contact with skin or if inhaled. Causes serious eye damage. May cause damage to organs (Respiratory Tract) through prolonged or repeated exposure. Harmful to aquatic life with long lasting effects.</p>	<p>Avoid release to the environment. IF ON SKIN: Wash with plenty of water. Call a POISON CENTER/ doctor if you feel unwell. Wear protective gloves/ protective clothing/ eye protection/ face protection. IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a POISON CENTER/ doctor if you feel unwell. IF IN EYES: Rinse cautiously with water for</p>	



				several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Get medical advice/attention if you feel unwell.	
--	--	--	--	---	--

Human Aging Panel 1 Detection Antibodies

Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information
Millipore - HAGE1-1020		Warning	Causes serious eye irritation.	Wash skin thoroughly after handling. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Wear eye protection/face protection. If eye irritation persists: Get medical advice/attention.	

15-CTACK Beads Magnetic

Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information
Millipore - HACTCK-MAG					No hazard pictogram, no signal word, no hazard statement(s), no precautionary statement(s) required.

18-FGF21 Beads Magnetic

Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information
Millipore - HAFGF21-MAG					No hazard pictogram, no signal word, no hazard statement(s), no precautionary



					statement(s) required.
--	--	--	--	--	------------------------

20-GDF11 Beads Magnetic

Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information
Millipore - HAGDF11-MAG					No hazard pictogram, no signal word, no hazard statement(s), no precautionary statement(s) required.

22-GDF15 Beads Magnetic

Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information
Millipore - HAGDF15-MAG					No hazard pictogram, no signal word, no hazard statement(s), no precautionary statement(s) required.

26-GnRH Beads Magnetic

Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information
Millipore - HAGNRH-MAG					Not a hazardous substance or mixture according to Regulation (EC) No 1272/2008.

34-anti IL-6 MAGNETIC Beads MILLIPLEX

Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information
Millipore - HIL6-MAG					No hazard pictogram, no signal word, no hazard statement(s), no precautionary statement(s) required.

35-IL-10 Beads Magnetic

Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information
Millipore - HAIL10-MAG					No hazard pictogram, no



					signal word, no hazard statement(s), no precautionary statement(s) required.
--	--	--	--	--	--

36-IL-18 Beads Magnetic

Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information
Millipore - HAIL18-MAG					Not a hazardous substance or mixture according to Regulation (EC) No 1272/2008.

37-Jag1 Beads Magnetic

Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information
Millipore - HAJAG1-MAG					Not a hazardous substance or mixture according to Regulation (EC) No 1272/2008.

39-anti Leptin MAGNETIC Beads MILL

Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information
Millipore - HLPTN-MAG					No hazard pictogram, no signal word, no hazard statement(s), no precautionary statement(s) required.

