

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Version 8.4 Revision Date 24.04.2025 Print Date 27.04.2025

GENERIC EU MSDS - NO COUNTRY SPECIFIC DATA - NO OEL DATA

SECTION 1: Identification of the substance/mixture and of the company/undertaking

Product identifiers

Product name : Canine Cytokine Magnetic Bead Panel

: CCYTOMAG-90K Product Number

: 667593 Catalogue No. Brand : Millipore

REACH No. : This product is a mixture. REACH Registration Number see

section 3.

1.2 Relevant identified uses of the substance or mixture and uses advised against

Identified uses : Biochemical research/analysis

1.3 Details of the supplier of the safety data sheet

Company

1.4 Emergency telephone

Emergency Phone # : +(44)-870-8200418 (CHEMTREC (GB))

> +(353)-19014670 (CHEMTREC Ireland) 001-803-017-9114 (CHEMTREC India)

This is a summary SDS for a kit. For the full SDS for each of the components listed in section 16 please visit our website.

SECTION 2: Hazards identification

2.1 Classification of the substance or mixture

Acute toxicity, (Category 4) H302: Harmful if swallowed.

H332: Harmful if inhaled. Acute toxicity, (Category 4)

Acute toxicity, (Category 4) H312: Harmful in contact with skin.

Serious eye damage, (Category H318: Causes serious eye damage.

Skin sensitization, (Category 1) H317: May cause an allergic skin reaction.

Specific target organ toxicity -H373: May cause damage to organs repeated exposure, (Category 2),

Respiratory Tract

through prolonged or repeated exposure.

Millipore- CCYTOMAG-90K Page 1 of 14



Long-term (chronic) aquatic hazard, (Category 2)

H411: Toxic to aquatic life with long lasting effects.

2.2 Label elements

Labelling according Regulation (EC) No 1272/2008

Pictogram

Signal Word Danger

Hazard Statements

H302 + H312 + H332 Harmful if swallowed, in contact with skin or if inhaled.

H317 May cause an allergic skin reaction.

H318 Causes serious eye damage.

H373 May cause damage to organs (Respiratory Tract) through

prolonged or repeated exposure.

H411 Toxic to aquatic life with long lasting effects.

Precautionary Statements

P273 Avoid release to the environment.

P280 Wear protective gloves/ protective clothing/ eye protection/ face

protection.

P302 + P352 + P312 IF ON SKIN: Wash with plenty of water.Call a POISON CENTER/

doctor if you feel unwell.

P304 + P340 + P312 IF INHALED: Remove person to fresh air and keep comfortable

for breathing. Call a POISON CENTER/ doctor if you feel unwell.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue

rinsing.

P314 Get medical advice/ attention if you feel unwell.

Supplemental Hazard

Statements

none

Reduced Labeling (<= 125 ml)

Pictogram



Signal Word Danger

Hazard Statements

H317 May cause an allergic skin reaction. H318 Causes serious eye damage.

Precautionary Statements

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue

rinsing.

Supplemental Hazard

Statements

none

2.3 Other hazards

This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher.

Ecological information:

Millipore- CCYTOMAG-90K Page 2 of 14



The substance/mixture does not contain components considered to have endocrine disrupting properties according to REACH Article 57(f) or Commission Delegated regulation (EU) 2017/2100 or Commission Regulation (EU) 2018/605 at levels of 0.1% or higher. Toxicological information:

The substance/mixture does not contain components considered to have endocrine disrupting properties according to REACH Article 57(f) or Commission Delegated regulation (EU) 2017/2100 or Commission Regulation (EU) 2018/605 at levels of 0.1% or higher.

SECTION 3: Composition/information on ingredients

Refer to component SDS

SECTION 4: First aid measures

Refer to component SDS

SECTION 5: Firefighting measures

Refer to component SDS

SECTION 6: Accidental release measures

Refer to component SDS

SECTION 7: Handling and storage

7.1 Precautions for safe handling

Advice on safe handling

Work under hood. Do not inhale substance/mixture. Avoid generation of vapours/aerosols.

Hygiene measures

Immediately change contaminated clothing. Apply preventive skin protection. Wash hands and face after working with substance.

For precautions see section 2.2.

7.2 Conditions for safe storage, including any incompatibilities

Storage conditions

Tightly closed.

Recommended storage temperature see product label.

Storage class

Storage class (TRGS 510): 10: Combustible liquids

7.3 Specific end use(s)

Apart from the uses mentioned in section 1.2 no other specific uses are stipulated

SECTION 8: Exposure controls/personal protection

Refer to component SDS

Millipore- CCYTOMAG-90K Page 3 of 14



SECTION 9: Physical and chemical properties

Refer to component SDS

SECTION 10: Stability and reactivity

Refer to component SDS

SECTION 11: Toxicological information

Refer to component SDS

SECTION 12: Ecological information

Refer to component SDS

SECTION 13: Disposal considerations

Refer to component SDS

SECTION 14: Transport information

14.1 UN number

ADR/RID: - IMDG: - IATA: -

14.2 UN proper shipping name

ADR/RID: Not dangerous goods IMDG: Not dangerous goods IATA: Not dangerous goods

14.3 Transport hazard class(es)

ADR/RID: - IMDG: - IATA: -

14.4 Packaging group

ADR/RID: - IMDG: - IATA: -

14.5 Environmental hazards

ADR/RID: no IMDG Marine pollutant: no IATA: no

14.6 Special precautions for user

Further information

Not classified as dangerous in the meaning of transport regulations.

14.7 Maritime transport in bulk according to IMO instruments

Not applicable for product as supplied.

Millipore- CCYTOMAG-90K Page 4 of 14



SECTION 15: Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

This material safety data sheet complies with the requirements of Regulation (EC) No. 1907/2006.

Refer to component SDS

15.2 Chemical Safety Assessment

For this product a chemical safety assessment was not carried out

SECTION 1	L6: Other	information
------------------	-----------	-------------

Kit Components: Canine Cytokine Stndrd CocktailMILLIPLEX	Millipore	LCC-8090	Acute Tox. 4; Eye Dam. 1; STOT RE 2; Aquatic Chronic 3; H302 + H312 + H332, H318, H373, H412
Canine Cytokine Quality Controls 1 and 2 - MILLIPLEX.	Millipore	LCC-6090	Acute Tox. 4; Eye Dam. 1; STOT RE 2; Aquatic Chronic 3; H302 + H312 + H332, H318, H373, H412
Serum Matrix - MILLIPLEX - COLD PAC	Millipore	LMC-SD	Aquatic Chronic 3; H412
Milliplex Assay Buffer	Millipore	L-AB	
10X Wash Buffer - MILLIPLEX	Millipore	L-WB	Skin Sens. 1; Aquatic Chronic 2; H317, H411
Canine Cytokine Detection Antibodies	Millipore	LCC-1090	Eye Irrit. 2; H319
Streptavidin- Phycoerythrin	Millipore	L-SAPE3	Eye Irrit. 2; H319
14-GM-CSF Beads - Magnetic - MILLIPLEX	Millipore	CGMCSF-MAG	
25-IFNy Beads - Magnetic - MILLIPLEX	Millipore	CIFNG-MAG	
61-anti KC MAGNETIC Beads MILLIPLEX	Millipore	MKC-MAG	
57-anti IP-10 MAGNETIC Beads MILLIPLEX	Millipore	MIP10-MAG	
29-IL-2 Beads - Magnetic - MILLIPLEX	Millipore	CIL2-MAG	
35-IL-6 Beads - Magnetic - MILLIPLEX	Millipore	CIL6-MAG	
36-anti IL-7 MAGNETIC Beads MILLIPLEX	Millipore	MIL7-MAG	
48-IL-8 Beads - Magnetic - MILLIPLEX	Millipore	CIL8-MAG	
63-IL-10 Beads - Magnetic - MILLIPLEX	Millipore	CIL10-MAG	
54-anti IL-15 MAGNETIC Beads MILLIPLEX	Millipore	MIL15-MAG	
67-IL-18 Beads - Magnetic - MILLIPLEX	Millipore	CIL18-MAG	
74-MCP-1 Beads - Magnetic - MILLIPLEX	Millipore	CMCP1-MAG	

Millipore- CCYTOMAG-90K Page 5 of 14



Magnetic - MILLIPLEX

Full text of other abbreviations

ADN - European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways; ADR - Agreement concerning the International Carriage of Dangerous Goods by Road; AIIC - Australian Inventory of Industrial Chemicals; ASTM -American Society for the Testing of Materials; bw - Body weight; CMR - Carcinogen, Mutagen or Reproductive Toxicant; DIN - Standard of the German Institute for Standardisation; DSL - Domestic Substances List (Canada); ECx - Concentration associated with x% response; ELx - Loading rate associated with x% response; EmS -Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx -Concentration associated with x% growth rate response; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan); ISO - International Organisation for Standardization; KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; n.o.s. -Not Otherwise Specified; NO(A)EC - No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention; PBT - Persistent, Bioaccumulative and Toxic substance; PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; RID - Regulations concerning the International Carriage of Dangerous Goods by Rail; SADT - Self-Accelerating Decomposition Temperature; SDS -Safety Data Sheet; TCSI - Taiwan Chemical Substance Inventory; TECI - Thailand Existing Chemicals Inventory; TSCA - Toxic Substances Control Act (United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative

Classification of the mixture		Classification procedure:
Acute Tox.4	H302	Calculation method
Acute Tox.4	H332	Calculation method
Acute Tox.4	H312	Calculation method
Eye Dam.1	H318	Calculation method
Skin Sens.1	H317	Calculation method
STOT RE2	H373	Calculation method
Aquatic Chronic2	H411	Calculation method

Further information

The information is believed to be correct but is not exhaustive and will be used solely as a quideline, which is based on current knowledge of the chemical substance or mixture and is applicable to appropriate safety precautions for the product. It does not represent any quarantee of the properties of the product. Sigma-Aldrich Corporation and its Affiliates shall not be held liable for any damage resulting from handling or from contact

Millipore- CCYTOMAG-90K Page 6 of 14



with the above product. See www.sigma-aldrich.com and/or the reverse side of invoice or packing slip for additional terms and conditions of sale.

Copyright 2020 Sigma-Aldrich Co. LLC. License granted to make unlimited paper copies for internal use only.

The branding on the header and/or footer of this document may temporarily not visually match the product purchased as we transition our branding. However, all of the information in the document regarding the product remains unchanged and matches the product ordered. For further information please contact mlsbranding@sial.com.

Millipore- CCYTOMAG-90K Page 7 of 14





Kit Components:

Canine Cyto	kine Stndrd C	ocktailMIL	LIPLEX		
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Millipore -		Danger	Harmful if	Avoid release to	
LCC-8090			swallowed, in	the	
	X		contact with skin	environment. IF	
	不多		or if inhaled.	ON SKIN: Wash	
	X		Causes serious	with plenty of	
	(!)		eye damage.	water.Call a	
			May cause	POISON	
			damage to	CENTER/ doctor	
			organs	if you feel	
			(Respiratory	unwell. Wear	
			Tract) through	protective	
			prolonged or	gloves/	
			repeated	protective	
			exposure.	clothing/ eye	
			Harmful to	protection/ face	
			aquatic life with	protection. IF	
			long lasting	INHALED:	
			effects.	Remove person	
				to fresh air and	
				keep	
				comfortable for	
				breathing. Call a POISON	
				CENTER/ doctor	
				if you feel	
				unwell. IF IN	
				EYES: Rinse	
				cautiously with	
				water for	
				several minutes.	
				Remove contact	
				lenses, if	
				present and	
				easy to do.	
				Continue	
				rinsing. Get	
				medical advice/	
				attention if you	
				feel unwell.	

Canine Cytokine Quality Controls 1 and 2 - MILLIPLEX.						
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information	
Millipore -		Danger	Harmful if swallowed, in	Avoid release to the		

Millipore- CCYTOMAG-90K Page 8 of 14



1.00.000	_			
LCC-6090	JE JE	contact with skin	environment. IF	
		or if inhaled.	ON SKIN: Wash	
		Causes serious	with plenty of	
		eye damage.	water.Call a	
		May cause	POISON	
		damage to	CENTER/ doctor	
		organs	if you feel	
		(Respiratory	unwell. Wear	
		Tract) through	protective	
		prolonged or	gloves/	
		repeated	protective	
		exposure.	clothing/ eye	
		Harmful to	protection/ face	
			protection, lace	
		aquatic life with	INHALED:	
		long lasting		
		effects.	Remove person	
			to fresh air and	
			keep	
			comfortable for	
			breathing. Call a	
			POISON	
			CENTER/ doctor	
			if you feel	
			unwell. IF IN	
			EYES: Rinse	
			cautiously with	
			water for	
			several minutes.	
			Remove contact	
			lenses, if	
			present and	
			easy to do.	
			Continue	
			rinsing. Get	
			medical advice/	
			attention if you	
			feel unwell.	

Serum Matrix	Serum Matrix - MILLIPLEX - COLD PAC						
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental		
identifier		Word	Statements	Statements	information		
Millipore -			Harmful to	Avoid release to			
LMC-SD			aquatic life with	the			
			long lasting	environment.			
			effects.	Dispose of			
				contents/			
				container to an			
				approved waste			
				disposal plant.			

Milliplex Assay Buffer						
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental	
identifier		Word	Statements	Statements	information	
Millipore - L-					No hazard	
					pictogram, no	

Millipore- CCYTOMAG-90K Page 9 of 14



AB			signal word, no
			hazard
			statement(s), no
			precautionary
			statement(s)
			required.

10X Wash Buffer - MILLIPLEX						
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information	
Millipore - L- WB	1	Warning	May cause an allergic skin reaction. Toxic to aquatic life with long lasting effects.	Avoid breathing mist or vapors. IF ON SKIN: Wash with plenty of water. Contaminated work clothing should not be allowed out of the workplace. If skin irritation or rash occurs: Get medical advice/ attention. Avoid release to the environment. Wear protective gloves.		

Canine Cytokine Detection Antibodies							
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental		
identifier		Word	Statements	Statements	information		
Millipore - LCC-1090		Warning	Causes serious eye irritation.	Wash skin thoroughly after handling. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Wear eye protection/ face protection. If eye irritation persists: Get medical advice/ attention.			

Millipore- CCYTOMAG-90K Page 10 of 14



The life science business of Merck operates as Millipore Sigma in the US and Canada $\,$

Streptavidin-Phycoerythrin							
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental		
identifier		Word	Statements	Statements	information		
Millipore - L-		Warning	Causes serious	Wash skin			
SAPE3	•		eye irritation.	thoroughly after			
	•			handling. IF IN			
				EYES: Rinse			
				cautiously with			
				water for			
				several minutes.			
				Remove contact			
				lenses, if			
				present and			
				easy to do.			
				Continue			
				rinsing. Wear			
				eye protection/			
				face protection.			
				If eye irritation			
				persists: Get			
				medical advice/			
				attention.			

14-GM-CSF B	14-GM-CSF Beads - Magnetic - MILLIPLEX								
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental				
identifier		Word	Statements	Statements	information				
Millipore -					No hazard				
CGMCSF-MAG					pictogram, no				
					signal word, no				
					hazard				
					statement(s), no				
					precautionary				
					statement(s)				
					required.				

25-IFNy Bea	ds - Magnetic	- MILLIP	LEX		
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Millipore -					No hazard
CIFNG-MAG					pictogram, no
					signal word, no
					hazard
					statement(s), no
					precautionary
					statement(s)
					required.

61-anti KC M	61-anti KC MAGNETIC Beads MILLIPLEX							
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental			
identifier		Word	Statements	Statements	information			
Millipore -					No hazard			
MKC-MAG					pictogram, no			
					signal word, no			
					hazard			
					statement(s), no			

Millipore- CCYTOMAG-90K Page 11 of 14



						precautionary statement(s) required.
--	--	--	--	--	--	--

57-anti IP-10	57-anti IP-10 MAGNETIC Beads MILLIPLEX									
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental					
identifier		Word	Statements	Statements	information					
Millipore -					No hazard					
MIP10-MAG					pictogram, no					
					signal word, no					
					hazard					
					statement(s), no					
					precautionary					
					statement(s)					
					required.					

29-IL-2 Bead	29-IL-2 Beads - Magnetic - MILLIPLEX								
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental				
identifier		Word	Statements	Statements	information				
Millipore -					Not a hazardous				
CIL2-MAG					substance or				
					mixture				
					according to				
					Regulation (EC)				
					No 1272/2008.				

35-IL-6 Beads - Magnetic - MILLIPLEX								
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental			
identifier		Word	Statements	Statements	information			
Millipore -					Not a hazardous			
CIL6-MAG					substance or			
					mixture			
					according to			
					Regulation (EC)			
					No 1272/2008.			

36-anti IL-7	36-anti IL-7 MAGNETIC Beads MILLIPLEX							
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information			
Millipore - MIL7-MAG					No hazard pictogram, no signal word, no hazard statement(s), no precautionary statement(s) required.			

48-IL-8 Beads - Magnetic - MILLIPLEX								
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental			
identifier		Word	Statements	Statements	information			
Millipore -					Not a hazardous			
CIL8-MAG					substance or			
					mixture			

Millipore- CCYTOMAG-90K Page 12 of 14

A

		according to
		Regulation (EC)
		No 1272/2008.

63-IL-10 Beads - Magnetic - MILLIPLEX								
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental			
identifier		Word	Statements	Statements	information			
Millipore -					Not a hazardous			
CIL10-MAG					substance or			
					mixture			
					according to			
					Regulation (EC)			
					No 1272/2008.			

54-anti IL-15	54-anti IL-15 MAGNETIC Beads MILLIPLEX									
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental					
identifier		Word	Statements	Statements	information					
Millipore -					No hazard					
MIL15-MAG					pictogram, no					
					signal word, no					
					hazard					
					statement(s), no					
					precautionary					
					statement(s)					
					required.					

67-IL-18 Beads - Magnetic - MILLIPLEX					
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Millipore -					Not a hazardous
CIL18-MAG					substance or
					mixture
					according to
					Regulation (EC)
					No 1272/2008.

74-MCP-1 Beads - Magnetic - MILLIPLEX					
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information
Millipore - CMCP1-MAG					No hazard pictogram, no signal word, no hazard statement(s), no precautionary statement(s) required.

76-TNF-a Beads - Magnetic - MILLIPLEX					
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Millipore -					No hazard
CTNFA-MAG					pictogram, no
					signal word, no

Millipore- CCYTOMAG-90K Page 13 of 14

A

		hazard
		statement(s), no
		precautionary
		statement(s)
		required.

Millipore- CCYTOMAG-90K Page 14 of 14

The life science business of Merck operates as MilliporeSigma in the US and Canada