

Glass Filtration Apparatus



Features

- ▶ The screwthread connector helps remove the silicone tube quickly.
- ▶ The silicone fixing sucker prevents the receiver flask from tipping over accidentally.


Applications

Solvent purification and general filtration such as microbiology, suspended solid analysis etc.

Specifications

Material	VF 7
Funnel	Borosilicate glass
Support base	Borosilicate glass
Membrane support	SS screen with PTFE gasket
Clamp	Aluminum
Stopper	Silicone rubber
Receiver flask	Borosilicate glass
Specifications	
Funnel capacity	300 ml
Flask capacity	1000 ml
Filter diameter	47 mm
Hose barb	ID8

*The flask withstands nearly full vacuum for use in vacuum filtration.



VF 7

Ordering Information

VF 7, Glass Filtration Set		167200-07
300 ml Glass filter holder		167120-33
1	Glass funnel, 300 ml	167210-03
2	Aluminum clamp	167240-01
3	SS support screen	167230-31
4	PTFE gasket	167230-41
5	Glass support base (Including No.8 silicone stopper)	167230-07
6	Receiver flask, 1000 ml (with screwthread connector)	167250-03
7	Silicone fixing sucker	167200-45

>> **11,000 THB** <<
(Exclude VAT)