



PCR 96 attachment

/// Data Sheet

The Matrix PCR 96 attachment is suitbale for 96-well-PCR plates. The ideal attachment design ensures an optimal temperature transfer and thus guarantees valid results. The thermoshakers Matrix Orbital Delta plus and Matrix Delta plus automatically recognize, which attachment is connected to the device via an interface in the locking mechanism. The attachment can be changed quickly and easily with an ergonomic lever.

Speed max.: 2000 rpm Temperature max.: 100 °C



Suitable for: Matrix Orbital Delta plus Matrix Delta plus



Technical Data Dimensions (W x H x D) [mm] 141 x 42 x 150 Weight [kg] 0.56