SAFETY DATA SHEET

Millipore

according to Regulation (EC) No. 1907/2006

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1	Product identifiers Product name	:	Human Cyto/Chemo MAG Premix 20 Plex Kit
	Product Number Catalogue No. Brand REACH No.	:	HCY4MG-64K-PX20 6F1256 Millipore This product is a mixture. REACH Registration Number see section 3.
1.2	Relevant identified u	ses	of the substance or mixture and uses advised against

Identified uses : Biochemical research/analysis

1.3 Details of the supplier of the safety data sheet

:

Company

1.4 Emergency telephone

Emergency Phone #	:	+(44)-870-8200418 (CHEMTREC (GB))
		+(353)-19014670 (CHEMTREC Ireland)
		001-803-017-9114 (CHEMTREC India)

This is a summary SDS for a kit. For the full SDS for each of the components listed in section 16 please visit our website.

SECTION 2: Hazards identification

2.1 Classification of the substance or mixture

Classification according to Regulation (EC) No 1272/2008 Acute toxicity, Oral (Category 4), H302 Acute toxicity, Inhalation (Category 4), H332 Acute toxicity, Dermal (Category 4), H312 Serious eye damage (Category 1), H318 Skin sensitization (Category 1), H317 Specific target organ toxicity - repeated exposure (Category 2), Respiratory Tract, H373 Long-term (chronic) aquatic hazard (Category 2), H411

For the full text of the H-Statements mentioned in this Section, see Section 16.

2.2 Label elements

Labelling according Regulation (EC) No 1272/2008

Millipore- HCY4MG-64K-PX20

The life science business of Merck operates as MilliporeSigma in the US and Canada



Page 1 of 11

Pictogram	
Signal Word	Danger
Hazard statement(s) H302 + H312 + H332 H317 H318 H373	Harmful if swallowed, in contact with skin or if inhaled. May cause an allergic skin reaction. Causes serious eye damage. May cause damage to organs (Respiratory Tract) through prolonged or repeated exposure.
H411	Toxic to aquatic life with long lasting effects.
Precautionary statement(s) P273 P280 P302 + P352 + P312 P304 + P340 + P312 P305 + P351 + P338 P314	Avoid release to the environment. Wear protective gloves/ protective clothing/ eye protection/ face protection. IF ON SKIN: Wash with plenty of water.Call a POISON CENTER/ doctor if you feel unwell. IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a POISON CENTER/ doctor if you feel unwell. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Get medical advice/ attention if you feel unwell.
-	
Supplemental Hazard Statements	none
Reduced Labeling (<= 1 Pictogram	

Signal Word	Danger
Hazard statement(s) H317 H318	May cause an allergic skin reaction. Causes serious eye damage.
Precautionary statement(s) P305 + P351 + P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
Supplemental Hazard Statements	none
Other hazards - none	

SECTION 3: Composition/information on ingredients

Refer to component SDS

Millipore- HCY4MG-64K-PX20

2.3

The life science business of Merck operates as MilliporeSigma in the US and Canada

Page 2 of 11

SECTION 4: First aid measures

Refer to component SDS

SECTION 5: Firefighting measures

Refer to component SDS

SECTION 6: Accidental release measures

Refer to component SDS

SECTION 7: Handling and storage

Precautions for safe handling 7.1

Advice on safe handling

Work under hood. Do not inhale substance/mixture. Avoid generation of vapours/aerosols.

Hygiene measures

Immediately change contaminated clothing. Apply preventive skin protection. Wash hands and face after working with substance. For precautions see section 2.2.

7.2 Conditions for safe storage, including any incompatibilities

Storage conditions

Tightly closed.

Recommended storage temperature see product label.

Storage class

Storage class (TRGS 510): 10: Combustible liquids

7.3 Specific end use(s)

Apart from the uses mentioned in section 1.2 no other specific uses are stipulated

SECTION 8: Exposure controls/personal protection

Refer to component SDS

SECTION 9: Physical and chemical properties

Refer to component SDS

SECTION 10: Stability and reactivity

Refer to component SDS

Millipore- HCY4MG-64K-PX20

The life science business of Merck operates as MilliporeSigma in the US and Canada

Page 3 of 11



SECTION 11: Toxicological information

Refer to component SDS

SECTION 12: Ecological information

Refer to component SDS

SECTION 13: Disposal considerations

Refer to component SDS

SECTION 14: Transport informat	ion	
14.1 UN number ADR/RID: 3316	IMDG: 3316	IATA: 3316
14.2 UN proper shipping name ADR/RID: CHEMICAL KIT IMDG: CHEMICAL KIT IATA: Chemical kit		
14.3 Transport hazard class(es) ADR/RID: 9	IMDG: 9	IATA: 9
14.4 Packaging group ADR/RID: III	IMDG: III	IATA: III
14.5 Environmental hazards ADR/RID: no	IMDG Marine pollutant: no	IATA: no
14.6 Special precautions for use	r	
Further information :	No data available	

SECTION 15: Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

This material safety data sheet complies with the requirements of Regulation (EC) No. 1907/2006.

Refer to component SDS

15.2 Chemical Safety Assessment

For this product a chemical safety assessment was not carried out

SECTION 16: Other information

Kit Components: 10X Wash Buffer -

Millipore

L-WB

Millipore- HCY4MG-64K-PX20

The life science business of Merck operates as MilliporeSigma in the US and Canada



Page 4 of 11

Skin Sens. 1; Aquatic Chronic 2;

MILLIPLEX Streptavidin- Phycoerythrin	Millipore	L-SAPE10	H317, H411 Eye Irrit. 2; H319
Milliplex Assay Buffer	Millipore	L-AB	
Human Cytokine/Chemokine Panel IV Standard	Millipore	MXH8064	Acute Tox. 4; Eye Dam. 1; STOT RE 2; Aquatic Chronic 3; H302 + H312 + H332, H318, H373, H412
Human Cytokine/Chemokine Panel IV Quality Control 1 & 2	Millipore	MXH6064	Acute Tox. 4; Eye Dam. 1; STOT RE 2; Aquatic Chronic 3; H302 + H312 + H332, H318, H373, H412
Human Cytokine/Chemokine Panel IV Detection Antibodies	Millipore	MXH1064	Eye Irrit. 2; H319
Serum Matrix HUMAN CYTOKINE IV 20-PLEX PREMIXED BEADS	Millipore Sigma-Aldrich	MXHSM-12 HCP4PX20-MAG	Aquatic Chronic 3; H412

The life science business of Merck operates as MilliporeSigma in the US and Canada

Page 5 of 11



Full text of other abbreviations

ADN - European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways; ADR - Agreement concerning the International Carriage of Dangerous Goods by Road; AIIC - Australian Inventory of Industrial Chemicals; ASTM -American Society for the Testing of Materials; bw - Body weight; CMR - Carcinogen, Mutagen or Reproductive Toxicant; DIN - Standard of the German Institute for Standardisation; DSL - Domestic Substances List (Canada); ECx - Concentration associated with x% response; ELx - Loading rate associated with x% response; EmS -Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx -Concentration associated with x% growth rate response; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan); ISO - International Organisation for Standardization; KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; n.o.s. -Not Otherwise Specified; NO(A)EC - No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NZIOC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention; PBT - Persistent, Bioaccumulative and Toxic substance; PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; RID - Regulations concerning the International Carriage of Dangerous Goods by Rail; SADT - Self-Accelerating Decomposition Temperature; SDS -Safety Data Sheet; TCSI - Taiwan Chemical Substance Inventory; TECI - Thailand Existing Chemicals Inventory; TSCA - Toxic Substances Control Act (United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative

Classification of the	Classification procedure:	
Acute Tox.4	H302	Calculation method
Acute Tox.4	H332	Calculation method
Acute Tox.4	H312	Calculation method
Eye Dam.1	H318	Calculation method
Skin Sens.1	H317	Calculation method
STOT RE2	H373	Calculation method
Aquatic Chronic2	H411	Calculation method

Further information

The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. The information in this document is based on the present state of our knowledge and is applicable to the product with regard to

Millipore- HCY4MG-64K-PX20

The life science business of Merck operates as MilliporeSigma in the US and Canada



Page 6 of 11

appropriate safety precautions. It does not represent any guarantee of the properties of the product. Sigma-Aldrich Corporation and its Affiliates shall not be held liable for any damage resulting from handling or from contact with the above product. See www.sigma-aldrich.com and/or the reverse side of invoice or packing slip for additional terms and conditions of sale.

Copyright 2020 Sigma-Aldrich Co. LLC. License granted to make unlimited paper copies for internal use only.

The branding on the header and/or footer of this document may temporarily not visually match the product purchased as we transition our branding. However, all of the information in the document regarding the product remains unchanged and matches the product ordered. For further information please contact mlsbranding@sial.com.

Millipore- HCY4MG-64K-PX20

The life science business of Merck operates as MilliporeSigma in the US and Canada

Page 7 of 11



Kit Components:

10X Wash Buf	10X Wash Buffer - MILLIPLEX							
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental			
identifier		Word	Statements	Statements	information			
Millipore - L- WB		Warning	May cause an allergic skin reaction. Toxic to aquatic life with long lasting effects.	Avoid breathing mist or vapors. IF ON SKIN: Wash with plenty of water. Contaminated work clothing should not be allowed out of the workplace. If skin irritation or rash occurs: Get medical advice/ attention. Avoid release to the environment. Wear protective gloves.				

Streptavidin-Phycoerythrin						
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information	
Millipore - L- SAPE10	!>	Warning	Causes serious eye irritation.	Wash skin thoroughly after handling. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Wear eye protection/ face protection. If eye irritation persists: Get		

Millipore- HCY4MG-64K-PX20

The life science business of Merck operates as MilliporeSigma in the US and Canada




Page 8 of 11

		medical advice/	
		attention.	

Milliplex Assay Buffer					
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information
Millipore - L- AB					No hazard pictogram, no signal word, no hazard statement(s), no precautionary statement(s) required

Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier	5	Word	Statements	Statements	information
Millipore - MXH8064		Danger	Harmful if swallowed, in contact with skin or if inhaled. Causes serious eye damage. May cause damage to organs (Respiratory Tract) through prolonged or repeated exposure. Harmful to aquatic life with long lasting effects.	Avoid release to the environment. IF ON SKIN: Wash with plenty of water.Call a POISON CENTER/ doctor if you feel unwell. Wear protective gloves/ protective clothing/ eye protection/ face protection. IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a POISON CENTER/ doctor if you feel unwell. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and	

The life science business of Merck operates as MilliporeSigma in the US and Canada



medical a	Get dvice/
attention	if you
feel unwe	II.

			V Quality Control 1		
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Millipore -		Danger	Harmful if	Avoid release to	
MXH6064			swallowed, in	the	
			contact with skin	environment. IF	
	* *		or if inhaled.	ON SKIN: Wash	
	Ă I		Causes serious	with plenty of	
			eye damage.	water.Call a	
			May cause	POISON	
			damage to	CENTER/ doctor	
			organs	if you feel	
			(Respiratory	unwell. Wear	
			Tract) through	protective	
			prolonged or	gloves/	
			repeated	protective	
			exposure.	clothing/ eye	
			Harmful to	protection/ face	
			aquatic life with	protection. IF	
			long lasting	INHALED:	
			effects.	Remove person	
				to fresh air and	
				keep	
				comfortable for	
				breathing. Call a	
				POISON	
				CENTER/ doctor	
				if you feel	
				unwell. IF IN	
				EYES: Rinse	
				cautiously with	
				water for	
				several minutes.	
				Remove contact	
				lenses, if	
				present and	
				easy to do.	
				Continue	
				rinsing. Get	
				medical advice/	
				attention if you	
				feel unwell.	1

The life science business of Merck operates as MilliporeSigma in the US and Canada

Page $10~{\rm of}~11$



Human Cytokine/Chemokine Panel IV Detection Antibodies					
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Millipore -		Warning	Causes serious	Wash skin	
MXH1064			eye irritation.	thoroughly after	
	•			handling. IF IN	
				EYES: Rinse	
				cautiously with	
				water for	
				several minutes.	
				Remove contact	
				lenses, if	
				present and	
				easy to do.	
				Continue	
				rinsing. Wear	
				eye protection/	
				face protection.	
				If eye irritation	
				persists: Get	
				medical advice/	
				attention.	

Serum Matrix					
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information
Millipore - MXHSM-12			Harmful to aquatic life with long lasting effects.	Avoid release to the environment. Dispose of contents/ container to an approved waste disposal plant.	

HUMAN CYTOKINE IV 20-PLEX PREMIXED BEADS					
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information
Sigma-Aldrich - HCP4PX20- MAG					Not a hazardous substance or mixture according to Regulation (EC) No 1272/2008.

The life science business of Merck operates as MilliporeSigma in the US and Canada

Page 11 of 11

