

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Version 8.6 Revision Date 11.04.2025 Print Date 21.04.2025

GENERIC EU MSDS - NO COUNTRY SPECIFIC DATA - NO OEL DATA

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1 Product identifiers

Product name : Human Cytokine/Chemokine MAG Panel 3

Product Number : HCYP3MAG-63K

Catalogue No. : 667603 Brand : Millipore

REACH No. :

1.2 Relevant identified uses of the substance or mixture and uses advised against

Identified uses : Biochemical research/analysis

1.3 Details of the supplier of the safety data sheet

Company :

1.4 Emergency telephone

Emergency Phone #: +(44)-870-8200418 (CHEMTREC (GB))

+(353)-19014670 (CHEMTREC Ireland) 001-803-017-9114 (CHEMTREC India)

This is a summary SDS for a kit. For the full SDS for each of the components listed in section 16 please visit our website.

SECTION 2: Hazards identification

2.1 Classification of the substance or mixture

Acute toxicity, (Category 4) H302: Harmful if swallowed.

Acute toxicity, (Category 4) H332: Harmful if inhaled.

Acute toxicity, (Category 4) H312: Harmful in contact with skin.

Skin irritation, (Category 2) H315: Causes skin irritation.

Serious eye damage, (Category H318: Causes serious eye damage.

1)

Skin sensitization, (Category 1) H317: May cause an allergic skin reaction.

Specific target organ toxicity - H373: May cause damage to organs repeated exposure, (Category 2), through prolonged or repeated exposure.

Respiratory Tract

Millipore- HCYP3MAG-63K Page 1 of 14



Long-term (chronic) aquatic hazard, (Category 2)

H411: Toxic to aquatic life with long lasting

effects.

2.2 Label elements

Labelling according Regulation (EC) No 1272/2008

Pictogram

Signal Word Danger

Hazard Statements

H302 + H312 + H332 Harmful if swallowed, in contact with skin or if inhaled.

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H318 Causes serious eye damage.

H373 May cause damage to organs (Respiratory Tract) through

prolonged or repeated exposure.

H411 Toxic to aquatic life with long lasting effects.

Precautionary Statements

P273 Avoid release to the environment.

P280 Wear protective gloves/ protective clothing/ eye protection/ face

protection.

P302 + P352 + P312 IF ON SKIN: Wash with plenty of water.Call a POISON CENTER/

doctor if you feel unwell.

P304 + P340 + P312 IF INHALED: Remove person to fresh air and keep comfortable

for breathing. Call a POISON CENTER/ doctor if you feel unwell.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue

rinsing.

P314 Get medical advice/ attention if you feel unwell.

Supplemental Hazard

Statements

none

Reduced Labeling (<= 125 ml)

Pictogram



Signal Word Danger

Hazard Statements

H317 May cause an allergic skin reaction.

H318 Causes serious eye damage.

Precautionary Statements

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue

rinsing.

Supplemental Hazard

Statements

none

2.3 Other hazards

This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher.

Millipore- HCYP3MAG-63K Page 2 of 14



Ecological information:

The substance/mixture does not contain components considered to have endocrine disrupting properties according to REACH Article 57(f) or Commission Delegated regulation (EU) 2017/2100 or Commission Regulation (EU) 2018/605 at levels of 0.1% or higher. Toxicological information:

The substance/mixture does not contain components considered to have endocrine disrupting properties according to REACH Article 57(f) or Commission Delegated regulation (EU) 2017/2100 or Commission Regulation (EU) 2018/605 at levels of 0.1% or higher.

SECTION 3: Composition/information on ingredients

Refer to component SDS

SECTION 4: First aid measures

Refer to component SDS

SECTION 5: Firefighting measures

Refer to component SDS

SECTION 6: Accidental release measures

Refer to component SDS

SECTION 7: Handling and storage

7.1 Precautions for safe handling

Advice on safe handling

Work under hood. Do not inhale substance/mixture. Avoid generation of vapours/aerosols.

Hygiene measures

Immediately change contaminated clothing. Apply preventive skin protection. Wash hands and face after working with substance.

For precautions see section 2.2.

7.2 Conditions for safe storage, including any incompatibilities

Storage conditions

Tightly closed.

Recommended storage temperature see product label.

Storage class

Storage class (TRGS 510): 10: Combustible liquids

7.3 Specific end use(s)

Apart from the uses mentioned in section 1.2 no other specific uses are stipulated

Millipore- HCYP3MAG-63K Page 3 of 14



SECTION 8: Exposure controls/personal protection

Refer to component SDS

SECTION 9: Physical and chemical properties

Refer to component SDS

SECTION 10: Stability and reactivity

Refer to component SDS

SECTION 11: Toxicological information

Refer to component SDS

SECTION 12: Ecological information

Refer to component SDS

SECTION 13: Disposal considerations

Refer to component SDS

SECTION 14: Transport information

14.1 UN number

ADR/RID: - IMDG: - IATA: -

14.2 UN proper shipping name

ADR/RID: Not dangerous goods IMDG: Not dangerous goods IATA: Not dangerous goods

14.3 Transport hazard class(es)

ADR/RID: - IMDG: - IATA: -

14.4 Packaging group

ADR/RID: - IMDG: - IATA: -

14.5 Environmental hazards

ADR/RID: no IMDG Marine pollutant: no IATA: no

14.6 Special precautions for user

Further information

Not classified as dangerous in the meaning of transport regulations.

Millipore- HCYP3MAG-63K Page 4 of 14



14.7 Maritime transport in bulk according to IMO instruments

Not applicable for product as supplied.

SECTION 15: Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

This material safety data sheet complies with the requirements of Regulation (EC) No. 1907/2006.

Refer to component SDS

15.2 Chemical Safety Assessment

For this product a chemical safety assessment was not carried out

SECTIO	N 16:	Other	inform	ation

Kit Components: Human Cytokine Panel III standard	Millipore	MXH8063	Acute Tox. 4; Skin Irrit. 2; Eye Dam. 1; STOT RE 2; Aquatic Chronic 3; H302 + H312 + H332, H315, H318, H373, H412
Human Cytokine Panel III Quality Control	Millipore	MXH6063	Acute Tox. 4; Skin Irrit. 2; Eye Dam. 1; STOT RE 2; Aquatic Chronic 3; H302 + H312 + H332, H315, H318, H373, H412
Serum Matrix - MILLIPLEX (w/ antisera)	Millipore	MXHSM-3	Acute Tox. 4; Aquatic Chronic 3; H312, H412
Milliplex Assay Buffer	Millipore	L-AB	11312, 11112
	•		
Bead Diluent COLD PACKS	Millipore	LBD	
10X Wash Buffer - MILLIPLEX	Millipore	L-WB	Skin Sens. 1; Aquatic Chronic 2; H317, H411
Human Cytokine Panel III Detection Abs	Millipore	MXH1063-1	Eye Irrit. 2; H319
Human Cytokine Panel III Detection Abs	Millipore	MXH1063-2	Eye Irrit. 2; H319
Streptavidin- Phycoerythrin (3.5 ml) MILLI	Millipore	L-SAPE6	Eye Irrit. 2; H319
51-M-CSF Beads - Magnetic MILLIPLEX	Millipore	HMCSF-MAG	
47-CXCL9/MIG Beads - Magnetic	Millipore	HMIG-MAG	
13-NAP-2/CXCL7 Beads - Magnetic	Millipore	HNAP2-MAG	
15-GCP-2/CXCL6 Beads - Magnetic	Millipore	HGCP2-MAG	
19-I-TAC/CXCL11 Beads - Magnetic	Millipore	HITAC-MAG	
21-HCC-1/CCL14a Beads - Magnetic	Millipore	HHCC1-MAG	
26-MIP3B/CCL19 Beads - Magnetic	Millipore	HMIP3B-MAG	
28-MIP3a/CCL20 Beads -	Sigma-Aldrich	HMIP3A-MAG	

Millipore- HCYP3MAG-63K Page 5 of 14



Magnetic

30-Lymphotactin/XCL1 Millipore HLTCTN-MAG

Beads-Magnetic

MILLIPLEX

34-IL-11 Beads - Millipore HIL11-MAG

Magnetic MILLIPLEX

36-IL-29/IFN-Lambda1 Millipore HIL29-MAG

Beads-Magnetic MILLIPLEX

Full text of other abbreviations

ADN - European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways; ADR - Agreement concerning the International Carriage of Dangerous Goods by Road; AIIC - Australian Inventory of Industrial Chemicals; ASTM -American Society for the Testing of Materials; bw - Body weight; CMR - Carcinogen, Mutagen or Reproductive Toxicant; DIN - Standard of the German Institute for Standardisation; DSL - Domestic Substances List (Canada); ECx - Concentration associated with x% response; ELx - Loading rate associated with x% response; EmS -Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx -Concentration associated with x% growth rate response; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan); ISO - International Organisation for Standardization; KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; n.o.s. -Not Otherwise Specified; NO(A)EC - No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention; PBT - Persistent, Bioaccumulative and Toxic substance; PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; RID - Regulations concerning the International Carriage of Dangerous Goods by Rail; SADT - Self-Accelerating Decomposition Temperature; SDS -Safety Data Sheet; TCSI - Taiwan Chemical Substance Inventory; TECI - Thailand Existing Chemicals Inventory; TSCA - Toxic Substances Control Act (United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative

Classification of the mixture		Classification procedure:
Acute Tox.4	H302	Calculation method
Acute Tox.4	H332	Calculation method
Acute Tox.4	H312	Calculation method
Skin Irrit.2	H315	Calculation method
Eye Dam.1	H318	Calculation method
Skin Sens.1	H317	Calculation method
STOT RE2	H373	Calculation method

Millipore- HCYP3MAG-63K Page 6 of 14



Further information

The information is believed to be correct but is not exhaustive and will be used solely as a guideline, which is based on current knowledge of the chemical substance or mixture and is applicable to appropriate safety precautions for the product. It does not represent any guarantee of the properties of the product. Sigma-Aldrich Corporation and its Affiliates shall not be held liable for any damage resulting from handling or from contact with the above product. See www.sigma-aldrich.com and/or the reverse side of invoice or packing slip for additional terms and conditions of sale.

Copyright 2020 Sigma-Aldrich Co. LLC. License granted to make unlimited paper copies for internal use only.

The branding on the header and/or footer of this document may temporarily not visually match the product purchased as we transition our branding. However, all of the information in the document regarding the product remains unchanged and matches the product ordered. For further information please contact mlsbranding@sial.com.

Millipore- HCYP3MAG-63K Page 7 of 14





Kit Components:

Human Cyto	kine Panel II	[standard			
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Millipore -		Danger	Harmful if	Avoid release to	
MXH8063			swallowed, in	the	
	X		contact with skin	environment. IF	
	(T. S.)		or if inhaled.	ON SKIN: Wash	
	X		Causes skin	with plenty of	
	! >		irritation.	water.Call a	
			Causes serious	POISON	
			eye damage.	CENTER/ doctor	
			May cause	if you feel	
			damage to	unwell. Wear	
			organs	protective	
			(Respiratory	gloves/	
			Tract) through	protective	
			prolonged or	clothing/ eye	
			repeated	protection/ face	
			exposure.	protection. IF	
			Harmful to	INHALED:	
			aquatic life with	Remove person	
			long lasting effects.	to fresh air and	
			errects.	keep comfortable for	
				breathing. Call a	
				POISON	
				CENTER/ doctor	
				if you feel	
				unwell. IF IN	
				EYES: Rinse	
				cautiously with	
				water for	
				several minutes.	
				Remove contact	
				lenses, if	
				present and	
				easy to do.	
				Continue	
				rinsing. Get	
				medical advice/	
				attention if you	
				feel unwell.	

Human Cytokine Panel III Quality Control						
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information	
Millipore -		Danger	Harmful if swallowed, in	Avoid release to the		

Millipore- HCYP3MAG-63K Page 8 of 14



MANUELLE	_			
MXH6063		contact with skin	environment. IF	
		or if inhaled.	ON SKIN: Wash	
		Causes skin	with plenty of	
		irritation.	water.Call a	
	~	Causes serious	POISON	
		eye damage.	CENTER/ doctor	
		May cause	if you feel	
		damage to	unwell. Wear	
		organs	protective	
		(Respiratory	gloves/	
		Tract) through	protective	
		prolonged or	clothing/ eye	
		repeated	protection/ face	
		exposure.	protection. IF	
		Harmful to	INHALED:	
		aquatic life with	Remove person	
		long lasting	to fresh air and	
		effects.	keep	
		Circus.	comfortable for	
			breathing. Call a	
			POISON	
			CENTER/ doctor	
			if you feel	
			unwell. IF IN	
			EYES: Rinse	
			cautiously with	
			water for	
			several minutes.	
			Remove contact	
			lenses, if	
			present and	
			easy to do.	
			Continue	
			rinsing. Get	
			medical advice/	
			attention if you	
			feel unwell.	

Serum Matrix - MILLIPLEX (w/ antisera)						
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information	
Millipore - MXHSM-3	•	Warning	Harmful in contact with skin. Harmful to aquatic life with long lasting effects.	Avoid release to the environment. IF ON SKIN: Wash with plenty of water.Call a POISON CENTER/ doctor if you feel unwell. Dispose of contents/ container to an approved waste disposal plant.		

Millipore- HCYP3MAG-63K Page 9 of 14



	Wear protective
	gloves/
	protective
	clothing. Take
	off contaminated
	clothing and
	wash it before
	reuse.

Milliplex Assay Buffer						
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental	
identifier		Word	Statements	Statements	information	
Millipore - L-					No hazard	
AB					pictogram, no	
					signal word, no	
					hazard	
					statement(s), no	
					precautionary	
1					statement(s)	
					required.	

Bead Diluent	Bead Diluent COLD PACKS					
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information	
Millipore - LBD			State.mento	o careionico	No hazard pictogram, no signal word, no hazard statement(s), no precautionary statement(s) required.	

10Y Wash Ru	10X Wash Buffer - MILLIPLEX							
Product identifier Millipore - L- WB	Pictogram	Signal Word Warning	Hazard Statements May cause an allergic skin reaction. Toxic to aquatic life with long	Precautionary Statements Avoid breathing mist or vapors. IF ON SKIN: Wash with plenty of water.	Supplemental information			
			lasting effects.	Contaminated work clothing should not be allowed out of the workplace. If skin irritation or rash occurs: Get medical				
				advice/ attention. Avoid release to the environment. Wear protective				

Millipore- HCYP3MAG-63K Page 10 of 14



	LUIUVES.	

Human Cytol	Human Cytokine Panel III Detection Abs						
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental		
identifier		Word	Statements	Statements	information		
Millipore -		Warning	Causes serious	Wash skin			
MXH1063-1			eye irritation.	thoroughly after			
	•			handling. IF IN			
				EYES: Rinse			
				cautiously with			
				water for			
				several minutes.			
				Remove contact			
				lenses, if			
				present and			
				easy to do.			
				Continue			
				rinsing. Wear			
				eye protection/			
				face protection.			
				If eye irritation			
				persists: Get			
				medical advice/			
				attention.			

Human Cytokine Panel III Detection Abs							
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental		
identifier		Word	Statements	Statements	information		
Millipore -		Warning	Causes serious	Wash skin			
MXH1063-2			eye irritation.	thoroughly after			
	·			handling. IF IN			
				EYES: Rinse			
				cautiously with			
				water for			
				several minutes.			
				Remove contact			
				lenses, if			
				present and			
				easy to do.			
				Continue			
				rinsing. Wear			
				eye protection/			
				face protection.			
				If eye irritation			
				persists: Get			
				medical advice/			
				attention.			

Streptavidin-Phycoerythrin (3.5 ml) MILLI						
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental	
identifier		Word	Statements	Statements	information	
Millipore - L- SAPE6	1	Warning	Causes serious eye irritation.	Wash skin thoroughly after handling. IF IN EYES: Rinse		

Millipore- HCYP3MAG-63K Page 11 of 14



	cautiously with	
	water for	
	several minutes.	
	Remove contact	
	lenses, if	
	present and	
	easy to do.	
	Continue	
	rinsing. Wear	
	eye protection/	
	face protection.	
	If eye irritation	
	persists: Get	
	medical advice/	
	attention.	

51-M-CSF Bea	51-M-CSF Beads - Magnetic MILLIPLEX						
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental		
identifier		Word	Statements	Statements	information		
Millipore -					No hazard		
HMCSF-MAG					pictogram, no		
					signal word, no		
					hazard		
					statement(s), no		
					precautionary		
					statement(s)		
					required.		

47-CXCL9/M	47-CXCL9/MIG Beads - Magnetic						
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information		
Millipore - HMIG-MAG					No hazard pictogram, no signal word, no hazard statement(s), no		
					precautionary statement(s) required.		

13-NAP-2/CXCL7 Beads - Magnetic							
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental		
identifier		Word	Statements	Statements	information		
Millipore -					No hazard		
HNAP2-MAG					pictogram, no		
					signal word, no		
					hazard		
					statement(s), no		
					precautionary		
					statement(s)		
					required.		

Millipore- HCYP3MAG-63K Page 12 of 14



15-GCP-2/CX	15-GCP-2/CXCL6 Beads - Magnetic							
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental			
identifier		Word	Statements	Statements	information			
Millipore -					No hazard			
HGCP2-MAG					pictogram, no			
					signal word, no			
					hazard			
					statement(s), no			
					precautionary			
					statement(s)			
					required.			

19-I-TAC/CX	19-I-TAC/CXCL11 Beads - Magnetic							
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental			
identifier		Word	Statements	Statements	information			
Millipore -					No hazard			
HITAC-MAG					pictogram, no			
					signal word, no			
					hazard			
					statement(s), no			
					precautionary			
					statement(s)			
					required.			

21-HCC-1/CC	21-HCC-1/CCL14a Beads - Magnetic						
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental		
identifier		Word	Statements	Statements	information		
Millipore -					No hazard		
HHCC1-MAG					pictogram, no		
					signal word, no		
					hazard		
					statement(s), no		
					precautionary		
					statement(s)		
					required.		

26-MIP3B/C0	26-MIP3B/CCL19 Beads - Magnetic						
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental		
identifier		Word	Statements	Statements	information		
Millipore - HMIP3B-MAG					No hazard pictogram, no		
THE SE PIAC					signal word, no hazard statement(s), no precautionary		
					statement(s) required.		
					required.		

28-MIP3a/CCL20 Beads - Magnetic							
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information		
Sigma-Aldrich - HMIP3A- MAG					No hazard pictogram, no signal word, no hazard		

Millipore- HCYP3MAG-63K Page 13 of 14



		statement(s), no
		precautionary
		statement(s)
		required.

30-Lymphotactin/XCL1 Beads-Magnetic MILLIPLEX					
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Millipore -					No hazard
HLTCTN-MAG					pictogram, no
					signal word, no
					hazard
					statement(s), no
					precautionary
					statement(s)
					required.

34-IL-11 Beads - Magnetic MILLIPLEX					
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Millipore -					Not a hazardous
HIL11-MAG					substance or
					mixture
					according to
					Regulation (EC)
					No 1272/2008.

36-IL-29/IFN-Lambda1 Beads-Magnetic MILLIPLEX					
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information
Millipore - HIL29-MAG					No hazard pictogram, no signal word, no hazard statement(s), no precautionary statement(s) required.

Millipore- HCYP3MAG-63K Page 14 of 14

