

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Version 8.4 Revision Date 18.10.2024 Print Date 23.04.2025

GENERIC EU MSDS - NO COUNTRY SPECIFIC DATA - NO OEL DATA

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1 Product identifiers

Product name : HUMAN CIRCULATING CANCER BIOMARKER

MAGNETIC BEAD PANEL 3

Product Number : HCCBP3MAG-58K

Catalogue No. : 667588 Brand : Millipore

REACH No. : This product is a mixture. REACH Registration Number see

section 3.

1.2 Relevant identified uses of the substance or mixture and uses advised against

Identified uses : Biochemical research/analysis

1.3 Details of the supplier of the safety data sheet

Company :

1.4 Emergency telephone

Emergency Phone # : +(44)-870-8200418 (CHEMTREC (GB))

+(353)-19014670 (CHEMTREC Ireland) 001-803-017-9114 (CHEMTREC India)

This is a summary SDS for a kit. For the full SDS for each of the components listed in section 16 please visit our website.

SECTION 2: Hazards identification

2.1 Classification of the substance or mixture

Acute toxicity, (Category 4) H302: Harmful if swallowed.

Acute toxicity, (Category 4) H332: Harmful if inhaled.

Acute toxicity, (Category 3) H311: Toxic in contact with skin.

Eye irritation, (Category 2) H319: Causes serious eye irritation.

Skin sensitization, (Category 1) H317: May cause an allergic skin reaction.

Long-term (chronic) aquatic H411: Toxic to aquatic life with long lasting

hazard, (Category 2) effects.

Millipore- HCCBP3MAG-58K Page 1 of 12



2.2 Label elements

Labelling according Regulation (EC) No 1272/2008

Pictogram

Signal Word Danger

Hazard Statements

H302 + H332 Harmful if swallowed or if inhaled.

H311 Toxic in contact with skin.

H317 May cause an allergic skin reaction. H319 Causes serious eye irritation.

H411 Toxic to aquatic life with long lasting effects.

Precautionary Statements

P273 Avoid release to the environment.

P280 Wear protective gloves/ protective clothing/ eye protection/ face

protection.

P301 + P312 IF SWALLOWED: Call a POISON CENTER/ doctor if you feel

unwell.

P302 + P352 + P312 IF ON SKIN: Wash with plenty of water. Call a POISON CENTER/

doctor if you feel unwell.

P304 + P340 + P312 IF INHALED: Remove person to fresh air and keep comfortable

for breathing. Call a POISON CENTER/ doctor if you feel unwell.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue

rinsing.

Supplemental Hazard

Statements

none

Reduced Labeling (<= 125 ml)

Pictogram

Signal Word Danger

Hazard Statements

H311 Toxic in contact with skin.

H317 May cause an allergic skin reaction.

Precautionary Statements

P302 + P352 + P312 IF ON SKIN: Wash with plenty of water.Call a POISON CENTER/

doctor if you feel unwell.

Supplemental Hazard

Statements

none

2.3 Other hazards - none

SECTION 3: Composition/information on ingredients

Refer to component SDS

SECTION 4: First aid measures

Refer to component SDS

Millipore- HCCBP3MAG-58K Page 2 of 12

A

SECTION 5: Firefighting measures

Refer to component SDS

SECTION 6: Accidental release measures

Refer to component SDS

SECTION 7: Handling and storage

7.1 Precautions for safe handling

Advice on safe handling

Work under hood. Do not inhale substance/mixture. Avoid generation of vapours/aerosols.

Hygiene measures

Immediately change contaminated clothing. Apply preventive skin protection. Wash hands and face after working with substance.

For precautions see section 2.2.

7.2 Conditions for safe storage, including any incompatibilities

Storage conditions

Tightly closed. Keep in a well-ventilated place. Keep locked up or in an area accessible only to qualified or authorized persons.

Recommended storage temperature see product label.

Storage class

Storage class (TRGS 510): 6.1C: Combustible, acute toxic Cat.3 / toxic compounds or compounds which causing chronic effects

7.3 Specific end use(s)

Apart from the uses mentioned in section 1.2 no other specific uses are stipulated

SECTION 8: Exposure controls/personal protection

Refer to component SDS

SECTION 9: Physical and chemical properties

Refer to component SDS

SECTION 10: Stability and reactivity

Refer to component SDS

SECTION 11: Toxicological information

Refer to component SDS

Millipore- HCCBP3MAG-58K Page 3 of 12



SECTION 12: Ecological information

Refer to component SDS

SECTION 13: Disposal considerations

Refer to component SDS

SECTION 14: Transport information

14.1 UN number

ADR/RID: - IMDG: - IATA: -

14.2 UN proper shipping name

ADR/RID: Not dangerous goods IMDG: Not dangerous goods IATA: Not dangerous goods

14.3 Transport hazard class(es)

ADR/RID: - IMDG: - IATA: -

14.4 Packaging group

ADR/RID: - IMDG: - IATA: -

14.5 Environmental hazards

ADR/RID: no IMDG Marine pollutant: no IATA: no

14.6 Special precautions for user

Further information

Not classified as dangerous in the meaning of transport regulations.

SECTION 15: Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

This material safety data sheet complies with the requirements of Regulation (EC) No. 1907/2006.

Refer to component SDS

15.2 Chemical Safety Assessment

For this product a chemical safety assessment was not carried out

SECTION 16: Other information

Kit Components:

Human Circulating Millipore HCC3-8058-3 Acute Tox. 4; Acute Tox. 3;

Cancer Biomarker Panel Aquatic Chronic 3; H302 +

3 Standard H332, H311, H412

Human Circulating Millipore HCC3-6058-3 Acute Tox. 4; Acute Tox. 3;

Millipore- HCCBP3MAG-58K Page 4 of 12



Cancer Biomarker Panel 3 Quality Control 1,2.			Aquatic Chronic 3; H302 + H332, H311, H412
Milliplex Assay Buffer Bead Diluent MILLIPLEX	Millipore Millipore	L-AB LBD-4	, ,
10X Wash Buffer -	Millipore	L-WB	Skin Sens. 1; Aquatic Chronic 2;
MILLIPLEX Human Circulating Cancer Biomarker Panel 3 Detection Antibodies	Millipore	HCC3-1058-3	H317, H411
Streptavidin- Phycoerythrin (3.5 ml) MILLI	Millipore	L-SAPE6	Eye Irrit. 2; H319
19-Galectin-3 Beads - Magnetic	Millipore	HGLCTN3-MAG	
22-Cathepsin D Beads - Magnetic	Millipore	HCTHPSND- MAG	
26-Fibroblast Activation Protein a Beads -	Millipore	HFAP-MAG	
Magnetic 34-MIA (Melanoma Inhibitory Activity Protein) Beads - Magnetic	Millipore	HMIA-MAG	
37- Myeloperoxidase Beads - Magnetic	Millipore	HMP0-MAG	
42-HSHBG (Sex Hormone Binding Globulin) Beads -	Millipore	HSHBG-MAG	
Magnetic 72-IGFBP3 Beads - Magnetic	Millipore	HIGFBP3-MAG	
75-Ferritin Beads - Magnetic	Millipore	HFRRTN-MAG	

Millipore- HCCBP3MAG-58K Page 5 of 12



Full text of other abbreviations

ADN - European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways; ADR - Agreement concerning the International Carriage of Dangerous Goods by Road; AIIC - Australian Inventory of Industrial Chemicals; ASTM -American Society for the Testing of Materials; bw - Body weight; CMR - Carcinogen, Mutagen or Reproductive Toxicant; DIN - Standard of the German Institute for Standardisation; DSL - Domestic Substances List (Canada); ECx - Concentration associated with x% response; ELx - Loading rate associated with x% response; EmS -Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx -Concentration associated with x% growth rate response; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan); ISO - International Organisation for Standardization; KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; n.o.s. -Not Otherwise Specified; NO(A)EC - No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention; PBT - Persistent, Bioaccumulative and Toxic substance; PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; RID - Regulations concerning the International Carriage of Dangerous Goods by Rail; SADT - Self-Accelerating Decomposition Temperature; SDS -Safety Data Sheet; TCSI - Taiwan Chemical Substance Inventory; TECI - Thailand Existing Chemicals Inventory; TSCA - Toxic Substances Control Act (United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative

Classification of the	Classification procedure:	
Acute Tox.4	H302	Calculation method
Acute Tox.4	H332	Calculation method
Acute Tox.3	H311	Calculation method
Eye Irrit.2	H319	Calculation method
Skin Sens.1	H317	Calculation method
Aquatic Chronic2	H411	Calculation method

Millipore- HCCBP3MAG-58K Page 6 of 12



Further information

The information is believed to be correct but is not exhaustive and will be used solely as a guideline, which is based on current knowledge of the chemical substance or mixture and is applicable to appropriate safety precautions for the product. It does not represent any guarantee of the properties of the product. Sigma-Aldrich Corporation and its Affiliates shall not be held liable for any damage resulting from handling or from contact with the above product. See www.sigma-aldrich.com and/or the reverse side of invoice or packing slip for additional terms and conditions of sale.

Copyright 2020 Sigma-Aldrich Co. LLC. License granted to make unlimited paper copies for internal use only.

The branding on the header and/or footer of this document may temporarily not visually match the product purchased as we transition our branding. However, all of the information in the document regarding the product remains unchanged and matches the product ordered. For further information please contact mlsbranding@sial.com.

Millipore- HCCBP3MAG-58K Page 7 of 12





Kit Components:

Human Circula	ating Cancer	Biomarker	Panel 3 Standard		
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Product		Signal	Hazard	Precautionary Statements Avoid breathing dust. IF SWALLOWED: Call a POISON CENTER/ doctor if you feel unwell. Avoid release to the environment. IF ON SKIN: Wash with plenty of water.Call a POISON CENTER/ doctor if you feel unwell. Wear protective gloves/ protective clothing. IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a	
				INHALED: Remove person to fresh air and keep comfortable for	

Human Circulating Cancer Biomarker Panel 3 Quality Control 1,2.					
Product	Pictogram	Signal Word	Hazard	Precautionary Statements	Supplemental information
identifier Millipore - HCC3-6058-3		Danger	Statements Harmful if swallowed or if inhaled. Toxic in contact	Statements Avoid breathing dust. IF SWALLOWED: Call a POISON	Information
			with skin. Harmful to aquatic life with long lasting	CENTER/ doctor if you feel unwell. Avoid release to the	
			effects.	environment. IF ON SKIN: Wash with plenty of	

Millipore- HCCBP3MAG-58K Page 8 of 12



water.Call a POISON CENTER/ doctor if you feel unwell. Wear protective gloves/ protective clothing. IF INHALED: Remove person
to fresh air and keep comfortable for breathing. Call a POISON CENTER/ doctor if you feel unwell.

Milliplex Assay Buffer					
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Millipore - L-					No hazard
AB					pictogram, no
					signal word, no
					hazard
					statement(s), no
					precautionary
					statement(s)
					required.

Bead Diluent MILLIPLEX					
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information
Millipore - LBD-4					No hazard pictogram, no signal word, no hazard statement(s), no precautionary statement(s) required.

10X Wash Buffer - MILLIPLEX					
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Millipore - L-		Warning	May cause an	Avoid breathing	
WB			allergic skin	mist or vapors.	
	N.		reaction.	IF ON SKIN:	
	***		Toxic to aquatic	Wash with	
			life with long	plenty of water.	
			lasting effects.	Contaminated	
				work clothing	
				should not be	

Millipore- HCCBP3MAG-58K Page 9 of 12



allowed out of the workplace. If skin irritation or rash occurs: Get medical
advice/ attention. Avoid release to the environment.
Wear protective gloves.

Human Circul	Human Circulating Cancer Biomarker Panel 3 Detection Antibodies				
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Millipore -					Not a hazardous
HCC3-1058-3					substance or
					mixture
					according to
					Regulation (EC)
					No 1272/2008.

Streptavidin-	Streptavidin-Phycoerythrin (3.5 ml) MILLI					
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental	
identifier		Word	Statements	Statements	information	
Millipore - L-		Warning	Causes serious	Wash skin		
SAPE6	•		eye irritation.	thoroughly after		
				handling. IF IN		
				EYES: Rinse		
				cautiously with		
				water for several minutes.		
				Remove contact		
				lenses, if		
				present and		
				easy to do.		
				Continue		
				rinsing. Wear		
				eye protection/		
				face protection.		
				If eye irritation		
				persists: Get		
				medical advice/		
				attention.		

19-Galectin-3 Beads - Magnetic					
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Millipore -					No hazard
HGLCTN3-					pictogram, no
MAG					signal word, no
					hazard
					statement(s), no
					precautionary

Millipore- HCCBP3MAG-58K

Page 10 of 12



	statement(s) required.
--	------------------------

22-Cathepsin	22-Cathepsin D Beads - Magnetic								
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental				
identifier		Word	Statements	Statements	information				
Millipore -					No hazard				
HCTHPSND-					pictogram, no				
MAG					signal word, no				
					hazard				
					statement(s), no				
					precautionary				
					statement(s)				
					required.				

26-Fibroblast	26-Fibroblast Activation Protein a Beads - Magnetic							
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental			
identifier		Word	Statements	Statements	information			
Millipore -					Not a hazardous			
HFAP-MAG					substance or			
					mixture			
					according to			
					Regulation (EC)			
					No 1272/2008.			

34-MIA (Mela	34-MIA (Melanoma Inhibitory Activity Protein) Beads - Magnetic							
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental			
identifier		Word	Statements	Statements	information			
Millipore -					Not a hazardous			
HMIA-MAG					substance or			
					mixture			
					according to			
					Regulation (EC)			
					No 1272/2008.			

37- Myeloperoxidase Beads - Magnetic							
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information		
Millipore - HMP0-MAG					No hazard pictogram, no signal word, no hazard statement(s), no precautionary statement(s) required.		

42-HSHBG (Sex Hormone Binding Globulin) Beads - Magnetic						
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental	
identifier		Word	Statements	Statements	information	
Millipore -					Not a hazardous	
HSHBG-MAG					substance or	
					mixture	
					according to	

Millipore- HCCBP3MAG-58K Page 11 of 12



			Regulation (EC)	Ī
			No 1272/2008.	

72-IGFBP3 Be	72-IGFBP3 Beads - Magnetic								
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental				
identifier		Word	Statements	Statements	information				
Millipore -					No hazard				
HIGFBP3-MAG					pictogram, no				
					signal word, no				
					hazard				
					statement(s), no				
					precautionary				
					statement(s)				
					required.				

75-Ferritin Beads - Magnetic							
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information		
Millipore - HFRRTN-MAG					No hazard pictogram, no		
					signal word, no hazard		
					statement(s), no precautionary statement(s)		
					required.		

Millipore- HCCBP3MAG-58K Page 12 of 12

