

# **SAFETY DATA SHEET**

according to Regulation (EC) No. 1907/2006

Version 8.3 Revision Date 28.08.2023 Print Date 24.04.2025

GENERIC EU MSDS - NO COUNTRY SPECIFIC DATA - NO OEL DATA

#### SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1 Product identifiers

Product name : Milliplex Mouse Cytokine/Chemokine

MAGNETIC BEAD Premixed 32 Plex Kit.

Product Number : MCYTMAG-70K-PX32

Catalogue No. : 667478 Brand : Millipore

REACH No. : This product is a mixture. REACH Registration Number see

section 3.

1.2 Relevant identified uses of the substance or mixture and uses advised against

Identified uses : Biochemical research/analysis

1.3 Details of the supplier of the safety data sheet

Company :

1.4 Emergency telephone

Emergency Phone # : +(44)-870-8200418 (CHEMTREC (GB))

+(353)-19014670 (CHEMTREC Ireland) 001-803-017-9114 (CHEMTREC India)

This is a summary SDS for a kit. For the full SDS for each of the components listed in section 16 please visit our website.

#### **SECTION 2: Hazards identification**

#### 2.1 Classification of the substance or mixture

# Classification according to Regulation (EC) No 1272/2008

Acute toxicity, Oral (Category 4), H302 Acute toxicity, Inhalation (Category 4), H332 Acute toxicity, Dermal (Category 4), H312 Serious eye damage (Category 1), H318 Skin sensitization (Category 1), H317

Specific target organ toxicity - repeated exposure (Category 2), Respiratory Tract, H373 Long-term (chronic) aquatic hazard (Category 2), H411

For the full text of the H-Statements mentioned in this Section, see Section 16.

Millipore- MCYTMAG-70K-PX32

Page 1 of 14



#### 2.2 Label elements

# Labelling according Regulation (EC) No 1272/2008

Pictogram

Signal Word Danger

Hazard statement(s)

H302 + H312 + H332 Harmful if swallowed, in contact with skin or if inhaled.

H317 May cause an allergic skin reaction.

H318 Causes serious eye damage.

H373 May cause damage to organs (Respiratory Tract) through

prolonged or repeated exposure.

H411 Toxic to aquatic life with long lasting effects.

Precautionary statement(s)

P273 Avoid release to the environment.

P280 Wear protective gloves/ protective clothing/ eye protection/ face

protection.

P302 + P352 + P312 IF ON SKIN: Wash with plenty of water.Call a POISON CENTER/

doctor if you feel unwell.

P304 + P340 + P312 IF INHALED: Remove person to fresh air and keep comfortable

for breathing. Call a POISON CENTER/ doctor if you feel unwell.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue

rinsing.

P314 Get medical advice/ attention if you feel unwell.

Supplemental Hazard

Statements

none

Reduced Labeling (<= 125 ml)

Pictogram

Signal Word Danger

Hazard statement(s)

H317 May cause an allergic skin reaction. H318 Causes serious eye damage.

Precautionary statement(s)

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue

rinsing.

Supplemental Hazard

Statements

none

2.3 Other hazards - none

Millipore- MCYTMAG-70K-PX32

The life science business of Merck operates as MilliporeSigma in the US and Canada



Page 2 of 14

# **SECTION 3: Composition/information on ingredients**

Refer to component SDS

#### **SECTION 4: First aid measures**

Refer to component SDS

# **SECTION 5: Firefighting measures**

Refer to component SDS

#### **SECTION 6: Accidental release measures**

Refer to component SDS

# **SECTION 7: Handling and storage**

#### 7.1 Precautions for safe handling

#### Advice on safe handling

Work under hood. Do not inhale substance/mixture. Avoid generation of vapours/aerosols.

#### **Hygiene measures**

Immediately change contaminated clothing. Apply preventive skin protection. Wash hands and face after working with substance.

For precautions see section 2.2.

# 7.2 Conditions for safe storage, including any incompatibilities

#### Storage conditions

Tightly closed.

Recommended storage temperature see product label.

#### Storage class

Storage class (TRGS 510): 10: Combustible liquids

#### 7.3 Specific end use(s)

Apart from the uses mentioned in section 1.2 no other specific uses are stipulated

# **SECTION 8: Exposure controls/personal protection**

Refer to component SDS

#### **SECTION 9: Physical and chemical properties**

Refer to component SDS

#### **SECTION 10: Stability and reactivity**

Refer to component SDS

Millipore- MCYTMAG-70K-PX32

Page 3 of 14



# **SECTION 11: Toxicological information**

Refer to component SDS

# **SECTION 12: Ecological information**

Refer to component SDS

# **SECTION 13: Disposal considerations**

Refer to component SDS

# **SECTION 14: Transport information**

14.1 UN number

ADR/RID: - IMDG: - IATA: -

14.2 UN proper shipping name

ADR/RID: Not dangerous goods IMDG: Not dangerous goods IATA: Not dangerous goods

14.3 Transport hazard class(es)

ADR/RID: - IMDG: - IATA: -

14.4 Packaging group

ADR/RID: - IMDG: - IATA: -

14.5 Environmental hazards

ADR/RID: no IMDG Marine pollutant: no IATA: no

# 14.6 Special precautions for user

#### **Further information**

Not classified as dangerous in the meaning of transport regulations.

# **SECTION 15: Regulatory information**

# 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

This material safety data sheet complies with the requirements of Regulation (EC) No. 1907/2006.

Refer to component SDS

# 15.2 Chemical Safety Assessment

For this product a chemical safety assessment was not carried out

Millipore- MCYTMAG-70K-PX32 Page 4 of 14



# **SECTION 16: Other information**

Kit Components:			
Mouse Cytokine/Chemokine Std	Millipore	MXM8070	Acute Tox. 4; Eye Dam. 1; STOT RE 2; Aquatic Chronic 3; H302 + H312 + H332, H318, H373,
			H412
Mouse Cytokine/Chemokine Standard	Millipore	MXM8070-2	Acute Tox. 4; Eye Dam. 1; STOT RE 2; Aquatic Chronic 3; H302 + H312 + H332, H318, H373, H412
Mouse Cytokine/Chemokine Quality Control	Millipore	MXM6070	Acute Tox. 4; Eye Dam. 1; STOT RE 2; Aquatic Chronic 3; H302 + H312 + H332, H318, H373, H412
Mouse Cytokine/Chemokine Quality Control	Millipore	MXM6070-2	Acute Tox. 4; Eye Dam. 1; STOT RE 2; Aquatic Chronic 3; H302 + H312 + H332, H318, H373, H412
Serum MatrixMILLIPLEX COLDPKS w antisera	Millipore	MXMSM	Aquatic Chronic 3; H412
Milliplex Assay Buffer	Millipore	L-AB	
10X Wash Buffer - MILLIPLEX	Millipore	L-WB	Skin Sens. 1; Aquatic Chronic 2; H317, H411
Mouse Cytokine Detection Antibodies	Millipore	MXM1070-1	
Mouse Cytokine Detection Antibodies	Millipore	MXM1070-2	
Mouse Cytokine Detection Antibodies	Millipore	MXM1070-3	
Streptavidin- Phycoerythrin	Millipore	L-SAPE3	
Streptavidin- Phycoerythrin	Millipore	L-SAPE4	Eye Irrit. 2; H319
Streptavidin- Phycoerythrin	Millipore	L-SAPE10	Eye Irrit. 2; H319
Mouse Cytokine/Chemokine MAG Pmx Beads32	Millipore	MCYPMX32- MAG	

 $Millipore\hbox{--} MCYTMAG-70K-PX32$ 

A

Page 5 of 14

#### Full text of other abbreviations

ADN - European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways; ADR - Agreement concerning the International Carriage of Dangerous Goods by Road; AIIC - Australian Inventory of Industrial Chemicals; ASTM -American Society for the Testing of Materials; bw - Body weight; CMR - Carcinogen, Mutagen or Reproductive Toxicant; DIN - Standard of the German Institute for Standardisation; DSL - Domestic Substances List (Canada); ECx - Concentration associated with x% response; ELx - Loading rate associated with x% response; EmS -Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx -Concentration associated with x% growth rate response; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan); ISO - International Organisation for Standardization; KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; n.o.s. Not Otherwise Specified; NO(A)EC - No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention; PBT - Persistent, Bioaccumulative and Toxic substance; PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; RID - Regulations concerning the International Carriage of Dangerous Goods by Rail; SADT - Self-Accelerating Decomposition Temperature; SDS -Safety Data Sheet; TCSI - Taiwan Chemical Substance Inventory; TECI - Thailand Existing Chemicals Inventory; TSCA - Toxic Substances Control Act (United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative

Classification of the	Classification procedure:	
Acute Tox.4	H302	Calculation method
Acute Tox.4	H332	Calculation method
Acute Tox.4	H312	Calculation method
Eye Dam.1	H318	Calculation method
Skin Sens.1	H317	Calculation method
STOT RE2	H373	Calculation method
Aquatic Chronic2	H411	Calculation method

#### **Further information**

The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. The information in this document is based on the present state of our knowledge and is applicable to the product with regard to

Millipore- MCYTMAG-70K-PX32 Page 6 of 14



appropriate safety precautions. It does not represent any guarantee of the properties of the product. Sigma-Aldrich Corporation and its Affiliates shall not be held liable for any damage resulting from handling or from contact with the above product. See www.sigma-aldrich.com and/or the reverse side of invoice or packing slip for additional terms and conditions of sale.

Copyright 2020 Sigma-Aldrich Co. LLC. License granted to make unlimited paper copies for internal use only.

The branding on the header and/or footer of this document may temporarily not visually match the product purchased as we transition our branding. However, all of the information in the document regarding the product remains unchanged and matches the product ordered. For further information please contact mlsbranding@sial.com.

 $Millipore\hbox{--} MCYTMAG\hbox{--}70K\hbox{--}PX32$ 

The life science business of Merck operates as MilliporeSigma in the US and Canada  $\,$ 



Page 7 of 14



# **Kit Components:**

Mouse Cytoki	ne/Chemoki	ne Std			
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information
Millipore - MXM8070		Danger	Harmful if swallowed, in contact with skin or if inhaled. Causes serious eye damage. May cause damage to organs (Respiratory Tract) through prolonged or repeated exposure. Harmful to aquatic life with long lasting effects.	Avoid release to the environment. IF ON SKIN: Wash with plenty of water.Call a POISON CENTER/ doctor if you feel unwell. Wear protective gloves/ protective clothing/ eye protection. IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a POISON CENTER/ doctor if you feel unwell. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Get medical advice/ attention if you feel unwell.	IIIIOIIIIauoii


Millipore- MCYTMAG-70K-PX32



Mouse Cytok Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier	rictogram	Word	Statements	Statements	information
Millipore -		Danger	Harmful if	Avoid release to	
MXM8070-2			swallowed, in	the	
	X		contact with skin	environment. IF	
	(下至)		or if inhaled.	ON SKIN: Wash	
	X		Causes serious	with plenty of	
			eye damage.	water.Call a	
			May cause	POISON	
			damage to	CENTER/ doctor	
			organs	if you feel	
			(Respiratory	unwell. Wear	
			Tract) through	protective	
			prolonged or	gloves/	
			repeated	protective	
			exposure.	clothing/ eye	
			Harmful to	protection/ face	
			aquatic life with	protection. IF	
			long lasting	INHALED:	
			effects.	Remove person	
				to fresh air and	
				keep comfortable for	
				breathing. Call a	
				POISON	
				CENTER/ doctor	
				if you feel	
				unwell. IF IN	
				EYES: Rinse	
				cautiously with	
				water for	
				several minutes.	
				Remove contact	
				lenses, if	
				present and	
				easy to do.	
				Continue	
				rinsing. Get	
				medical advice/	
				attention if you	
				feel unwell.	

Mouse Cytoki	Mouse Cytokine/Chemokine Quality Control						
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental		
identifier		Word	Statements	Statements	information		
Millipore -		Danger	Harmful if	Avoid release to			
MXM6070			swallowed, in	the			
			contact with skin	environment. IF			
	<u> </u>		or if inhaled.	ON SKIN: Wash			
			Causes serious	with plenty of			

Page 9 of 14





ovo damaca	water.Call a	
eye damage.	POISON	
May cause		
damage to	CENTER/ doctor	
organs	if you feel	
(Respiratory	unwell. Wear	
Tract) through	protective	
prolonged or	gloves/	
repeated	protective	
exposure.	clothing/ eye	
Harmful to	protection/ face	
aquatic life with	protection. IF	
long lasting	INHALED:	
effects.	Remove person	
	to fresh air and	
	keep	
	comfortable for	
	breathing. Call a	
	POISON	
	CENTER/ doctor	
	if you feel	
	unwell. IF IN	
	EYES: Rinse	
	cautiously with	
	water for	
	several minutes.	
	Remove contact	
	lenses, if	
	present and	
	easy to do.	
	Continue	
	rinsing. Get	
	medical advice/	
	attention if you	
	feel unwell.	

<b>Mouse Cytok</b>	ine/Chemoki	ne Quality	Control		
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information
Millipore - MXM6070-2	!	Danger	Harmful if swallowed, in contact with skin or if inhaled. Causes serious eye damage. May cause damage to organs (Respiratory Tract) through prolonged or repeated	Avoid release to the environment. IF ON SKIN: Wash with plenty of water.Call a POISON CENTER/ doctor if you feel unwell. Wear protective gloves/protective	

Page 10 of 14



ovnosuro	clothing/ eye	
exposure.		
Harmful to	protection/ face	
aquatic life with	protection. IF	
long lasting	INHALED:	
effects.	Remove person	
	to fresh air and	
	keep	
	comfortable for	
	breathing. Call a	
	POISON	
	CENTER/ doctor	
	if you feel	
	unwell. IF IN	
	EYES: Rinse	
	cautiously with	
	water for	
	several minutes.	
	Remove contact	
	lenses, if	
	present and	
	easy to do.	
	Continue	
	rinsing. Get	
	medical advice/	
	-	
	attention if you	
	feel unwell.	

Serum Matri	Serum MatrixMILLIPLEX COLDPKS w antisera						
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental		
identifier		Word	Statements	Statements	information		
Millipore -			Harmful to	Avoid release to			
MXMSM			aquatic life with	the			
			long lasting	environment.			
			effects.	Dispose of			
				contents/			
				container to an			
				approved waste			
				disposal plant.			

Milliplex Assa	Milliplex Assay Buffer						
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental		
identifier		Word	Statements	Statements	information		
Millipore - L-					No hazard		
AB					pictogram, no		
					signal word, no		
					hazard		
					statement(s), no		
					precautionary		
					statement(s)		
					required		

The life science business of Merck operates as MilliporeSigma in the US and Canada  $\,$ 



Page 11 of 14

10X Wash But	tter - MILLIF	LEX	1	T	1
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Millipore - L- WB	<b>! ★</b>	Warning	May cause an allergic skin reaction. Toxic to aquatic life with long lasting effects.	Avoid breathing mist or vapors. IF ON SKIN: Wash with plenty of water. Contaminated work clothing should not be allowed out of the workplace. If skin irritation or rash occurs: Get medical advice/ attention. Avoid release to the environment. Wear protective gloves.	

Mouse Cytoki	Mouse Cytokine Detection Antibodies						
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information		
Millipore - MXM1070-1					Not a hazardous substance or mixture according to Regulation (EC) No 1272/2008.		

Mouse Cytokine Detection Antibodies							
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental		
identifier		Word	Statements	Statements	information		
Millipore -					Not a hazardous		
MXM1070-2					substance or		
					mixture		
					according to		
					Regulation (EC)		
					No 1272/2008.		

Mouse Cytoki	Mouse Cytokine Detection Antibodies							
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information			
Millipore - MXM1070-3					Not a hazardous substance or mixture according to			

Page 12 of 14



		Regulation (EC)
		No 1272/2008.

Streptavidin-Phycoerythrin						
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental	
identifier		Word	Statements	Statements	information	
Millipore - L-					Not a hazardous	
SAPE3					substance or	
					mixture	
					according to	
					Regulation (EC)	
					No 1272/2008.	

Streptavidin-Phycoerythrin						
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental	
identifier		Word	Statements	Statements	information	
Millipore - L- SAPE4	•	Warning	Causes serious eye irritation.	Wash skin thoroughly after handling. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Wear eye protection/ face protection. If eye irritation persists: Get medical advice/ attention.		

Streptavidin-Phycoerythrin						
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information	
Millipore - L- SAPE10	1	Warning	Causes serious eye irritation.	Wash skin thoroughly after handling. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue		

Page 13 of 14



rinsing. Wear
eye protection/
face protection.
If eye irritation
persists: Get
medical advice/
attention.

Mouse Cytok	Mouse Cytokine/Chemokine MAG Pmx Beads32						
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information		
Millipore - MCYPMX32- MAG					Not a hazardous substance or mixture according to Regulation (EC) No 1272/2008.		

