

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Version 8.0 Revision Date 17.05.2024 Print Date 27.04.2025

GENERIC EU MSDS - NO COUNTRY SPECIFIC DATA - NO OEL DATA

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1 Product identifiers

Product name : Feline Cytokine MAG Premix 19 Plex BULK

Product Number : FCYTMAG20KPX19BK

Catalogue No. : 6D3102 Brand : Millipore

REACH No. : This product is a mixture. REACH Registration Number see

section 3.

1.2 Relevant identified uses of the substance or mixture and uses advised against

Identified uses : Use restricted under TSCA to research and development or as

analytical reagent. Uses regulated under FDA or FIFRA are not

affected., Biochemical research/analysis

1.3 Details of the supplier of the safety data sheet

Company :

1.4 Emergency telephone

Emergency Phone #: +(44)-870-8200418 (CHEMTREC (GB))

+(353)-19014670 (CHEMTREC Ireland) 001-803-017-9114 (CHEMTREC India)

This is a summary SDS for a kit. For the full SDS for each of the components listed in section 16 please visit our website.

SECTION 2: Hazards identification

2.1 Classification of the substance or mixture

Acute toxicity, (Category 4) H302: Harmful if swallowed.

Acute toxicity, (Category 4) H332: Harmful if inhaled.

Acute toxicity, (Category 3) H311: Toxic in contact with skin.

Eye irritation, (Category 2) H319: Causes serious eye irritation.

Skin sensitization, (Category 1) H317: May cause an allergic skin reaction.

Millipore- FCYTMAG20KPX19BK

Page 1 of 11



Long-term (chronic) aquatic hazard, (Category 2)

H411: Toxic to aquatic life with long lasting

effects.

2.2 Label elements

Labelling according Regulation (EC) No 1272/2008

Pictogram

V V

Signal Word Danger

Hazard Statements

H302 + H332 Harmful if swallowed or if inhaled.

H311 Toxic in contact with skin.

H317 May cause an allergic skin reaction.

H319 Causes serious eye irritation.

H411 Toxic to aquatic life with long lasting effects.

Precautionary Statements

P273 Avoid release to the environment.

P280 Wear protective gloves/ protective clothing/ eye protection/ face

protection.

P301 + P312 IF SWALLOWED: Call a POISON CENTER/ doctor if you feel

ınwell.

P302 + P352 + P312 IF ON SKIN: Wash with plenty of water.Call a POISON CENTER/

doctor if you feel unwell.

P304 + P340 + P312 IF INHALED: Remove person to fresh air and keep comfortable

for breathing. Call a POISON CENTER/ doctor if you feel unwell.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue

rinsing.

Supplemental Hazard

Statements

none

Reduced Labeling (<= 125 ml)

Pictogram

Signal Word Danger

Hazard Statements

H311 Toxic in contact with skin.

H317 May cause an allergic skin reaction.

Precautionary Statements

P302 + P352 + P312 IF ON SKIN: Wash with plenty of water. Call a POISON CENTER/

doctor if you feel unwell.

Supplemental Hazard

Statements

none

2.3 Other hazards - none

Millipore- FCYTMAG20KPX19BK

Page 2 of 11



SECTION 3: Composition/information on ingredients

Refer to component SDS

SECTION 4: First aid measures

Refer to component SDS

SECTION 5: Firefighting measures

Refer to component SDS

SECTION 6: Accidental release measures

Refer to component SDS

SECTION 7: Handling and storage

7.1 Precautions for safe handling

Advice on safe handling

Work under hood. Do not inhale substance/mixture. Avoid generation of vapours/aerosols.

Hygiene measures

Immediately change contaminated clothing. Apply preventive skin protection. Wash hands and face after working with substance.

For precautions see section 2.2.

7.2 Conditions for safe storage, including any incompatibilities

Storage conditions

Tightly closed. Keep in a well-ventilated place. Keep locked up or in an area accessible only to qualified or authorized persons.

Storage stabilityRecommended storage temperature

2 - 8 °C

Storage class

Storage class (TRGS 510): 6.1C: Combustible, acute toxic Cat.3 / toxic compounds or compounds which causing chronic effects

7.3 Specific end use(s)

Apart from the uses mentioned in section 1.2 no other specific uses are stipulated

SECTION 8: Exposure controls/personal protection

Refer to component SDS

Millipore- FCYTMAG20KPX19BK

Canada

A

Page 3 of 11

SECTION 9: Physical and chemical properties

Refer to component SDS

SECTION 10: Stability and reactivity

Refer to component SDS

SECTION 11: Toxicological information

Refer to component SDS

SECTION 12: Ecological information

Refer to component SDS

SECTION 13: Disposal considerations

Refer to component SDS

SECTION 14: Transport information

14.1 UN number

ADR/RID: 3316 IMDG: 3316 IATA: 3316

14.2 UN proper shipping name

ADR/RID: CHEMICAL KIT IMDG: CHEMICAL KIT Chemical kit

14.3 Transport hazard class(es)

ADR/RID: 9 IMDG: 9 IATA: 9

14.4 Packaging group

ADR/RID: III IMDG: III IATA: III

14.5 Environmental hazards

ADR/RID: no IMDG Marine pollutant: no IATA: no

14.6 Special precautions for user

Further information : No data available

 $Millipore\hbox{--} FCYTMAG20KPX19BK$

A

Page 4 of 11

SECTION 15: Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

This material safety data sheet complies with the requirements of Regulation (EC) No. 1907/2006.

Refer to component SDS

15.2 Chemical Safety Assessment

For this product a chemical safety assessment was not carried out

SECTION 16: Other information

Kit Components:			
Feline Cytokine Standard	Millipore	FCY-8020	Acute Tox. 4; Acute Tox. 3; Aquatic Chronic 2; H302 + H332, H311, H411
Feline Cytokine QC 1 & 2	Millipore	FCY-6020	Acute Tox. 4; Acute Tox. 3; Aquatic Chronic 2; H302 + H332, H311, H411
Serum Matrix - MILLIPLEX - COLD PAC	Millipore	LMC-SD	Aquatic Chronic 3; H412
Milliplex Assay Buffer	Millipore	L-AB	
10X Wash Buffer - MILLIPLEX	Millipore	L-WB	Skin Sens. 1; Aquatic Chronic 2; H317, H411
Feline Cytokine Detection Antibodies	Millipore	FCY-1020	Eye Irrit. 2; H319
Streptavidin- Phycoerythrin (3.5 ml) MILLI	Millipore	L-SAPE6	Eye Irrit. 2; H319
Feline Cytokine Premixed Beads	Millipore	FCYPMX19-MAG	

 $Millipore\hbox{--} FCYTMAG20KPX19BK$

M

Page 5 of 11

Full text of other abbreviations

ADN - European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways; ADR - Agreement concerning the International Carriage of Dangerous Goods by Road; AIIC - Australian Inventory of Industrial Chemicals; ASTM -American Society for the Testing of Materials; bw - Body weight; CMR - Carcinogen, Mutagen or Reproductive Toxicant; DIN - Standard of the German Institute for Standardisation; DSL - Domestic Substances List (Canada); ECx - Concentration associated with x% response; ELx - Loading rate associated with x% response; EmS -Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx -Concentration associated with x% growth rate response; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan); ISO - International Organisation for Standardization; KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; n.o.s. -Not Otherwise Specified; NO(A)EC - No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention; PBT - Persistent, Bioaccumulative and Toxic substance; PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; RID - Regulations concerning the International Carriage of Dangerous Goods by Rail; SADT - Self-Accelerating Decomposition Temperature; SDS -Safety Data Sheet; TCSI - Taiwan Chemical Substance Inventory; TECI - Thailand Existing Chemicals Inventory; TSCA - Toxic Substances Control Act (United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative

Classification of the	Classification procedure:	
Acute Tox.4	H302	Calculation method
Acute Tox.4	H332	Calculation method
Acute Tox.3	H311	Calculation method
Eye Irrit.2	H319	Calculation method
Skin Sens.1	H317	Calculation method
Aquatic Chronic2	H411	Calculation method

Further information

The information is believed to be correct but is not exhaustive and will be used solely as a guideline, which is based on current knowledge of the chemical substance or mixture and is applicable to appropriate safety precautions for the product. It does not represent any guarantee of the properties of the product. Sigma-Aldrich Corporation and its Affiliates shall not be held liable for any damage resulting from handling or from contact

Millipore- FCYTMAG20KPX19BK

A

with the above product. See www.sigma-aldrich.com and/or the reverse side of invoice or packing slip for additional terms and conditions of sale.

Copyright 2020 Sigma-Aldrich Co. LLC. License granted to make unlimited paper copies for internal use only.

The branding on the header and/or footer of this document may temporarily not visually match the product purchased as we transition our branding. However, all of the information in the document regarding the product remains unchanged and matches the product ordered. For further information please contact mlsbranding@sial.com.

 $Millipore\hbox{--} FCYTMAG20KPX19BK$

The life science business of Merck operates as MilliporeSigma in the US and Canada



Page 7 of 11




Kit Components:

Feline Cytokin	e Standard				
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
		_		-	

Feline Cytokine QC 1 & 2					
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Millipore -		Danger	Harmful if	Avoid breathing	
FCY-6020	335		swallowed or if	dust. IF	
	AV.		inhaled.	SWALLOWED:	
	***		Toxic in contact	Call a POISON	
	~		with skin.	CENTER/ doctor	
			Toxic to aquatic	if you feel	
			life with long	unwell. Avoid	

 $Millipore\hbox{--} FCYTMAG20KPX19BK$

Page 8 of 11





lasting effects.	release to the	
	environment. IF	
	ON SKIN: Wash	
	with plenty of	
	water.Call a	
	POISON	
	CENTER/ doctor	
	if you feel	
	unwell. Wear	
	protective	
	gloves/	
	protective	
	clothing. IF	
	INHALED:	
	Remove person	
	to fresh air and	
	keep	
	comfortable for	
	breathing. Call a	
	POISON	
	CENTER/ doctor	
	if you feel	
	unwell.	

Serum Matrix - MILLIPLEX - COLD PAC					
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Millipore -			Harmful to	Avoid release to	
LMC-SD			aquatic life with	the	
			long lasting	environment.	
			effects.	Dispose of	
				contents/	
				container to an	
				approved waste	
				disposal plant.	

Milliplex Assay Buffer					
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Millipore - L-					No hazard
AB					pictogram, no
					signal word, no
					hazard
					statement(s), no
					precautionary
					statement(s)
					required

 $Millipore\hbox{--} FCYTMAG20KPX19BK$

A

10X Wash Buffer - MILLIPLEX					
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information
Millipore - L- WB	**************************************	Warning	May cause an allergic skin reaction. Toxic to aquatic life with long lasting effects.	Avoid breathing mist or vapors. IF ON SKIN: Wash with plenty of water. Contaminated work clothing should not be allowed out of the workplace. If skin irritation or rash occurs: Get medical advice/ attention. Avoid release to the environment. Wear protective gloves.	

Feline Cytok	Feline Cytokine Detection Antibodies					
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental	
identifier		Word	Statements	Statements	information	
Millipore -		Warning	Causes serious	Wash skin		
FCY-1020			eye irritation.	thoroughly after		
	•			handling. IF IN		
				EYES: Rinse		
				cautiously with		
				water for		
				several minutes.		
				Remove contact		
				lenses, if		
				present and		
				easy to do.		
				Continue		
				rinsing. Wear		
				eye protection/		
				face protection.		
				If eye irritation		
				persists: Get		
				medical advice/		
				attention.		

 $Millipore\hbox{--} FCYTMAG20KPX19BK$

M

Streptavidin-I	Streptavidin-Phycoerythrin (3.5 ml) MILLI					
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental	
identifier		Word	Statements	Statements	information	
Millipore - L-		Warning	Causes serious	Wash skin		
SAPE6	•		eye irritation.	thoroughly after		
	•			handling. IF IN		
				EYES: Rinse		
				cautiously with		
				water for		
				several minutes.		
				Remove contact		
				lenses, if		
				present and		
				easy to do.		
				Continue		
				rinsing. Wear		
				eye protection/		
				face protection.		
				If eye irritation		
				persists: Get		
				medical advice/		
				attention.		

Feline Cytok	Feline Cytokine Premixed Beads				
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Millipore -					No hazard
FCYPMX19-					pictogram, no
MAG					signal word, no
					hazard
					statement(s), no
					precautionary
					statement(s)
					required

 $Millipore\hbox{--} FCYTMAG20KPX19BK$

A