

SAFETY DATA SHEET

according to Regulation (EC) No. 1907/2006

Version 8.6 Revision Date 20.02.2025 Print Date 27.04.2025

GENERIC EU MSDS - NO COUNTRY SPECIFIC DATA - NO OEL DATA

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1 Product identifiers

Product name : Canine Pituitary Expanded Panel

Product Number : CANPIT-96K Catalogue No. : 6E9654 Brand : Millipore

REACH No. :

1.2 Relevant identified uses of the substance or mixture and uses advised against

Identified uses : Reagent for development and research

1.3 Details of the supplier of the safety data sheet

Company :

1.4 Emergency telephone

Emergency Phone #: +(44)-870-8200418 (CHEMTREC (GB))

+(353)-19014670 (CHEMTREC Ireland) 001-803-017-9114 (CHEMTREC India)

This is a summary SDS for a kit. For the full SDS for each of the components listed in section 16 please visit our website.

SECTION 2: Hazards identification

2.1 Classification of the substance or mixture

Acute toxicity, (Category 4) H302: Harmful if swallowed.

Acute toxicity, (Category 4) H332: Harmful if inhaled.

Acute toxicity, (Category 4) H312: Harmful in contact with skin.

Serious eye damage, (Category H318: Causes serious eye damage.

1)

Skin sensitization, (Category 1) H317: May cause an allergic skin reaction.

Specific target organ toxicity - H373: May cause damage to organs repeated exposure, (Category 2), through prolonged or repeated exposure.

Respiratory Tract

Long-term (chronic) aquatic H411: Toxic to aquatic life with long lasting

Millipore- CANPIT-96K Page 1 of 13



2.2 Label elements

Labelling according Regulation (EC) No 1272/2008

Pictogram

Signal Word Danger

Hazard Statements

H302 + H312 + H332 Harmful if swallowed, in contact with skin or if inhaled.

H317 May cause an allergic skin reaction.

H318 Causes serious eye damage.

H373 May cause damage to organs (Respiratory Tract) through

prolonged or repeated exposure.

H411 Toxic to aquatic life with long lasting effects.

Precautionary Statements

P273 Avoid release to the environment.

P280 Wear protective gloves/ protective clothing/ eye protection/ face

protection.

P302 + P352 + P312 IF ON SKIN: Wash with plenty of water.Call a POISON CENTER/

doctor if you feel unwell.

P304 + P340 + P312 IF INHALED: Remove person to fresh air and keep comfortable

for breathing. Call a POISON CENTER/ doctor if you feel unwell.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue

rinsing.

P314 Get medical advice/ attention if you feel unwell.

Supplemental Hazard

Statements

none

Reduced Labeling (<= 125 ml)

Pictogram

Signal Word Danger

Hazard Statements

H317 May cause an allergic skin reaction.

H318 Causes serious eye damage.

Precautionary Statements

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue

rinsing.

Supplemental Hazard

Statements

none

2.3 Other hazards - none

SECTION 3: Composition/information on ingredients

Refer to component SDS

Millipore- CANPIT-96K Page 2 of 13

SECTION 4: First aid measures

Refer to component SDS

SECTION 5: Firefighting measures

Refer to component SDS

SECTION 6: Accidental release measures

Refer to component SDS

SECTION 7: Handling and storage

7.1 Precautions for safe handling

Advice on safe handling

Work under hood. Do not inhale substance/mixture. Avoid generation of vapours/aerosols.

Hygiene measures

Immediately change contaminated clothing. Apply preventive skin protection. Wash hands and face after working with substance.

For precautions see section 2.2.

7.2 Conditions for safe storage, including any incompatibilities

Storage conditions

Tightly closed.

Recommended storage temperature see product label.

Storage class

Storage class (TRGS 510): 10: Combustible liquids

7.3 Specific end use(s)

Apart from the uses mentioned in section 1.2 no other specific uses are stipulated

SECTION 8: Exposure controls/personal protection

Refer to component SDS

SECTION 9: Physical and chemical properties

Refer to component SDS

SECTION 10: Stability and reactivity

Refer to component SDS

SECTION 11: Toxicological information

Refer to component SDS

Millipore- CANPIT-96K Page 3 of 13

A

SECTION 12: Ecological information

Refer to component SDS

SECTION 13: Disposal considerations

Refer to component SDS

SECTION 14: Transport information

14.1 UN number

ADR/RID: - IMDG: - IATA: -

14.2 UN proper shipping name

ADR/RID: Not dangerous goods IMDG: Not dangerous goods IATA: Not dangerous goods

14.3 Transport hazard class(es)

ADR/RID: - IMDG: - IATA: -

14.4 Packaging group

ADR/RID: - IMDG: - IATA: -

14.5 Environmental hazards

ADR/RID: no IMDG Marine pollutant: no IATA: no

14.6 Special precautions for user

Further information

Not classified as dangerous in the meaning of transport regulations.

14.7 Maritime transport in bulk according to IMO instruments

Not applicable for product as supplied.

SECTION 15: Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

This material safety data sheet complies with the requirements of Regulation (EC) No. 1907/2006.

Refer to component SDS

15.2 Chemical Safety Assessment

For this product a chemical safety assessment was not carried out

A

Millipore- CANPIT-96K Page 4 of 13

SECTION 16: Other information

Kit Components:			
10X Wash Buffer - MILLIPLEX	Millipore	L-WB	Skin Sens. 1; Aquatic Chronic 2; H317, H411
Milliplex Assay Buffer	Millipore	L-AB	,
Bead Diluent COLD PACKS	Millipore	LBD	
Streptavidin- Phycoerythrin	Millipore	L-SAPE	Eye Irrit. 2; H319
Canine Pituitary Expanded Panel Standard	Millipore	CANP-8096	Acute Tox. 4; Eye Dam. 1; STOT RE 2; Aquatic Chronic 3; H302 + H312 + H332, H318, H373, H412
Canine Pituitary Expanded Panel Quality Control 1	Millipore	CANP-6096	Acute Tox. 4; Eye Dam. 1; STOT RE 2; Aquatic Chronic 3; H302 + H312 + H332, H318, H373, H412
Canine Pituitary Expanded Panel Quality Control 2	Millipore	CANP-7096	Acute Tox. 4; Eye Dam. 1; STOT RE 2; Aquatic Chronic 3; H302 + H312 + H332, H318, H373, H412
Canine Pituitary Expanded Panel Detection Antibodies	Millipore	CANP-1096	Eye Irrit. 2; H319
Serum Matrix	Millipore	CANP-SM	Aquatic Chronic 3; H412
13-ACTH MAGNETIC Beads - MILLIPLEX	Millipore	RACTH-MAG	
26-FSH MAGNETIC Beads - MILLIPLEX	Millipore	RFSH-MAG	
28-GH Beads Magnetic	Millipore	CANPGH-MAG	
62-TSH MAGNETIC Beads - MILLIPLEX	Millipore	RTSH-MAG	
64-LH MAGNETIC Beads - MILLIPLEX	Millipore	CLH-MAG	
66-Prolactin Beads Magnetic	Millipore	CANPHPL-MAG	

Millipore- CANPIT-96K Page 5 of 13



Full text of other abbreviations

ADN - European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways; ADR - Agreement concerning the International Carriage of Dangerous Goods by Road; AIIC - Australian Inventory of Industrial Chemicals; ASTM -American Society for the Testing of Materials; bw - Body weight; CMR - Carcinogen, Mutagen or Reproductive Toxicant; DIN - Standard of the German Institute for Standardisation; DSL - Domestic Substances List (Canada); ECx - Concentration associated with x% response; ELx - Loading rate associated with x% response; EmS -Emergency Schedule; ENCS - Existing and New Chemical Substances (Japan); ErCx -Concentration associated with x% growth rate response; GHS - Globally Harmonized System; GLP - Good Laboratory Practice; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; ICAO - International Civil Aviation Organization; IECSC - Inventory of Existing Chemical Substances in China; IMDG - International Maritime Dangerous Goods; IMO - International Maritime Organization; ISHL - Industrial Safety and Health Law (Japan); ISO - International Organisation for Standardization; KECI - Korea Existing Chemicals Inventory; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); MARPOL - International Convention for the Prevention of Pollution from Ships; n.o.s. -Not Otherwise Specified; NO(A)EC - No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NOELR - No Observable Effect Loading Rate; NZIoC - New Zealand Inventory of Chemicals; OECD - Organization for Economic Co-operation and Development; OPPTS - Office of Chemical Safety and Pollution Prevention; PBT - Persistent, Bioaccumulative and Toxic substance; PICCS - Philippines Inventory of Chemicals and Chemical Substances; (Q)SAR - (Quantitative) Structure Activity Relationship; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; RID - Regulations concerning the International Carriage of Dangerous Goods by Rail; SADT - Self-Accelerating Decomposition Temperature; SDS -Safety Data Sheet; TCSI - Taiwan Chemical Substance Inventory; TECI - Thailand Existing Chemicals Inventory; TSCA - Toxic Substances Control Act (United States); UN - United Nations; UNRTDG - United Nations Recommendations on the Transport of Dangerous Goods; vPvB - Very Persistent and Very Bioaccumulative

Classification of the	e mixture	Classification procedure:
Acute Tox.4	H302	Calculation method
Acute Tox.4	H332	Calculation method
Acute Tox.4	H312	Calculation method
Eye Dam.1	H318	Calculation method
Skin Sens.1	H317	Calculation method
STOT RE2	H373	Calculation method
Aquatic Chronic2	H411	Calculation method

Millipore- CANPIT-96K Page 6 of 13



Further information

The information is believed to be correct but is not exhaustive and will be used solely as a guideline, which is based on current knowledge of the chemical substance or mixture and is applicable to appropriate safety precautions for the product. It does not represent any guarantee of the properties of the product. Sigma-Aldrich Corporation and its Affiliates shall not be held liable for any damage resulting from handling or from contact with the above product. See www.sigma-aldrich.com and/or the reverse side of invoice or packing slip for additional terms and conditions of sale.

Copyright 2020 Sigma-Aldrich Co. LLC. License granted to make unlimited paper copies for internal use only.

The branding on the header and/or footer of this document may temporarily not visually match the product purchased as we transition our branding. However, all of the information in the document regarding the product remains unchanged and matches the product ordered. For further information please contact mlsbranding@sial.com.

Millipore- CANPIT-96K Page 7 of 13





Kit Components:

10X Wash Buf	fer - MILLIP	LEX			
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Millipore - L- WB		Warning	May cause an allergic skin reaction. Toxic to aquatic life with long lasting effects.	Avoid breathing mist or vapors. IF ON SKIN: Wash with plenty of water. Contaminated work clothing should not be allowed out of the workplace. If skin irritation or rash occurs: Get medical advice/ attention. Avoid release to the environment. Wear protective gloves.	

Milliplex Assay Buffer					
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Millipore - L-					No hazard
AB					pictogram, no
					signal word, no
					hazard
					statement(s), no
					precautionary
					statement(s)
					required.

Bead Diluent	Bead Diluent COLD PACKS				
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information
Millipore - LBD					No hazard pictogram, no signal word, no hazard statement(s), no precautionary statement(s) required.

Millipore- CANPIT-96K Page 8 of 13



Streptavidin-F	Streptavidin-Phycoerythrin				
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Millipore - L- SAPE	<u>(1)</u>	Warning	Causes serious eye irritation.	Wash skin thoroughly after handling. IF IN EYES: Rinse cautiously with water for several minutes.	
				Remove contact lenses, if present and easy to do. Continue rinsing. Wear eye protection/ face protection. If eye irritation persists: Get medical advice/ attention.	

Canine Pituit	ary Expanded	d Panel Sta	<u>in</u> dard		
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Millipore -		Danger	Harmful if	Avoid release to	
CANP-8096			swallowed, in	the	
			contact with skin	environment. IF	
			or if inhaled.	ON SKIN: Wash	
			Causes serious	with plenty of	
	! >		eye damage.	water.Call a	
	~		May cause	POISON	
			damage to	CENTER/ doctor	
			organs	if you feel	
			(Respiratory	unwell. Wear	
			Tract) through	protective	
			prolonged or	gloves/	
			repeated	protective	
			exposure. Harmful to	clothing/ eye protection/ face	
			aquatic life with	protection, IF	
			long lasting	INHALED:	
			effects.	Remove person	
			circus.	to fresh air and	
				keep	
				comfortable for	
				breathing. Call a	
				POISON	
				CENTER/ doctor	
				if you feel	
				unwell. IF IN	
				EYES: Rinse	
				cautiously with	
				water for	

Millipore- CANPIT-96K Page 9 of 13



several minutes.
Remove contact
lenses, if
present and
easy to do.
Continue
rinsing. Get
medical advice/
attention if you
feel unwell.

Canine Pituita	Canine Pituitary Expanded Panel Quality Control 1					
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental	
identifier		Word	Statements	Statements	information	
Millipore -		Danger	Harmful if	Avoid release to		
CANP-6096			swallowed, in	the		
			contact with skin	environment. IF		
	<u> </u>		or if inhaled.	ON SKIN: Wash		
	X		Causes serious	with plenty of		
	(!)		eye damage.	water.Call a		
	_		May cause	POISON		
			damage to	CENTER/ doctor		
			organs	if you feel		
			(Respiratory	unwell. Wear		
			Tract) through	protective		
			prolonged or	gloves/		
			repeated	protective		
			exposure. Harmful to	clothing/ eye		
			aquatic life with	protection/ face protection. IF		
			long lasting	INHALED:		
			effects.	Remove person		
			errects.	to fresh air and		
				keep		
				comfortable for		
				breathing. Call a		
				POISON		
				CENTER/ doctor		
				if you feel		
				unwell. IF IN		
				EYES: Rinse		
				cautiously with		
				water for		
				several minutes.		
				Remove contact		
				lenses, if		
				present and		
				easy to do.		
				Continue		
				rinsing. Get		
				medical advice/		
				attention if you		
				feel unwell.		

Millipore- CANPIT-96K Page 10 of 13



			ality Control 2	1	1
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Millipore -		Danger	Harmful if	Avoid release to	
CANP-7096			swallowed, in	the	
			contact with skin	environment. IF	
	<u> </u>		or if inhaled.	ON SKIN: Wash	
	X		Causes serious	with plenty of	
	(!)		eye damage.	water.Call a	
			May cause	POISON	
			damage to	CENTER/ doctor	
			organs	if you feel	
			(Respiratory	unwell. Wear	
			Tract) through	protective	
			prolonged or	gloves/	
			repeated	protective	
			exposure.	clothing/ eye	
			Harmful to	protection/ face	
			aquatic life with	protection. IF	
			long lasting	INHALED:	
			effects.	Remove person	
				to fresh air and	
				keep	
				comfortable for	
				breathing. Call a	
				POISON	
				CENTER/ doctor	
				if you feel	
				unwell. IF IN	
				EYES: Rinse	
				cautiously with	
				water for	
				several minutes.	
				Remove contact	
				lenses, if	
				present and	
				easy to do.	
				Continue	
				rinsing. Get	
				medical advice/	
				attention if you	
				feel unwell.	

Canine Pitui	tary Expande	d Panel Det	ection Antibodies	}	
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information
Millipore - CANP-1096	!	Warning	Causes serious eye irritation.	Wash skin thoroughly after handling. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if	

Millipore- CANPIT-96K Page 11 of 13



	present and easy to do. Continue	
	rinsing. Wear eye protection/ face protection. If eye irritation persists: Get	
	medical advice/ attention.	

Serum Matrix					
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Millipore -			Harmful to	Avoid release to	
CANP-SM			aquatic life with	the	
			long lasting	environment.	
			effects.	Dispose of	
				contents/	
				container to an	
				approved waste	
				disposal plant.	

13-ACTH MAGNETIC Beads - MILLIPLEX					
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental
identifier		Word	Statements	Statements	information
Millipore -					No hazard
RACTH-MAG					pictogram, no
					signal word, no
					hazard
					statement(s), no
					precautionary
					statement(s)
					required.

26-FSH MAGI	26-FSH MAGNETIC Beads - MILLIPLEX					
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information	
Millipore - RFSH-MAG					No hazard pictogram, no signal word, no hazard statement(s), no precautionary statement(s) required.	

28-GH Beads Magnetic						
Product identifier	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information	
Millipore - CANPGH-MAG					No hazard pictogram, no signal word, no hazard statement(s), no	

Millipore- CANPIT-96K Page 12 of 13

A

		precautionary
		statement(s)
		required.

62-TSH MAGNETIC Beads - MILLIPLEX						
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental	
identifier		Word	Statements	Statements	information	
Millipore -					No hazard	
RTSH-MAG					pictogram, no	
					signal word, no	
					hazard	
					statement(s), no	
					precautionary	
					statement(s)	
					required.	

64-LH MAGNETIC Beads - MILLIPLEX						
Product identifier Millipore - CLH-MAG	Pictogram	Signal Word	Hazard Statements	Precautionary Statements	Supplemental information No hazard pictogram, no signal word, no hazard statement(s), no	
					precautionary statement(s) required.	

66-Prolactin Beads Magnetic						
Product	Pictogram	Signal	Hazard	Precautionary	Supplemental	
identifier		Word	Statements	Statements	information	
Millipore -					No hazard	
CANPHPL-					pictogram, no	
MAG					signal word, no	
					hazard	
					statement(s), no	
					precautionary	
					statement(s)	
					required.	

Millipore- CANPIT-96K Page 13 of 13

