

100 μL Gastight Syringe Model 1710 RN, Small Removable Needle, 22s gauge, 2 in, point style 2

Specification Sheet (Part/REF # 81030)

Story of the story

Split/splitless injection consists of a heated chamber with a glass liner into which the sample is injected through the septum. A microsyringe is used to inject the sample through a rubber septum into a flash vaporizer chamber at the head

Product Specifications

of the column.

Volume	100 μL
Syringe Type	Gastight
Syringe Series	1700
Termination	Removable Needle (RN)
Needle	Small RN
Gauge	22s gauge
Point Style	2
Needle Length	2 inch (51 mm)
Fluid Path	SST, borosilicate glass, PTFE
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	1000 psig (68.9 bar)
Maximum Temperature	115°C / 239°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.057 inches (1.46 mm)



US: +1-800-648-5950 EU: +40-356-635-055



Barrel Outer Diameter	0.305 inches (7.75 mm)	
Needle Inner Diameter	0.0059 inches (0.15 mm)	
Needle Outer Diameter	0.0283 inches (0.72 mm)	
Std Needle Dead Volume	1.13 μL	
Sales Quantity	1 syringe per pack	
Specifications are subject to change without notice		Spec. Version B

https://www.hamiltoncompany.com/laboratory-products/syringes/81030



US: +1-800-648-5950 EU: +40-356-635-055