

DRY BATH INCUBATOR



MINI DRY BATH INCUBATOR





FOUR E'S



Features:

- Temperature is up to 100°C
- LED display for easy reading
- Various applications with different blocks
- \bullet High temperature precision is up to $\pm~0.5^{\circ}\text{C}$
- Programmable operation can store up t o 5 programs
- 3D digital wheel adjustment of speed and temperature
- Simultaneous display of set /actual time and temperature
- Low operation voltage (12v) can be used directly wit h car charger

SPECIFICATIONS:

Model	TC0401005		
Temperature control range [°C]	RT + 5°C~100°C		
Timing range	0~99h59min		
Temp. control accuracy [°C]	±0.5		
Display accuracy [°C]	±0.1		
Heating time (40°C to100°C temp.)	≤15min		
Dimension (mm)	164 x 116 x 98		
Net weight	0.9kg		

MINI DRY BATH INCUBATOR WITH HEATING LID





Features:

- Heating lid
- LED display for easy reading
- Temperature is up to 100°C
- Various applications with different blocks
- High temperature precision is up to ± 0.5°C
- Programmable operation can store up t o 5 programs
- 3D digital wheel adjustment of speed and temperature
- Simultaneous display of set /actual time and temperature
- Low operation voltage (12v) can be used directly wit h car charger

SPECIFICATIONS:

Of Edit Ida Hollo.				
Model	TC0401006			
Temperature control range [°C]	RT + 5°C~100°C			
Timing range	0~99h59min			
Temp. control accuracy [°C]	±0.5			
Display accuracy [°C]	±0.1			
Heating time (40°C to100°C temp.)	≤15min			
Heating lid				
Dimension (mm)	164 x 116 x 98			
Net weight	0.9kg			

ACCESSORIES:



TC0401005-30 24 x 0.5ml block



TC0401005-31 40 x 0.2ml block



TC0401005-32 15 x 1.5ml block



TC0401005-33 15 x 2.0ml block



TC0401005-45 8-12.5x12.5x32mm (used in cuvette)



4 x 15ml block

TC0401005-100 2 x 50ml block





ACCESSORIES

ACCESSORIES				
Part No.	Description			
TC0401001-35	single block, Eppendorf tubes, 20 x 1.5mL, Φ11.5mm x 36.9mm			
TC0401001-54	single block, Conical tubes, 12 x 15ml, Φ17.1mm x 44.5mm			
TC0401001-55	single block, vials, 6 x 50ml Φ29.0 x 47.6mm			
TC0401001-60	single block, PCR tube, 30 x 0.5mL, Φ7.9mm x 27.6mm			
TC0401001-61	single block, microtube combination, 6/10/5 x 0.5/1.5/2ml			
TC0401001-62	single block, Microtubes, 20 x 1.5mL, Φ11.1mm x 39.1mm			
TC0401001-63	single block, Eppendorf tubes, 20 x 2.0ml, Φ11.5mm x 38.1mm			
TC0401001-64	single block, Corning tubes, 20 x 2.0ml, Φ10.9mm x 38.1mm			
TC0401001-66	single block, centifuge tube combinat ion, 4/3/2 x 1.5/15/50ml			
TC0401001-67	single block, Round tubes, 30 x 6mm, Φ8.3mm x 48.4mm			
TC0401001-68	single block, Round tubes, 24 x 10mm, Φ10.7mm x 48.4mm			
TC0401001-69	single block, Round tubes, 16 x 12/13mm, Φ13.9mm x 48.4mm			
TC0401001-70	single block, Round tubes, 20 x 12/13mm, Φ13.9mm x 48.4mm			
TC0401001-71	single block, Round tubes, 12 x 15/16mm, Ф17.5mm x 48.4mm			
TC0401001-72	single block, Round tubes, 8 x 20mm, Φ21mm x 48.4mm			
TC0401001-73	single block, Round tubes, 6 x 25mm, Φ26.2mm x 48.4mm			
TC0401001-74	single block, Round tubes, 4 x 35mm, Φ35mm x 47.6mm			
TC0401001-75	single block, Round tubes, 12 x 17/18mm, Ф19.1mm x 48.4mm			
TC0401001-76	single block, vials, 20 x 12mm, Φ12.7mm x 30mm			
TC0401001-77	single block, vials, 20 x 15mm, Φ15.8mm x 35mm			
TC0401001-78	single block, vials, 12 x 17mm, Φ17.8mm x 45mm			
TC0401001-79	single block, vials, 12 x 19mm, Φ19.7mm x 45mm			
TC0401001-80	single block, vials, 9 x 21mm, Φ21.7mm x 45mm			
TC0401001-81	single block, vials, 8 x 23mm, Φ23.8mm x 45mm			
TC0401001-82	single block, vials, 8 x 25mm, Φ25.8mm x 45mm			
TC0401001-83	single block, vials, 6 x 28mm, Φ28.8mm x 45mm			
TC0401001-84	single block, vials, 15 x 16mm, Φ16.4mm x 45mm			
TC0401001-85	single block, 10x8 PCR tube st rips, 0.2ml tubes			
TC0401001-86	single block, 64 PCR t ube strips, 0.2ml tubes			
TC0401001-88	single block, 2x6 cuvettes, 12.5mm			

SPECIFICATIONS:

SPECIFICATIONS.			
Model	TC0401001	TC0401002	TC0401003
Number of blocks	1	2	4
Heat output [W]	165	250	430
Heating temperature	room temp.	room temp.	room temp.
range [°C]	+5°C-150°C	+5°C-150°C	+5°C-150°C
Temperature display	LED	LED	LED
Adjustment and display	0.1	0.1	0.1
resolution [K]	0.1	0.1	0.1
Connection for ext.	DIN 12 878	DIN 12 878	DIN 12 878
temperature sensor	BIIV 12 070	DII 12 07 0	BIIV 12 070
PT 1000 variat ion;	< =± (0.15	< =± (0.15	< = ± (0.15 +
DIN EN 60751 KI. A [K]	+ 0.002 x ITI)	+ 0.002 x ITI)	0.002 x ITI)
Temperature stability within	±0.2	±0.5	±0.5
the blocks < 60°C * [±°C]		_510	2010
Temperature stability within	±0.2	±1	±1
the blocks > 60°C * [±°C]	_0		
Temperature Uniformity	0.2	0.5	0.5
< 60°C * [K]		5.5	3.5
Temperature Uniformity	0.2	0.5	1
> 60°C * [K]		- 1 - 1	·
Heating rate / Heat up t ime	5	4.5	4
with external sensor * [K/min]			
Set-up plate material	Aluminium alloy	Aluminium alloy	Aluminium alloy
Set-up plate dimensions [mm]	96x76	96x152	96x304
Fixed safety circuit [°C]	180	180	180
Timer	√	√	√
Time setting range	1min - 99h59min	1min - 99h59min	1min - 99h59min
Dimensions (W x H x D) [mm]	152 x 86 x 190	152 x 86 x 300	152 x 86 x 465
Weight [kg]	1.5	2.5	7.5
Permissible ambient	5 - 40	5 - 40	5 - 40
temperature [°C]			
Permissible relative humidity	80%	80%	80%
Protection class according	IP 21	IP 21	IP 21
to DIN EN 60529			
Voltage [V]	200-240/115	200-240/115	200-240/115
Frequency [Hz]	50/60	50/60	50/60
Power input [W]	165	250	412